

## 4. Meals: Setting ISF\_weights in /Preferences

v.4.5



**Please note that with autoISF you are in an early-dev. environment,**

where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)

### 4.1 Getting started

- 4.1.1 Reminder of pre-requisites
- 4.1.2 autoISF factors overview in typical glucose curve
- 4.1.3 Getting ready to set your autoISF\_weights

### 4.2 Initial bg rise: bgAccel\_ISF

- 4.2.1 Mimicking a HCL bolus
- 4.2.2 Widened safety restrictions
- 4.2.3 Start value for bgAccel\_ISF tuning
- 4.2.4 \_weight influence on the result
- 4.2.5 Characteristics of a well tuned-ISF
- 4.2.6 Suitability for many types of meal
- 4.2.7 Summary
- 4.2.8 Acceleration at neg. delta

### 4.3 Strong bg rise: pp\_ISF

- 4.3.1 Main function of pp\_ISF
- 4.3.2 Tuning pp\_ISF\_weight
- 4.3.3 Loop states with very little insulin need
- 4.3.4 "Quality control"

### 4.4 Sluggish rise towards a bg peak: bgBrake\_ISF

### 4.5 Sluggish rise into a bg plateau, or late plateau-ing at high bg: dura\_ISF and bg\_ISF

- 4.5.1 dura\_ISF for sluggish rise into a bg plateau
- 4.5.2 dura\_ISF for late/high bg plateaus
- 4.5.3 "One size fits all" -dura\_ISF
- 4.5.4 Options to limit iob from dura\_ISF
- 4.5.5 How dura\_ISF works
- 4.5.6 Set your dura\_ISF
- 4.5.7 Set your bg\_ISF
- 4.5.8 "Quality control" on your tuning
- 4.5.9 How "UAM" concludes insulinRequ.

### 4.6 Tuning your initial settings

### 4.7 Covering more complex scenarios

### 4.8 Profile helper

#### [Available related case studies:](#)

[Case study 4.1: Pizza](#)

[Case study 4.2: Low carb meal](#)

[Case study 4.3: Hands-off FCL on Xmas](#)

## 4.1 Getting started

Caution: This entire e-book is about Full Closed Looping (FCL).

**In case you intend to work with giving boli**, many suggestions made - notably in this [section 4](#), and in [section 2](#) – should **not** be followed. Y

You should then primarily use the autoISF Quick Guide (from <https://github.com/ga-zelle/autoISF>), and **do extra research**, on your own data. (Look at the chart in [section 4.1.2](#) - your bolus very much would change things there!).

If you shy away, for now, from FCL, please have a look into sections discussing methods with “**Meal Announcement**”, [section 07](#), and [section 13.3](#)-

### 4.1.1 Reminder of pre-requisites

This section 4. is about the core FCL aspects of autoISF. Before doing anything with this section, please make sure you have studied the preceding [sections 1](#) and [2](#) on the general pre-requisites for FCL and the developers “Quick guide” (see [section 3](#)) on the principal workings of autoISF.

Core points are briefly summarized below.

#### Start with proper “safety” settings

Before you start tuning your autoISF for FCL, **make sure you have** appropriately:

- **widened the SMB size restrictions** ([section 2.1](#)),
- **elevated** the max allowed ISF amplification via your set **autoISFmax** ( [section 2.2](#))

Both of these points are extremely important: If you set (or keep in place) narrow restrictions, this will **not** allow to see effects from a more aggressively tuned ISF. Even worse, it would cover-up too aggressive settings (e.g. on the ...\_ISF\_weights that we get to in a moment), and invariably make your loop bounce against the restriction(s).

This could even work fine, if your meal spectrum isn't broad: If, in your HCL, the **same** bolus size pretty much fitted all your meals, it could now, in FCL, be replaced by rushing, with super-aggressively modulated ISFs, into the set restrictions, to produce - with only a brief delay – the required iob that would be about equivalent to what you formerly had bolussed in your HCL.

A system that is really fit **for the variance** we all like to enjoy in our daily lives, though, would be characterized by “tolerating” pretty wide open safety restrictions\*), while having cautiously calibrated other, notably ISF modulating, parameters (as described in [sections 4.2 – 4.5](#)).

\*) Still, for safety (as also suggested in section 2.1 and 2.2), start your tuning on a middle ground, and only gradually widen SMB size and autoISFmax during your initial tuning.

80 Also make sure you have

- 81 • **set your iobTH%** (refer to [section 2.4](#) and if available [4.8](#))

82

83 Furthermore, in your early test phase, it is recommended to:

- 84 • Run the system as dummy, not connected to your body (or, on own risk, connect only as long  
85 as you watch closely)
- 86 • In AAPS preferences, switch your autoISF FCL ( = **autoISF/”Enable adaptation of ISF to  
87 glucose behavior”**) ON only during daytime hours of a meal, e.g. *11-18h*, for fully automatic  
88 “full closed loop” management *of lunches*.

89 You can do this switching manually *at 11 h and 18 h every day, or* set up an  
90 Automation that does that (see [section 3.4](#) ).

91 Take **typical but not extreme** meals. Omit sweet drinks, or drink only slowly. You are going for a  
92 “good enough” compromise, that works with your range of usual meals.

93

94 **It is wasted time to do a lot of iterations to “optimize” settings based on just 1 type of meal.**

95 See [case study 8.2](#)

96

97 Occasionally, watch the time-pattern of bg, iob (SMBs given), and insulin activity after meal start.

98 Aside from serious “mathematical” attempts to tune settings based on data from the SMB tab (or

99 the Emulator, section 10), just watching the curves develop on your AAPS main screen can, over

100 time, give you “a feel” what settings, and eating behaviors, are benign or detrimental to good %TIR  
101 performance.

102

103 [Importance of proper profile ISFs.](#)

104

105 Starters on autoISF FCL who are coming from using HCL with **dynamicISF** must be aware of the  
106 following: It is absolutely essential to build your FCL on a properly set **profile** ISFs (likely a  
107 circadian pattern over 24 hrs).

108

109 It may not apply to you, but many dynamicISF users did never bother to determine their ISFs that  
110 would maximize their HCL performance, but employ dynamicISF so to speak for going  
111 „dynamically“ through a wide range of possible ISFs, until eventually hitting a sweet spot, and the  
112 whole thing works better than before, with what they had used as a profile ISF (often only one, e.g.  
113 coming from Autotune).

114

115 The following is important to understand, as it leads straight into the core idea behind FCL with  
116 autoISF, too: It is a good idea to establish a well-running hybrid closed loop with set (non-dynamic)  
117 **ISF** (set in **profile** for each hour of the day). That ISF must be **aggressive enough** that it gets you  
118 down from a high around 200 mg/dl to target. That is roughly also the way you experimentally  
119 determined it (so I hope. See [https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-](https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/blob/HCL--settings-main-repo-(pdf)/ISF%20determination%20V.3.33.pdf)  
120 [settings/blob/HCL--settings-main-repo-\(pdf\)/ISF%20determination V.3.33.pdf](https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/blob/HCL--settings-main-repo-(pdf)/ISF%20determination%20V.3.33.pdf) ).

- 121 • Using *that strong* value also at *lower bg*, (on the way “up” , after meal start), is very positive:  
122 We do *not* want to have a *softer* acting loop when at *lower bg* (which is what dynamicISF  
123 tends to do!). autoISF will, in contrast, temporarily sharpen your ISF when, at low bg,  
124 acceleration is detected..
- 125 • On the way down from peak, towards glucose target, a somewhat too strong ISF should not  
126 hurt because much of the time your loop (well supplied with insulin before, „on the way up“) is  
127 zero temping , or at least has only a small gap to correct, from predicted bg to target bg.
- 128 • You have no business to be much above 200 mg/dl where an *even stronger ISF* may or  
129 may not help. It sure does not help at an occlusion which is about the only reason to see  
130 super high values as an experienced looper.

131

132 **Pegging ISF strength to bg level** therefore **does not make sense in FCL**. You will use the  
133 autoISF toolbox to get strongest ISF **at low** but beginning-to-rise bg,

134

135 Note: There are very much refined versions of dynamicISF that can have beneficial applications,  
136 notably in HCL ...And, yes, I know, bg levels can also correlate with insulin sensitivity. But let us not  
137 get into “chicken or egg type” discussions.

138 Rather, focus on doing a good tuning job, and use superior approaches to account for sensitivity  
139 changes in a more pro-active manner (before running into sky-high bg (or into hypos)):

140

141 Going to autoISF FCL, you absolutely must **anchor on the proper profile\_ISF**.

142

143 The profile is not “set in stone”, though. To use above terminology again: **Pegging ISF strength to**  
144 **your current insulin sensitivity** – very much like you had done all along in HCL - **does make**  
145 **sense in your FCL...**

146 (...and the fact that autoISF afterwards “anyways” often strongly modifies ISF is not a reasonable  
147 counter-argument).

148

149 There are fully automated, as well as manual ways for sensitivity adaptations of the profile ISF:

150

151 Profile ISFs can get **fully automatically adapted**, e.g. by Autosens, or by the Activity Monitor,  
152 which in autoISF we rather use ([section 6.5](#)).

153 Which of your basic related settings (in AAPS/Preferences) produce exactly which adaptation can  
154 be seen right in the top lines of your SMB tab, at each loop decision. Likewise, it can be retrieved  
155 later in logfile analysis (see Emulator, [section 10](#))

156

157 Furthermore, when using autoISF you can – as you did in the past, e.g. around exercise, or in  
158 times of illness – temporarily **manually modify** your profile ISFs

159

160 Also these effects are quantified in SMB tab and logfiles \*).

161 \*) Furthermore, the results from autoISF are explained in the SMB tab, and multiplied with (original or adjusted)  
162 profile ISF to result in the ISF (called “sens”) used in the current insulinRequired calculation

163

164 All three top buttons in AAPS (%profile switch, exercise and TT) can be freely used to adapt to  
165 changes in sensitivity/resistance, turning into a yellow color to alert you to this. (More about your  
166 “FCL cockpit” see [section 5.2.2.2](#) ).

167

168 For a start, please spend a couple of days (if not weeks) to **get your key autoISF related settings**  
169 **right**, strictly **on/for days with your normal insulin sensitivity**. This is what this [section 4](#) is  
170 about.

171

172 [Importance of starting from a well-performing Hybrid Closed Loop](#)

173

174 A **satisfying performance in Hybrid Closed Loop** mode is a pre-requisite. Expect to reproduce  
175 about the same %TIR also in your FCL, but with less daily interaction, once established.

176 Note that this refers to prior use of „vanilla“ software, without fancy „dynamic add-ons“ (such as:  
177 Autotune determined factors, dynamicISF etc). that may have introduced bias into the profile  
178 settings you bring with you into FCL now.

179

180 To reach a satisfying performance you must start from a hybrid closed loop in which you did  
181 **master your meal management well** using theoref(1) algo SMB+UAM.

182 This is a pre-requisite **to be able to forget it** ... - because the initial tuning that we now turn to  
183 demands, that you analyze your **prior best practice as your blueprint** to find appropriate settings  
184 and „teach“ your FCL to come up with the necessary job.

185

186 This is the main subject of this section 4 (finding settings for automatic meal management).

187

188

189

190

191 Do not copy settings from other FCL loopers

192  
193 When setting your parameters, don't use any given numerical example (not even as “a starting  
194 point”). Instead, anchor on **data from your *successful* Hybrid Closed Loop!**

195  
196 Most *examples given in this paper* are from an adult diabetic (Lyumjev, G6) whose insulin sensitivity  
197 can be characterized as follows: approximately 37 U TDD, thereof 13 U profile basal, at about 200g  
198 daily carbs from mainly lunch and dinner; no couch snacks or sweet drinks. The user also  
199 participates in multiple instances of daily moderate exercise such as dog walking, biking and  
200 gardening. In Hybrid Closed Loop, a typical meal bolus was 8 U that was sometimes reduced such  
201 as when activity followed the meal.

202 After seeing some more inputs from a variety of users we might put together a profile helper for  
203 some rough orientation, and for plausibility cross-checking, in [section 4.8](#)

204  
205 Importance of going step-by-step

206  
207 [Section 5](#) will explore avenues to manage “disturbances”, i.e. time blocks or situations that might  
208 demand enhanced or reduced loop aggressiveness.

209 [Section 6](#) will focus on the exercise mode, and the activity monitor.

210 In case you have a strong interest in the Activity monitor ([section 6.6](#)), you can start with  
211 calibrating that, and have it run already in the weeks when you go through sections 4 and 5.  
212 In case you use an EatingSoonTT at meal start (the author recommends to try without), note that any  
213 active TT shuts activity monitor automatically off while that TT is active.

214  
215 Resist the temptation to make use of the tools presented in sections 5 and 6 too early.

216 On your **first** setting-up and tuning attempt, **it is strongly recommended that you not “play**  
217 **around” with all ultimately available features, but stick to the sequence of steps to take.**

218  
219 Yes, “playing around” with the many extra buttons often will help find an improvement. But you  
220 likely create an instable FCL that, already at fairly standard situations, uses up some of your FCL's  
221 principal capacity to correct for disturbances. This limits what will be left to manage extreme  
222 situations.

223  
224 Caution: **Once you created a maze of little errors and counter-strategies/counter-errors, it will be**  
225 **nearly impossible to find your way** out of this mess, **towards better settings**, at any later point of time.

226  
227 Note that it is principally not easy to conclude on suitability of tuning:

228 • AutoISF comes with very many (currently 18) extra parameters, and even when employing the  
229 emulator ([sections 10](#) and [11](#)), it is quite hard to analyze their interaction.

230 One principal reason why things are difficult to analyze is, that you really can only analyze one  
231 decision, and that will put you on another bg curve. So, you can never see the full effect, along more  
232 than half an hour or so, that *any* change would really result in.

233

234 • Understandably, many loopers rather “move forward” to an over-patch for identified problems, and  
235 not bother with a more “puristic” step-by-step approach to do things right from the ground up.

236 • Aware of above sketched conundrum, the AAPS autoISF developers offer the ultimate tool to  
237 investigate “what-if”, regarding a setting change you may contemplate: A nice lady voice on your  
238 smartphone can tell you, at each loop decision, where your contemplated change would make a  
239 difference (in SMB size). This offers an opportunity to watch closely, with or without implementing  
240 that change. (It is always your spontaneous choice, whether you want to “follow the lady’s  
241 suggestion and manually add to the SMB, as suggested). More see in [Section 11.4](#)

242 But, we are getting ahead too far here. You first must find a starting point for key settings, which works  
243 reasonably for not too-challenging meals in your personal spectrum.

244 Before getting into this, let’s first have a look on how autoISF basically works. (More see in Quick guide by  
245 the developer, referenced in [section 3.2](#); or directly at <https://github.com/ga-zelle/autoISF>).

246

#### 247 [4.1.2 autoISF factors overview in typical glucose chart](#)

248

249 The core challenge of your UAM Full Closed Loop is to recognize a meal start from the glucose  
250 trend, and ramping up iob.


251

252 When setting up your autoISF Full Closed Loop, **you must set several ISF\_weight parameters**  
253 **in AAPS Preferences/OpenAPS SMB/autoISF settings.**

254

255 They relate to different stages of the typical glucose curve after starting a meal:





Note: **bg\_ISF** is not used much in FCL, as it is rather late to act on high (or low) bg level that developed. But, feel free to experiment, e.g. in case you have indications, in your data, that in the past dynamicISF was useful to manage bg extremes in some situations.

The core advantage of using autoISF withoref(1) SMB+UAM (in FCL as well as in hybrid closed loop) is that it manages the glucose curve it sees developing, **no matter what the underlying reason** is.

42 potential factors were identified (see: <https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/blob/FCL-w/autoISF/42%20factors%20influence%20bg.pdf> ), so, no wonder, that loopers who meticulously input their carbs will often *not* see the expected result.

#### 4.1.3 Getting ready to set your autoISF\_weights

Before you progress, make sure you studied the flowcharts in [section 3](#) that describe how autoISF calculates the **effective**(ly used) ISF.

**Warning:** Any bolus you „sneak in“ will severely distort the glucose curve. That could render your tuning of weights (see below) useless, and could **make your loop act in unpredictable ways** (potentially also dangerous, however, your set iobTH ([section 2.4](#)) should help here, too).



In case you feel tempted to use boli, be ready for some own extra research, and refer to [section 7](#).

After doing the prep work as outlined in [section 2](#) **you now get to calibrate your FCL to your normal meal spectrum** by initially **setting and tuning the various \_ISF\_weights**, that dynamically change with bg curve characteristics as sketched in the chart on the previous page.

**Please stay away from extremes** (regarding both, meals and exercise) **when you go through this [section 4](#)**. It is about getting a first *roughly right* set of settings, as a basis.

**Researching your standard meal patterns, and finding settings for the various \_ISF\_weights is the core job in setting up your autoISF FCL.**

Depending how varied your diet and general lifestyle are (and your expectation of %TIR you like to reach), this could be the main job at hand. However, there is much more you *could do later*, and that will be outlined in later sections 5 and 6.

Consult sometimes your SMB tab, to see how the applied effective ISF (named **sens** there) is calculated. (Example given in [section 5.4.5](#)).

## 4.2 Meal detection and managing the initial bg rise: bgAccel\_ISF

### 4.2.1 Mimicking a HCL bolus in FCL using bgAccel\_ISF

When looping without carb inputs and without giving a bolus ourselves, the first crucial setting is to set the **bgAccel\_ISF\_weight** so that SMBs are requested immediately when the loop detects an acceleration in your blood glucose (bg) that is starting to rise.

Ideally **within about 20 minutes after acceleration detection, which would be the first up to 4 SMBs, as much iob should automatically be supplied as we would have given with our bolus in hybrid closed loop.**

As the biggest principal challenge for the FCL is big **high/fast carb** meals (from within your personal “spectrum”), we start with a focus to get sufficiently big SMBs going for those.

Note, though, that in a **low carb** meal scenario, the first 4 SMBs would have to automatically result much smaller (which, after careful tuning, is possible with the same parameter settings, see e.g. [case studies 4.2 vs 4.3](#)).

Rule of thumb: Two of the first three SMBs each (in this test based on a big meal) should be about  $\frac{1}{4}$  to  $\frac{1}{3}$  the size of a bolus in your HCL „career“ (for a similar meal).

#### **Going over 1/3 could be problematic**

- if your diet contains occasional **low carb** (or brief **snacking**), it is not helpful if your settings make your loop invariably “bounce” over your iobTH (and then you would need extra snacks to balance the auto-generated iob, to prevent hypos),
- also if your **CGM quality** is sometimes unreliable, and might produce an artefact that could be mistaken for a meal start.

Be vigilant about this topic! And please do not choose the supposedly easy way, to just set safety restrictions (allowed max SMB size, or autoISFmax) so low, that your loop never can exceed  $\frac{1}{3}$ . Try to really tune the \_ISF\_weights appropriately. (Only that way, your loop can “accommodate” the entire meal spectrum, and also states of adapted general insulin sensitivity).

### 4.2.2 Widened safety restrictions

Already when tuning the bgAccel\_ISF\_weight it can become evident that safety restrictions (as discussed in [section 2](#)) must be **widened** further:

- 334 • Especially if your *profile basal* rate is very small, the **smb\_max\_range\_extention** and/or  
335 the **autoISF\_max** "must" often be increased further.
- 336 • Pay attention also to the **iobTH%** and, potentially, iobMAX
- 337 • Note that the smb\_delivery\_ratio "only" portions the insulinReq differently over the next 15  
338 minutes (see also [section 2.3](#)), and therefore is **not** a prime tuning parameter.

339 In the end you should **not set these safety limits too tight**, so "nudging" aggressiveness by  
340 another 10 or 20% from your cockpit, later, will not bounce into restrictions.

341

342 On the other hand, setting **narrower** restrictions for max allowed SMB size can also become  
343 necessary:

- 344 • Poorer CGM quality demands narrower restrictions for safety reasons.
- 345 • If you use a 1-minute CGM, please observe [section 1.4.2](#)

346

#### 347 [4.2.3 Start value for your bgAccel\\_ISF tuning](#)

348

349 bgAccel\_ISF\_weight is set default to zero in AAPS Preferences/SMB/autoISF.

350 **To start**, I would try 0.05 or **max 0.1**, and keep trying in max 0.05 steps. Soon move to 0.02 steps.  
351 From my (very limited) overview, many use around 0.2, and possibly even higher if their hourly  
352 basal rate is 0.1U or lower. ([Consult section 4.8 when available](#)). Do not be tempted to rush this  
353 setting by using large jumps in adjustments.

354 **To monitor what is happening**, and start tuning, in search of appropriate settings, you must keep  
355 (real-time) track of how autoISF uses your set bgAccel\_ISF\_weight:

356

- 357 • To do this in the **SMB tab** is possible but not very practical. You would end up making a lot  
358 of screenshots (quickly in the crucial minutes after a SMB was given, or when *you thought it*  
359 *should be* given), for later analysis.

- 360 • The superior method is to just copy **logfiles** from your phone/internal  
361 memory/AAPS/logs ...

362     ○ all zip files there

363     ○ look up how many days of data are covered there on a rolling basis, and copy out  
364     onto your PC (see [section 10.1.1](#)) before the older ones get forever lost

365 ... and analyze them at your convenience later, using the **emulator** (see [section 10](#); used  
366 e.g. in last pages of [case study 4.2](#)).

- 367 • Some emulator-based analysis is also possible within AAPS on your phone ([section-11](#)).

368 In any case, it is worth the effort to tune the **bgAccel\_ISF\_weight** in such a way that high glucose  
369 increases are already nipped in the bud, so to speak.

370

371 To summarize: In FCL, the first 3 or 4 SMBs should not be much delayed, and amount to similar  
372 iob like your “former boli in HCL”.

373 Depending on details about the carb absorption characteristics of your meal, and the performance  
374 of your CGM, also **pp\_ISF** (see 4.3) might be a fairly early contributor to getting iob up.

375

#### 376 4.2.4 How changing the \_weights influences the resulting calculated insulinRequired

377

378 To get a feel for how changing the \_weights influences the resulting calculated insulinRequired, it is  
379 best to start cautiously and just do 10 to max 20% steps up, and watch out for the effects. Doing  
380 similar step sizes should yield about similar effects each time.

381 (You can skip reading the following example calculations, unless you want to know more  
382 quantitatively how things work).

383 *Example 1: Going from bgAccel\_ISF\_weight of 0.2 to 0.16 (20% less).*

384 *If your profile\_ISF is 40 mg/dl/U and with bgAccel\_ISF\_weight = 0.20 you saw acce\_ISF*  
385 *factor of 1.31, this would ((if the acce influence dominates and is used as effective ISF))*  
386 *lead to the effectively used ISF of  $40/1.31 = 30.53$  mg/dl/U. For an intended correction by –*  
387 *10 mg/dl the insulinRequired would calculate to  $10 / 30.53 = 0.328$  U.*

388 *Now, going with a 20% reduced bgAccel\_ISF\_weight of 0.16:*

389 *acce\_ISF =  $1 + \text{bgAccel\_ISF\_weight} * \text{internalFactor}$*

390 *before  $1,31 = 1 + 0.20 * iF \Rightarrow 0.31 = 0.20 * iF \Rightarrow iF = 1,55$*

391 *after  $? = 1 + 0.16 * iF \Rightarrow ? = 1 + 0.16 * 1.55 = 1.25$*

392 *New effective ISF would be  $40 / 1.25 = 32.05$  mg/dl/U. For an intended correction by – 10*  
393 *mg/dl the insulinRequired would calculate to  $10 / 32.05 = 0.312$  U, which is 4.9% less.*

394

395 *Example 2: Going from bgAccel\_ISF\_weight of 0.2 to 0.10 (50% less; or doubling in the*  
396 *other direction).*

397 *If your profile\_ISF is 40 mg/dl/U and with bgAccel\_ISF\_weight = 0.20 you saw acce\_ISF*  
398 *factor of 1.31, this would ((if the acce influence dominates and is used as effective ISF))*  
399 *lead to the effectively used ISF of  $40/1.31 = 30.53$  mg/dl/U. For an intended correction by –*  
400 *10 mg/dl the insulinRequired would calculate to  $10 / 30.53 = 0.328$  U.*

401 *Now, going with a 50% reduced bgAccel\_ISF\_weight of 0.10:*

402 *acce\_ISF =  $1 + \text{bgAccel\_ISF\_weight} * \text{internalFactor}$*

403 *before  $1,31 = 1 + 0.20 * iF \Rightarrow 0.31 = 0.20 * iF \Rightarrow iF = 1,55$*

404 *after  $? = 1 + 0.10 * iF \Rightarrow ? = 1 + 0.10 * 1.55 = 1.155$*

New effective ISF would be  $40 / 1.155 = 34.63$  mg/dl/U. For an intended correction by – 10 mg/dl the insulinRequired would calculate to  $10 / 34.63 = 0.289$  U, which is 12 % less (going the other way, 0.328 is 13.5 % more).

Example 2 (-50%) reduces `_weight` 2.5 times lower than example 1 (-20%), and the resulting effect (-12% vs. -4.9% insulin Required) is also factor 2.5 different.

Note: “Your” internal factor “iF” might differ; for sure it is very different between the various ...\_ISF components.

**Never forget to look into how other ..ISFs play into the effective ISF (named `sens` in the SMB tab) which overall results.**

#### 4.2.5 Characteristics of a well tuned-in `bgAccel_ISF_weight`

Your starting point was to set the `bgAccel_ISF_weight` so FCL works in a rather high carb meal.

Now you must check (and potentially fine tune) so it **will not “shoot iob too high”** with the first 3 or 4 SMBs **in other meals from your spectrum:**

- For meals that are in the **lower** (!) range of the “fast **carb load**” of your cluster, the necessary insulin supply for the first two hours or so might pretty much be provided already with the first 3 or 4 SMBs

The glucose curve, at such meals, begins to flatten early in this SMB phase, so a deceleration (**braking**) follows very soon (-> [section 4.4](#)). Clearly, the first 3 SMBs, in such cases, must remain below iobTH.

- Low carb** meals are principally easiest for the FCL. However, you must secure that your `bgAccel_ISF` driven **first SMBs** remain small. This is principally possible also with a fairly aggressive `bgAccel_ISF_weight` set, because both acceleration and initial deltas are small when eating low carb. (Regarding the detected acceleration, the stakes may be high for the CGM and smoothing method you chose).
- A stage where moderate amounts of carb absorption and of insulin usage/need hold a balance could protract – at moderate bg elevation -over hours. The **dura\_ISF** might play a bigger role, then, as e.g. in the low carb example in [case study 4.2](#).

In case you run into limitations, see next sub-chapter.

#### 440 4.2.6 Suitability for many types of meal

441

442 For a **hands-off FCL**, your settings have to fit

- 443 • **in each** of your **meal times**

444 What helps here is that, *between* your daily mealtime slots, your **circadian profile ISFs**  
445 (upon which the autoISF modulations build) automatically make a differentiation (as was  
446 the case in your HCL).

- 447 • for the whole **range of your meals**. All this is principally possible, but: ...

448

449 What if you still have meals that you cannot make fit?

450

451 In extreme cases you will have to balance too high running iob with additional carbs (a late  
452 additional snack against going too low), and in the opposite case, you will have to reckon with  
453 temporarily exceeding the glucose target range, and losing some %TIR for this day.

454

455 If your meals vary very strongly, there are **avenues to ease your initial tuning job, or to**  
456 **optimize overall resulting loop performance:**

- 457 • Automations allow you to differentiate. For instance it is possible to apply different  
458 iobTH\_percent and/or different bgAccel\_ISF\_weights for meals in different **time windows**  
459 or geo locations (details see [sections 3.4](#) and [5.1](#))

460 In case you use autoISF on the Trio or iAPS platform for i-phones, you may need to use a  
461 third party automation software, or “middleware” (! [call for a case study 4.X](#) )

- 462 • you can pre-program **custom buttons for special** meal (or snack) **types**, with different  
463 underlying FCL settings (see “cockpit”, [section 5.2.2.3](#))

- 464 • You can **modulate FCL aggressiveness manually** making use of the top 3 buttons in the  
465 AAPS home screen: These turn yellow during temporary switched %profile or glucose  
466 target ([section 5.2.2.2](#))

467

468 Experimenting with the three above mentioned “avenues”, the author found:

- 469 • the last point easiest to occasionally use, and the first one hardest.
- 470 • it worth investing some effort (also using the emulator a couple of times) to iterate through  
471 the typical meal spectrum a couple of times, for finding a “good enough” set  
472 of ..\_ISF\_weights and other settings (like autoISFmax, iobTH% etc), **and not do much**  
473 **extra differentiation**. (More see in [section 5](#)).

#### 474 4.2.7 Summary on tuning for the initial SMBs via bgAccel\_ISF

475

476 **Early strong iob** also will **ease the tuning task for the subsequent phases** of the meal, because  
477 there is, then, largely zero-temping (as well known from HCL-times after your administered bolus).  
478 Also, the lower and shorter lasting the glucose peak, the lesser the hypo danger from the activity  
479 tail of SMBs given *when* glucose was „stuck“ high.

480

481 However, it is important **not too super-aggressively** tune bgAccel\_ISF\_weight up, so, regardless  
482 of the type of meal, very big SMBs invariably would result.

483

484 Rather, the rough idea should be:

- 485 • SMBs driven by bgAccel\_ISF: initial iob for **all meals**. SMB sizes vary, because accelerations  
486 and deltas vary.

487 So, at high carb meals it depends on your settings, and on the evolving bg curve, whether the  
488 first few bgAccel\_ISF driven SMBs get you already up to iobTH in high carb meals, or whether  
489 this happens in the *overlapping* next stage.

490 So, looking a bit ahead to the next chapters:

- 491 • SMBs driven by pp\_ISF: to the extent there is strong (near-linear) bg rise (at **big meals rich in**  
492 **carbs**) with big or small deltas, iob is now driven towards (and potentially over) iobTH.

493 In low carb meals this period can be extremely short, with iob remaining under iobTH (example  
494 see [case study 4.2](#))

- 495 • SMBs driven by bgBrake\_ISF, bg\_ISF, or dura\_ISF:

496 Note that *all of these* can overlap with the pp\_ISF stage. Consult the csv table output  
497 from the Emulator (example given at end of [case study 4.2](#)) as to which of the \_ISF  
498 categories drives the effectively used ISF (and what change of the ...\_ISF\_weights  
499 would change this. Consult decision flowcharts for effective\_ISF in pages 1-6 of the  
500 Quick Guide.pdf in <https://github.com/ga-zelle/autolSF> ).

501 Depending on the shape of the bg curve after the initial strong rise, and depending on  
502 insulinReq. and on iob (> iobTH?), autoISF can provide more SMBs to bring bg to target. This  
503 case applies to **low carb** meals. The dura\_ISF is also useful to manage temporary insulin  
504 resistance often observed late in **fatty** meals.

505

506 **It is worth investing effort** (following the sequence of steps in sections 01-04 of this FCL e-book)  
507 **in your initial project to establish a good set of ISF\_weights** for your meal spectrum. This will  
508 keep interventions in daily life to a minimum.



509 Unless your lifestyle, or health and body weight change radically, this should be a **one-time effort**  
510 (in your initial weeks establishing your FCL), with *no need* to fine-tune much later (see [section 8](#)).  
511

#### 512 4.2.8 Note regarding acceleration “happening again” in late part of [dropping](#) glucose 513 (*Skip, unless interested*)

514 After the peak, in the late stage of *falling* bg, the glucose curve is like an accelerating  
515 parabola again. The algorithm tries to evaluate when and at which bg level complete  
516 digestion of the meal and a bg minimum will result. Insulin required to stabilize around  
517 target bg is usually very small, and the adaptation of ISF in that stage relatively  
518 unimportant. See in your SMB tab, how, at “already falling” bg, the ISF modulation is taken  
519 back.

In version 2.2.8.2 there was a potential deficiency in situations where glucose was falling and the glucose acceleration was already positive. That meant a minimum glucose level can be extrapolated. If that happens to be less than target and expected in less than 15 minutes then there should be no strengthening of ISF as it would lower glucose even more. Therefore bgBrake\_ISF\_weight is used now instead of bgAccel\_ISF\_weight. But those situations were rare and less critical than might be expected at first sight. The reason is that in most cases the predictions ended up even below their threshold meaning SMB were disabled.

520  
521

### 522 4.3 Managing strong bg rises: pp\_ISF

523

#### 524 4.3.1 Main function of pp\_ISF in autoISF FCL

525

526 Between acceleration and deceleration there is a more or less linear further increase of bg and of  
527 insulin need.

528 • With **higher carb load** meals, or meals that come with a sweet drink, the increase will be  
529 particularly strong, and (if not already driven there by bgAccel\_ISF) now reach, and with the  
530 last “allowed” SMB exceed, the valid iobTH.

531 • With **low carb** meals, there is only a very un-pronounced (short, with weak deltas) “pp\_ISF  
532 phase”. (Example see end of [case study 4.2](#)).

533

534 autoISF should now “fight” this with the help of the post-prandial ISF, set via **pp\_ISF\_weight**, after  
535 you have set your bgAccel\_ISF\_weight.

536

#### 537 4.3.2 Tuning pp\_ISF\_weight

538

539 To tune-in your **pp\_ISF\_weight**, please do this with a really high carb meal (from within your  
540 typical meal spectrum) *after* you have set a halfway suitable (not too aggressive)  
541 bgAccel\_ISF\_weight.

Note that if you rush into pp\_ISF tuning while “still having a too aggressive bgAccel\_ISF”, the latter is covering up the requirement you now really want to calibrate for in pp\_ISF!

So, at a meals in the upper spectrum of your carb load, carefully begin with a starting value for *pp\_ISF\_weight* of 0.005. Observe the reactions and check the SMB-tab before you increase it cautiously for the next days.

Best practice is to analyze the emulator tables (discussed in [section 10](#), and example given in the pizza [case study 4.1](#))

#### 4.3.3 Loop states with very little insulin need (*iob* > *iobTH*, or 0 %TBR)

Normally (except for very low carb meals) the SMBs triggered by *bgAccel\_ISF\_weight* and *pp\_ISF\_weight* should be sufficient to reach and slightly exceed the ***iobTH*** (see [section 2.4](#)) so all *the other* autoISF parameters are relatively unimportant for now.

A reason why this can work at all, also for quite a variety of meals, lies in the fact that there is an hourly carb absorption limit of about 30g/h

(Reference: Dana

Lewis: <https://github.com/danamlewis/artificialpancreasbook/blob/master/8.-tips-and-tricks-for-real-life-with-an-aps.md#heres-the-detailed-explanation-of-what-we-learned>. (That limit can be lower, e.g. with gastroparesis or certain medications, but that would make things even easier)

So while meals might wildly vary in composition and size: What is digested, and needs insulin in the first ~90 minutes (when FCL tries to catch up with insulin need and differs strongly from HCL, with *bgAccel\_ISF* and *pp\_ISF* in the leading role), will be relatively close...for meals with similar *initial* glucose acceleration and rises, anyways...

The others, **low carb** with much slower initial acceleration and rise, are easy recognized as different by the loop, see [section 4.4](#) that follows.

Depending on the type of meal and “aggressiveness” of your *bgAccel\_ISF\_weight* and *pp\_ISF\_weight* tuning, the *iob* will already be so high that, in the phase of decelerated glucose rise towards the peak (the “last part of the rise”), no more *insulinRequired* is seen by the loop.

Therefore the ***bgBrake\_ISF\_weight*** is often unimportant in meals with a relevant carb content.

For potential relevance in low carb meals, see [section 4.4](#).

#### 581 4.3.4 “Quality control” on how well your tuning can replace your former HCL bolussing

582

583 Warning: **Occasionally consult the SMB tab to see how your settings really work.**

584 A setting (...ISF\_weight) that is actually set too aggressive might be masked.

585 **Tuning only works if** the effects of the settings being tuned are **not** unintentionally **limited by**  
586 **other** (e.g. „safety“) **settings**.

587

588 Also, **always look at two or three *different* meals** before deciding whether a tuning “fits” („good  
589 enough“ for each of them). You probably will have to iterate back and forth doing this for two or  
590 three different kinds of meals ...

591 • [Case Study 4.1](#) (Pizza Meal) contains, towards the end, an example how you can go about  
592 tuning the \_weights for various \_ISF factors of autoISF.

593 • [Case Study 8.2](#) shows that it is **not** worth it to seek “optimized” settings based on just one  
594 meal.

595 ... until you find *one* good enough set of settings *for all* of them. Do not rush this, establishing a  
596 solid foundation will be well worth your time.

597

598 **The following sections will deal with similar issues like you were facing in HCL after your**  
599 **given bolus lost much of its power, and SMBs were needed for the “eCarbs”.**

600

#### 601 4.4 Sluggish rise towards a bg peak: bgBrake\_ISF

602

603 At a **low carb** meal, or an attempt at doing a **weight reduction diet**, (and probably also with  
604 gastroparesis, or if you take one of these novel GLP-1 drugs that slow meal absorption -

605 **Somebody, please supply a case study!** - the glucose goes up only sluggishly, and iobTH should  
606 not be reached at all.

607 In case you *exclusively* do very slow absorbing meals, you could of course also adjust your iobTH  
608 setting low enough to suit your *uniform* situation.

609

610 Acceleration, and the phase of strong glucose rise, are quickly over at slow-absorbing meals, and  
611 there can be:

612 • a decelerating bulge of insulin action that projects over an hour or longer. This is where the  
613 importance of the **bgBrake\_ISF** can come in.

614 • a bg curve that hovers for an hour or longer around an elevated bg level, because  
615 additionally absorbed carbs, and consumption of the moderate SMBs delivered, tend to  
616 keep a balance for a while. **Dura\_ISF** can deal with this (see next chapter). An example for  
617 this is given in [Case study 4.2](#).

618 Note that in some data outputs (e.g. the csv/xls tables coming from the Emulator, e.g. in Case study 4.2,  
619 big table at the end there), you will see only “**acce\_ISF**” results.

- 620 • In case of positive acceleration, these are driven by the bgAccel\_ISF\_weight setting, and results  
621 are >1.
- 622 • **In case of negative acceleration** (decelerating rise), **bgBrake\_ISF\_weight** is **applied**, , and  
623 results are < 1. (Example see in graph in [section 10.3.3.3](#) ).

624

625 In full closed loop, the bgBrake\_ISF\_weight is often only about half as large as the  
626 bgAccel\_ISF\_weight (but that would also depend on your personal diet pattern and  
627 eating/digestion speed). Also here, one should approach the tuning gradually, increasing the  
628 weight coming from small values.

629

630 Please observe that **tuning bgBrake\_ISF\_weight must strictly be done with types of meals for**  
631 **which there is insulin need at de-celerating but still rising bg.**

632 bgBrake\_ISF is totally irrelevant for hi carb meals where your loop shot over iobTH already  
633 by the time your rising towards the bg peak slows down!

634 Likewise, if your initial bgAccel\_weight is set so strong that your first SMBs catapult you  
635 over the iobTH, no matter what type of meal: Then you must **first** find a reasonable setting  
636 for this parameter, one that works “good enough” to control your carb loaded meals, and  
637 then see whether there is “room” (and need) for milder loop response at low carb meals.

638

639 In case you cannot quite get all the ISF\_weights “right” so the occasional low carb meal will not get  
640 over-treated: Avenues to adapt your loop aggressiveness are discussed in [section 5](#).

641 For instance you will be able to (if needed):

- 642
- 643 • use a temp. reduced %profile
- 644 • temp. lower iobTH or bgAccel\_ISF\_weight
- 645 • construct for yourself a “DIY cockpit” with an extra “snack” or “low carb” button with an  
646 underlying suitable Automation

647

648 In the **late stage of still rising (!) glucose**, the Full Closed Loop typically sharply reduces  
649 SMBs already because it is “painfully aware” of the following principal conflict:

650

- 651 • iob (like formerly given in HCL via your bolus) must go high quickly, in order to limit the high

- However, if there is too much insulin in the system, a **hypoglycemia can happen later** within the DIA time window, because the loop can, later, only correct to a very limited extent (namely, only to the extent that it can set basal to zero).

Therefore, the core problem is that the Full Closed Loop must build up iob very quickly, **but not too much**, in the initial phase of a meal, and high bg values (out of range, >180 mg/dl) can not always be avoided.

## 4.5 Sluggish rise into a bg plateau, or late plateauing at high bg: dura\_ISF and bg\_ISF

Depending how your personal diet spectrum looks, you need to tune-in your dura\_ISF primarily with large hi-FPU meals, and/or for meals at the low carb end of your diet

### 4.5.1 dura\_ISF for sluggish rise into a bg plateau

A (in that case, often not very high) plateau can form in **low carb meals**, when, basically, carb and insulin “burn rates” might keep a balance over an hour or longer, requiring occasional moderate size SMBs.(See an example in [case study 4.2](#))

### 4.5.2 dura\_ISF for late/high bg plateaus

With **large or high fat/protein meals**, often a long high bg plateau is encountered (sometimes associated with 2nd “late, long stretched hill” forming for this, in the bg curve).

For such situations, autoISF features the modulation of ISF depending on bg level and on duration of **plateau** formation.

### 4.5.3 “One size fits all” -dura\_ISF

Absolute “pros” could primarily calibrate their dura\_ISF for low carb.

Dura\_ISF has in-built amplification at higher bg levels. So, effects will automatically be boosted in case much higher plateaus develop after greasy feasts.

Should that not per se be sufficient, there is more the DIY “pro” can do:

- by adding an Automation that gives an extra boost “against” the temporary insulin resistance associated with fats (via increasing the baseline, in terms of a temp.130% profile switch, for instance. Compare at:

<https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#stagnation-at-high-bg-values> ),

- or by making additional use of the bg\_ISF (or dynamicISF) (-> Tune it in parallel.)

The author's preference would be to go via Automation, but only in case the in-built differentiation via bg level make it necessary.

#### 4.5.4 Options to limit iob delivered from dura\_ISF


Rather than relying on your initial tuning to keep you safe from hypos also in the future, there are some extra precautions you could take. Some were discussed in Discord or in dev circles, regarding what could be done:

- 1) To limit the danger of going low, it can make sense to design an **Automation** which pauses the delivery of more insulin.

This one was suggested by Alex999

If a glucose plateau built under 140 mg/dl, do not treat via dura\_ISF (because the defined Action is to set an elevated TT to a level that will not require more correction insulin.

An alternative Action would be to set, near the actual glucose target, an odd-numbered TT (which blocks any SMB be given, while valid).



- 2) In an autoISF update, the **duration** in which iob is added up could be **capped** after max. 1.5 hours of any “stubborn high”.
- 3) Instead of 2), or additionally, the total **iob accrued in that “dura phase”** could be **capped** by a new related safety setting. It would probably be anchored on iobTH, and could also become a tuneable setting, maybe even a new parameter useable in Automations, too.

#### 4.5.5 How dura\_ISF works

Conditions for dura\_ISF to become active:

- 1) glucose is varying within a **+/- 5% interval** only;

- 2) the average glucose (*dura\_ISF\_average*) within that interval is **above target**;
- 3) this situation lasted **at least for** the last **10 minutes** (*dura\_ISF\_minutes*)

Effect:

- 4) The strengthening of ISF is stronger the longer the situation lasts, and **the higher** the average glucose is above target:
- 5) This can be individually tuned by the **duraISF\_weight to automatically manage** high plateaus in bg values.

The formula looks like (this, and more, see page 6 of the Quick guide at [Github/ga-zelle/autoISF](https://github.com/ga-zelle/autoISF)):

$$\text{dura\_ISF} = 1 + \frac{\text{avg05-target\_bg}}{\text{target\_bg}} * \frac{\text{dura05}}{60} * \text{dura\_ISF\_weight}$$

where

avg05 = *dura\_ISF\_average*  
dura05 = *dura\_ISF\_minutes*

*Off topic:* *dura\_ISF* is also very useful in Hybrid Closed Loop. It can be used to elegantly manage, fully automatically, a temporary insulin resistance from fatty acids. Please refer to other papers for details (for instance, section „Late stage of meals“ of “Meal Management Basics”, available here: <https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings> ).

#### 4.5.6 Set your *dura\_ISF*

Set a **start value of 0.2** for your *dura\_ISF\_weight*, and increase only cautiously with an eye on hypo prevention 2-3 hours later.

Caution: Fine tuning this parameter only makes sense **after** you tuned your *bgAccel\_ISF* and *pp\_ISF* well (so your thin yellow insulin activity curve shifts as *far to the left*, towards meal start, as *possible*, which will lower bg peaks and ease the job for *dura\_ISF*).

#### 4.5.7 Set your *bg\_ISF*

Since in Full Closed Loop we make our loop give us the maximum SMB size it can, at the beginning of a rise, it is crucial to **resist the temptation to continue** with a particularly **strong ISF** in the meal phase with the **highest glucose** values .

This is a reason why in Full Closed Loop we do not make much use of the ***bg\_ISF*** component of *autoISF*.



- Wanting to get most of our insulin from SMBs delivered at fairly low (but beginning-to-rise) bg implies that we do **not** make ISF weaker at low bg. Under preferences/OpenAPS SMB/autoISF/bg\_ISF settings you could set **lower ISF\_range\_weight** = 0.0.

If you want to analyze in your data, whether you might benefit from a milder ISF at low bg values (e.g. if you often go below target after correction of only mildly elevated bg in the preceding hours), you may want to try lower ISF\_range\_weight = 0.1 or 0.2. Study the effects from bgISF, and increase, or decrease, the bgISF\_weight to fine tune the sought-after affect.

- The **higher\_ISF\_range\_weight** is used when bg is above target, It then strengthens ISF the more the higher the set weight is. 0 disables this contribution, i.e. ISF is constant in the whole range above target.

In FCL, this factor should be fairly irrelevant: Near glucose peak, zero-tempering usually prevails anyway, so the settings we try might often not be used really by the loop. Very likely, you can live with setting the weight to = 0.0 here, too.

If you want to analyze in your data, whether you might benefit from a stronger ISF at high bg values (e.g. if you often remain above target after correction of elevated bg in the preceding hours), you may want to try higher ISF\_range\_weight = 0.1 or 0.2. Study the effects from bg\_ISF, and increase, or decrease, the higher\_ISF\_range\_weight to fine tune the sought-after affect.

#### 4.5.8 “Quality control” on your tuning for the later half of your meal time

The later stages of meal management (both, in HCL and in FCL) struggle with the problem that there is a **hypo danger** from the “tail” of insulin activity from earlier SMBs that were needed to fight high bg or plateaus associated with temporary insulin resistance.

Once your bg sits high, neither you, nor a hybrid closed loop with all the carb info, nor your FCL can work wonders.

Very important:

- Iterate between **2 or 3 kinds of meals** (from your typical spectrum) to find **one** set of settings that works *good-enough for all*. That should be possible.
- If you can’t make it work for certain meal types, see [sections 4.7](#) and [5](#). what you can do then.

794 Observe hypo trends after meals, and

795 • resist the temptation to elevate the **dura\_ISF\_weight** very high.

796 • try to stay away from **bg\_ISF** or dynamicISF in Full Closed Loop:

797     ○ In FCL you probably can afford to shut bg\_ISF entirely off via setting both related

798     \_weights to 0.0.

799     ○ At least be careful, use small ISF\_range\_weights and check whether you are happy

800     with the contributions to the effectively used ISFs

801     ○ *Off topic:* If, coming from dynamicISF usage, you stay in **Hybrid** Closed Loop, but now with

802     autoISF, you probably can use the bg\_ISF parameter with higher \_weights to emulate what

803     you like to replicate from your dynamicISF experience.

804

805 bg highs will take time to resolve.

806 Interestingly, an after-dinner walk can work wonders sometimes (take glucose tablets along).

807 Zero-temping and too tightly set safety limits can be **stumbling blocks in your tuning** project

808 Investigating effects of set ISF\_weights is not really possible in periods of zero-temping.

809 Too aggressive settings might not come into play most of the time.

810 However, some *other* time they might come into play, and *then* produce a hypo 1-2 hours later.

811

812 Therefore, **carefully study the SMB tab** (or better yet, do an **emulator** based analysis, see

813 [sections 10-11](#)) to see

814 • what the selected weights **would** do, **if** there was **no zero-temping** at the time, and

815 • whether you bump into a **set limit** already (if your bgAccel\_ISF\_weight makes you exceed

816 allowed max. SMB size, then further tuning your settings only makes sense with either

817 allowing bigger SMBs, or limiting bgAccel\_ISF\_weight to a lower number at which you will

818 not frequently bounce into the SMB limit)

819 • at which **other** times (rather than the one you currently look at and try to improve) that

820 selected setting might backfire

821

822

#### 823 4.5.9 How your “UAM” concludes insulin need for your un-declared carbs former 4.5.3 quoted elsewhere!

824

825 The UAM Full Closed Loop doesn't get any information from you as to how many grams of carbs  
826 will be there, for absorption.

827

#### 828 Looking back

829

830 For the recent 5-minute segments, the UAM oref(1) loops can precisely calculate how many grams  
831 of carb “must have been digested” based on the CGM values seen, and your IC and ISF profile  
832 parameters

833 For more detail see chapter 1.2 on how dynamic carb absorption is calculated, in “IC (carb  
834 ratio)...pdf” at: <https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings>.

835

#### 836 Looking into the next minutes, hour

837

838 However, *here* we worry about the *late meal stage*, and our FCL has gotten no information from us  
839 about how many grams **in total** were eaten, and certainly we do not bother to give eCarbs with  
840 estimated **absorption times** (that are so essential in iOS Loop).

841

842 So, in FCL you leave your loop without knowledge when your steady-state max carb absorption  
843 phase...

844     ○ the earlier mentioned 30g/h, or

845     ○ with gastroparesis, or if on GLP-1 drug treatment, probably on a lower g/h level

846     ○ sometimes prolonged (“faked”) by a brief episode of insulin resistance to fats

847 ...might end. Nor, whether extra carbs were added, later, or “FPUs are lurking”.

848

849 The FCL now needs to provide desired amounts of insulin, while facing potentially induced hypo  
850 danger later (considering the DIA of all insulin in use).

851

852 Fortunately, the UAM Full Closed Loop is *not completely clueless* regarding how carb absorption  
853 *will continue*:

854

855 It will work with a **prediction** of *further* carb absorption, building on the **carb deviations**

856 (=calculation of how much got absorbed in the *past* 5 minute segments), and phase out further

857 *expected* carb decay in the course of the next 1 to max 3 hours.

858

859

860 For more detail see

- 861 • <https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand-determine-basal.html#understanding-the-basic-logic-written-version>
- 863 • *or do a real-time study* with (screenshots from) your SMB tab info.

864

## 865 Discussion

866

867 **This UAM prediction** about further carb absorption can be worse, but **can also be better than a**  
868 **prediction based on the user's „e-Carb“ input in Hybrid Closed Loop.**

869

870 In any case, and even when having perfect knowledge about how exactly the carbs fade out in the  
871 next hours, there would still be a principal problem for the loop: Heavy insulin „fire“ against highs  
872 will not work immediately (depending on the insulin's time-to-peak), and notably it comes with a  
873 significant hypo danger from the „tail“ of insulin activity.

874 A big bolus, or also a series of boli, will rarely work exactly for several hours matching the  
875 absorption of carbs (from what, how much and how fast the user ate).

876

877 *Off topic closing remark:* With meticulous attention to all carb-related profile parameters, and daily  
878 inputs on amounts and absorption times, plus some “mindfulness” as to which diet habits disturb the  
879 possible balance, there are “pro” loopers (notably on iOS Loop) who achieve impressive %TIR  
880 performance. – The author consciously chose the other way, to put substantial effort into a  
881 personalized upfront system calibration, and work with aoref(1) algo system that allows (nearly)  
882 every day hands-off FCL.

883

## 884 4.6 Tuning your initial settings

885

886 Be pro-active: **The earlier large SMBs come** (driven by bgAccel\_ISF and pp\_ISF) ...

887 Note: Also your CGM smoothing may play a role here, that you may want to look into !

888 ...the **less high** the overall increase in BG will be, and (provided you set a proper iobTH)  
889 the **lesser** the **risk** will be **for a hypo** after the meal.

890

891 Therefore, **put most of your FCL tuning effort into determining suitable weights for**  
892 **bgAccel\_ and for pp\_ISF, and for finding a suitable iobTH\_percent.**

893

894 Low carbers probably should pay more attention on **dura\_ISF**, besides seeing to it that  
895 bgAccel\_ISF is not too aggressive (see above, [section 4.2.5](#) and [case study 4.2](#)).

896

897 Later, your “FCL cockpit” will give you access to **temporarily modulate** these essential  
898 parameters (see [section 5.2.](#)), providing you an opportunity

- 899 • in your tuning phase, for more research on the fly, so to speak
- 900 • everyday, for temp. adaptations to altered insulin sensitivity, or to special  
901 disturbances (should you, occasionally, see a need).

902

903 After you tuned your **initial settings** well, there should rarely arise a need for “fine tuning” later,  
904 see [section 8](#) and [case study 8.2!](#)

905

906 The experience of the author is that it is possible to tune the above mentioned weights for very  
907 different meals in such a way that the glucose almost always remains acceptably in range.

908

909 However, if you come to the conclusion that **differentiated settings** (for different meals or meal  
910 time clusters) would be easier to establish, and/or work better for you, the following sections  
911 suggest many options you could try and use.

912

## 913 [4.7 Maneuvering through more complex scenarios](#)

914

915 You now can move on, to accommodate more complex scenarios.

916

917 To deal with **different disturbances than presented by the meal spectrum you were**  
918 **calibrating for**, there will be **temporary modulations** of your FCL possible. as was already  
919 touched on in [section 4.2.6](#)

920

921 Depending

- 922 ○ how satisfied you are with your initially reached result, or which **more extreme meals**  
923 (smaller? faster/slower carbs? totally different fat/protein content?) you would like your FCL to  
924 manage as well, or
- 925 ○ whether you seek **temporary** adjustments that **make your FCL act more aggressive, or**  
926 **softer**

927 you have a variety of options to deal with that, and this will be the topic in [section 5.](#)

928

929 It is suggested to do **major exercise** still *in your hybrid closed loop* setting, *until* you have your  
930 FCL up and running for meals on normal days with no or only moderate exercise.

931

932 Later, implement extras as discussed in [section 6](#) to fully implement your FCL.

933 While FCL is about **fully automatic cruising**, your **AAPS main screen** will serve you as your  
934 “FCL cockpit” to check how everything is running, and to aid your loop maneuvering through some  
935 special disturbances.

936

937 For a “deeper look” into a loop decision, consult sometimes the **SMB tab** (example given in [section](#)  
938 [5.4.5](#)). There you can see how the autoISF modulation of ISF is overall applied to arrive at the  
939 actually used **effective ISF (“sens”)**:

- 940 • In the SMB tab, *above the “start autoISF.” line*, the profile ISF is given (“ISF unchanged”),  
941 eventually with adaptation by activity monitor (“adjusting ...ISF from ... to ..”) or by a TT  
942 (“adjusting ...ISF from ... to ..”) or by a %temp. profile set (still called “ISF  
943 unchanged” then, meaning unchanged yet by autoISF).
- 944 • *Then* follows the autoISF section, explaining in detail how the recently encountered bg  
945 curve characteristics suggest adaptations, and what overall the conclusion is (“final ISF  
946 factor”, calculated following the flowcharts as explained in detail in section 03.).
- 947 • *Below the autoISF section*, the effective ISF (sens) results from dividing the (unchanged or  
948 adapted) ISF *prior to* “start autoISF”, with the determined “final ISF factor” at the end of the  
949 autoISF section of the SMB tab.

950

951

## 952 4.8 Profile helper

953

954

955 xls based tool has been under development

956

957

- needs more user data
- we have second thoughts, because of “Do not copy settings from other FCL loopers” (see  
959 section 4.1.1 )

960

- chapter with link to xls tool *might* follow later