

11. Emulator on your AAPS Smartphone

V15

The Emulator on the PC discussed above is very good for making your initial tuning for a meal spectrum, when weighing different effects over the entire course of time after each meal, for a variety of your meals.

A very useful additional tool, is the Emulator running on the AAPS loop smartphone

Installation guide and basic use instructions see here: https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.2_ai3.0_dev/Documentation%20in%20English/A%20%20README.md

11.1 (1) autoISF **tab** in AAPS main screen

or ALTERNATIVE: -> (2)

(Note that these features are not yet available.

For now, even the simplified emulator analysis as proposed in this [section 11.1](#) must be done on the PC, following instructions like e.g. https://github.com/ga-zelle/APS-what-if/blob/A3.1.0.3_ai2.2.8.1/Documentation%20in%20English/How-to-create-the-autoISF-factor-plot.pdf)

The **autoISF tab** in the **AAPS main screen** gives access for the past 3-4 hours, of how the various ISFs came into effect.

Hence the table gives clues about which of the weights to modify for any desired effect on SMB sizes

add ISF bg; delta and/or pp

| ISF | ISF | ISF | ISF | ISF | ISF | Ins.Requ | SMB | TBR |
|------|------|------|-------|------|------|----------|------|--------|
| prof | dura | high | delta | acce | emul | emul | emul | emul |
| 43 | 43 | 43 | 43 | 39,9 | 39,9 | -0,03 | 0 | 0,15 |
| 43 | 43 | 43 | 43 | 44,8 | 43 | 0 | 0 | 0,165 |
| 43 | 42,2 | 43 | 43 | 42,7 | 41 | -0,02 | 0 | 0,06 |
| 43 | 42,2 | 43 | 43 | 43 | 40,6 | 0 | 0 | 0,055 |
| 43 | 42,2 | 43 | 43 | 42,7 | 41 | 0,07 | 0 | 0,7175 |
| 43 | 42,2 | 43 | 43 | 42,3 | 40,6 | 0,12 | 0 | 0,715 |
| 43 | 41,7 | 43 | 43 | 39,6 | 39,6 | 0,2 | 0,1 | 0,983 |
| 43 | 42,2 | 43 | 43 | 41,5 | 40,6 | 0 | 0 | 0,483 |
| 43 | 41,7 | 43 | 43 | 41,4 | 40,6 | 0,02 | 0 | 0,623 |
| 43 | 41,7 | 43 | 26,9 | 41,2 | 26,9 | 0,63 | 0,4 | 1,843 |
| 43 | 43 | 43 | 14,3 | 23,3 | 21,5 | 4,09 | 2,6 | 5,5 |
| 43 | 43 | 43 | 12,6 | 24,8 | 21,5 | 2,33 | 1,5 | 0 |
| 41 | 41 | 41 | 8,9 | 25,1 | 20,5 | 3,85 | 2,5 | 0 |
| 41 | 41 | 41 | 12,1 | 33,7 | 20,5 | 0,83 | 0,5 | 0 |
| 41 | 41 | 41 | 8,2 | 32,1 | 20,5 | 1,654 | 1 | 0 |
| 41 | 41 | 41 | 18,6 | 39,9 | 20,5 | -1,12 | 0 | 0 |
| 41 | 41 | 41 | 13,7 | 49,3 | 24,7 | -1,46 | 0 | 0 |
| 41 | 41 | 41 | 12,8 | 29,4 | 20,5 | 1,682 | 1 | 0 |
| 41 | 41 | 41 | 17,1 | 46,4 | 23,2 | 0,3 | 0,1 | 0 |
| 41 | 41 | 41 | 14,6 | 44 | 22 | 1,422 | 0,9 | 0 |
| 41 | 41 | 41 | 14,6 | 44 | 22 | 1,447 | 0,9 | 0 |
| 41 | 41 | 41 | 15,8 | 38 | 20,5 | 0,975 | 0,6 | 0 |
| 41 | 36,6 | 41 | 22,8 | 51 | 28,4 | 0,852 | 0,5 | 0 |
| 41 | 34,7 | 41 | 41 | 60,3 | 51,1 | -0,8 | 0 | 0 |

<- add selected weights in middle columns

Table scrolls through display (landscape oriented)

Tbd DEV:


(1) give ISFs or factors on profile ISF in middle

| MEZ | bg | iob |
|-------|-----|------|
| 12:03 | 98 | 0,71 |
| 12:08 | 95 | 0,59 |
| 12:13 | 93 | 0,47 |
| 12:18 | 90 | 0,36 |
| 12:23 | 89 | 0,27 |
| 12:27 | 88 | 0,23 |
| 12:33 | 88 | 0,2 |
| 12:38 | 86 | 0,29 |
| 12:43 | 85 | 0,25 |
| 12:48 | 88 | 0,22 |
| 12:53 | 98 | 0,69 |
| 12:58 | 110 | 3,42 |
| 13:02 | 128 | 4,79 |
| 13:08 | 140 | 7,09 |
| 13:12 | 160 | 7,35 |
| 13:18 | 166 | 8,04 |
| 13:23 | 176 | 7,7 |
| 13:28 | 187 | 7,32 |
| 13:33 | 194 | 7,91 |
| 13:38 | 203 | 7,58 |
| 13:38 | 203 | 7,55 |
| 13:43 | 211 | 8,03 |
| 13:48 | 215 | 8,15 |
| 13:53 | 215 | 8,16 |
| 13:53 | 215 | 8,11 |
| 13:57 | 211 | 7,66 |
| 14:03 | 206 | 7,15 |
| 14:07 | 201 | 6,66 |

25

26

27 Additionally, key info from this table is presented in graphical form like this **example** (by ga-
28 zelle; LC - try to **get picture fitting table** above - **with 3 different _ISFs contributing** strongly
29 to SMB -over time)



11.1 (2) autoISF charts in AAPS main screen (ALTERNATIVE)

autoISF-related charts can be activated below the glucose chart in the AAPS home screen (*below IOB, SENS etc*) to give access to how the various ISFs came into effect in the past 3-4 hours,

Hence the graph / table / data similar to those suggested in the preceding [section 11.1 \(1\)](#) gives clues about which of the weights to modify for any desired effect on SMB sizes

- (can we add pictures from prototype ga-zelle ?)

11.2 Real-time checking a „what-if“ question using speech synthesis

The emulator on your smartphone can help clarify a "what if..." question.

In running the Emulator on the phone, one can define in the .vdf file of the Emulator, which setting one would like to be differently aggressive than in the active AAPS.

At times when this different setting would have resulted in smaller or greater SMB insulin delivery, the notification is reported **via speech synthesis**, and you can assess the situation in real-time yourself.

53 If for instance a suggested extra or bigger SMB makes sense, you can add this portion
54 manually*and in time observe, whether this bolus was OK and you should switch to the
55 different setting you were investigating.

56

57 **In Full Closed Loop, you don't need any buttons at the bottom of the AAPS main screen. But for such*
58 *test phases it is practical to reinstall the insulin button at the bottom of the AAPS main screen*
59 *(Preferences/Overview/Buttons/Insulin -> ON).*

60

61 After a few iterations, you'll get a feel for whether you want to incorporate this tightening
62 into the active AAPS.

63

64 Warning: Your settings must always work for a variety of meals. Do not put too
65 much effort into optimizing one situation!

66

67

68 *(add link to more info AND more info here, and / or a case study)*