

5. Time Blocks in which No Aggressive autoISF Should Run V16

Once the initial tuning according to [section 4](#) is done, you are ready to use autoISF for your automated meal management.

You will have two major other challenges to manage:

- recognize and manage (partial) occlusions, or other technical (CGM or BT related) obstacles
- deal with times when the loop should be set "milder" as a precaution.

How big this challenge is depends very much on your individual lifestyle. [Sections 3 and 4](#) discuss this in more detail.

In order to run the loop fully automatically around the clock, the times outside the meal blocks must also be precisely analyzed, and solutions to problems must be sought.

5.1 Manage Time Blocks of Different Aggressiveness w/ Automations and/or setting odd TT

Personalized Automations tailor the loop exactly to your data so **fully automated handling** of time blocks with different aggressiveness of the loop can be made.

In setting up your FCL, you therfore now have another difficult and time-consuming job at hand, to define **automated** solutions for any of your „other“ situations, outside of meal management.

5.1.1 Using Automations to reduce FCL aggressiveness outside of meal times

A widely used ACTION that strongly modifies how fast your FCL can add more iob is setting an odd-numbered temp. glucose target which makes the loop operate without giving any SMBs (%TBR modulation only).

Ensure the even/odd logic in the settings is toggled on in Preferences> openAPS SMB> autoISF settings> smb delivery settings>: "Enable alternative activation of SMB depending on TempTarget" ON as well as "Enable...depending on profile target" ON.

So, from patterns you find in YOUR data, at times where you want your loop act differently, you need to carve out CONDITIONS that describe the respective situations (and either for how long it typically lasts, or at which *other* CONDITIONS you want your loop get back to default FCL operation).


From, autoISF 3.0 onwards, also the following parameters are provided as CONDITION and/or as ACTION for defining YOUR Automations:

- Enable ISF weights / Disable ISF weights => Allows temp. ON/OFF for the key ISF modulation parts of autoISF
- Trigger/set iobTH percent => Keeps default aggressiveness, but only until a modified iob threshold is surpassed
- Trigger/set bgAccel_ISF_weight => Modifies the default aggressiveness

Recommended solution for nights

The recommended solution for nights is to set an odd profile target that prevents SMBs.

If you do have nights that would benefit from a couple of SMBs (to treat temp. highs from a late fatty pizza, raclette and such): Define suitable Automations like the two „night“ ones in this list:



Do not underestimate the „trickyness“ of getting these Automations „right“.

Night data (with your thought-out Automations in place) need to be analyzed to see

- whether the bg and iob limits defined in the given example work sensibly
- whether the TT duration is chosen appropriately
- Swapping the sequence in which the automations appear in the automation list would also lead to different SMB impacts.

5.1.2 Using Automations to focus (or differentiate) FCL aggressiveness to (in) meal time windows

If, aside from meal management, you were rather happy in hybrid closed loop, you could continue to run in that mode, and just focus your new autoISF FCL on management of meals (on all meals, or only on a sub-set of them, like only dinners – which might make sense especially in your initial transitioning phase).

For this, you define Automations

- that set **meal time windows** in which autoISF gets fully turned on
- or: that turn off autoISF (or just bgAccel_ISF off) in time windows in which surely no meal occurs.

72 Other early DEV AAPS variants (see [section 13.3](#)) all work with meal-time windows. The window is
73 either set by time of day in the settings, or it always must be „set“ by the user via giving a
74 mandatory small pre-bolus before any meal starts. **Outside** of these time windows, these loops
75 then runs with less aggressive SMBs like oref(1) SMB+UAM in AAPS Master.
76 This mode is not really FCL, but an advance over traditional HCL that often achieves satisfying
77 degrees of automation and performance.
78 The term **Meal Announcement** (MA) is often used to label this closed looping mode.

79

80 5.1.3 Activity Monitor

81 If you choose to make use of your smartphone's **stepcounter**, you can (automatically)
82 adjust insulin sensitivity ratio to activity level in the past minutes to one hour time frame.

83

84 This is another little tuning opportunity, in which you study your body's response to light
85 exercise (like walking) or to not moving at all (like desk, couch), and select appropriate
86 settings which, in the future, will automatically adjust insulin delivery to suit activity state of
87 the past minutes (up to 1 hour).(AAPS Preferences/OpenAPS SMB/Activity modifies
88 sensitivity/ -> set two scaling factors.)

89

90 This autoISF feature (new since V.3.0) is much quicker responding than Autosens or
91 dynamicISF to adjust insulin sensitivity to your current „lifestyle state“.

92

93 More see [sections 3.5](#) and [6.5](#)

94

95 5.1.4 Completely hands-off Full Closed Loop?

96

97 Remaining 24/7 in a „hands-off“ FCL can be a realistic goal with autoISF 3.0 if special challenges -
98 as discussed in [section 5.1.1](#) were analyzed and adressed.

99

100 Clearly it depends very much on your lifestyle, and how interested, willing, and capable you are to
101 recognize, deal with, (and in the future avoid?) situations that get you outside of your desired %TIR
102 on occasion.

103 So, this is also about what %TIR you are aiming at, and can accept, as it averages out for
104 the week, for instance.

105

106 Even if a principal capability for a fully automatic running FCL is given, this still
107 means that

108 • the user should be knowledgeable about what exactly is going on, and

- has a capability to „nudge“, or even to take over.

In [section 5.2](#) that immediately follows, we discuss the options the **FCL cockpit** on the AAPS home screen gives you, to occasionally tweak your FCL aggressiveness, should you see a need to do that.

5.2 FCL Cockpit: AAPS home screen as primary access point for modulating aggressiveness in daily use

Note: autoISF 3.0 is an early dev variant of AAPS, and as user you are participating in an on-going development. Of note, autoISF 3.0 is first launched without many of the described cockpit features that are presented below.

For the time being, multi-step work-arounds may become necessary. In many cases, going into AAPS Preferences and changing settings would be needed (...plus not forgetting to change these settings back, afterwards).

This is also an open invitation for you to contact us in case you could help program a module for one of the required user interface extras.

For future integration into AAPS Master, an eye should be kept also on the question which other modes (like FCL using Automations, or dynamicISF etc.) might benefit from some of the extra features.

The loop can run fully automatically without any user interaction (after the initial tuning phase, and related settings made in AAPS /preferences/SMB/autoISF).

But, just like in the airplane cockpit: Cruising in full auto mode should involve having an eye on the instruments, and on potential disturbances ahead in the environment.

E.g.: storm ahead => instruct your plane to climb to another flight height.

Analogy: exercise ahead => setting an exercise TT, or => pressing a button that activates a sequence of instructions (some of them probably hinging on conditions, like actual iob) how to manage through that exercise situation).

So, for the occasional „disturbance“ coming up, you should find an easy way to

- call up a pre-programmed routine for automatic management, with auto-adjusted aggressiveness, or:
- tweak a setting or two, to temporarily adjust the aggressiveness
- There may also arise a desire to just exit the FCL mode, and be your own captain for mastering a special situation.

All this is facilitated within seconds right from the AAPS home screen's **cockpit features** *(to the extent they are already incorporated):*

- The button that is integrated into the **violet FCL icon** serves as emergency off button, to quickly stop FCL, or to at least to immediately stop any more SMBs (...just for a couple of minutes, or for the remaining meal time: pick from the options offered with just one keystroke).

Via the violet FCL icon on your AAPS home screen, you also can access a temp. switch-off button for SMBs (see section that next follows below).

- The **three top fields** (%profile, exercise, TT) provide access to temp. tuning of core parameters, and/or to some pre-programmed routines.

Taken together with some **new indicator fields** about your loop state, this makes the AAPS home screen your **cockpit** for Full Closed Looping.

Let us look on each of these cockpit elements in some detail:

5.2.1 Violet FCL icon and underlying buttons

Novices to FCL, or really anyone running into a very special situation, may appreciate that the new closed loop icon on the AAPS home screen in pink (for FCL) has buttons to quickly shut off getting more SMBs (1st row), or to enter other loop modes (second row).


It functions very much as the other ones that you know from HCL already, and in fact you get offered some of the same options (for instance, to switch the (full) closed loop off for 15 minutes for going to take a shower)

Note that in FCL you leave all BG regulation, notably against meal spikes, to the loop. So, try not to disconnect in phases when your FCL must ramp up your iob.

The required insulin would still be supplied after you reconnect. However, without the user pre-bolussing, the delay would be more of an issue in FCL than it had been in HCL.

178

179 Just pressing on the FCL icon, a dialogue box comes up:



<-add-arrow, form-like-like-CL-Loop...¶

<-dotted+arrow, form-like-like-LGS¶



DEV: when in Hybrid-Closed-Loop = top headline + green circle), there is no Suspend SMB part but it starts right with „Closed-Loop“. In that screen. There, the 1st element must be option to pick¶

.....pink circle, FULL-CL-LOOP¶

where here the green HYBRID-CL-LOOP stands¶

For loopers who did not set up FCL, a feedback must come up ~ „FCL not installed“ if they press on that button¶

180

181

182 Pressing „**Suspend SMB**“ provides fast and easy „emergency braking“ regarding delivery of more
183 SMBs:

184 Select the one with the desired number of minutes: 5 or 12 for just blocking the potential next

185 SMB(s), and up to 3 hours to manage the entire rest of this meal with %TBR from then on.

186 Whenever, and whenever, your FCL is in „no SMBs allowed mode (e.g. automatically after

187 surpassing an iobTH also, or might be triggered by an odd TT), the FCL icon will turn into a dotted

188 one

189 Instead of remaining **duration to end time** it indicates in the middle „the condition“, „**iob**“ or „**TT**“

190 Add an indication **if** suspend SMB comes from an Automation, e.g. add an „**(A)**“, **underneath** the

191 #minutes, iob, or TT in the middle of the dotted violet field.

192 So, as in other (already in HCL existing) cases, those icons show in the middle the minutes left that

193 they will be running, or the condition which would have to go away for this temp. setting to stop.

194 It always auto-reverts into the FCL state and FCL icon, when time (or other condition) has elapsed.

195 Pressing „**HYBRID CL. LOOP**“ or other buttons from the 2nd row provides fast and easy
196 „emergency **exit**“ **into other modes**.

197 This enables beginners an easy „temp. escape“ into their well-known HCL (green) at any
198 point of time. bgAccel_ISF_weight is set to zero when going FCL->HCL. HCL can run with
199 autoISF (for instance dura_ISF) uninhibited otherwise. (check implications for HCL users of
200 autoISF ??).

201 Note: These options from row 2 have no time limit. Loop will **not** by itself go back to FCL. You see
202 the different loop icon as a reminder to manually revert, when ready.

203

204 5.2.2 Buttons „Insulin“, „Calculator“ etc at bottom of AAPS home screen

205

206 These buttons are **not useful any longer in FCL**, and automatically disappear whenever in FCL
207 mode (also in Suspend SMB state), and re-appear when leaving FCL. This applies also when an
208 Automation or technical system failure shut off FCL.

209 Users who, maybe in the beginning phase, feel better having those buttons, can override
210 the removal (of the insulin button, or any other) by going into /preferences/overview/buttons
211 and forcing them on. They only remain on until the next re-entry into FCL mode, when auto-
212 off happens again.

213 The reason why we do this: It really is important to let the loop loop, and not interfere more
214 than absolutely needed. Any bolus the user gives will sure distort the bg curve, on which
215 autoISF, especially when aggressively tuned for FCL, builds a lot of its decisions!

216

217 5.2.3. Three top fields (%profile, exercise, TT)

218

219 Depending on the variedness of lifestyle, the desired %TIR, and the initial tuning effort put in, the
220 user may want occasionally to „tweek“ the **aggressiveness of her/his FCL**.

221

222 The top 3 fields (grey in default mode, **yellow when temp. in mode with changed**

223 **agressiveness**) serve as quick and easy entry points to make temp. switches (as users will be
224 used to for %profile switches, or for setting an EatingSoonTT in HCL, .. which they still can do in
225 FCL ... but more:)

226 Expert FCL users might need this feature rarely, but probably at least to manage activity after
227 meals: Each require opposite aggressiveness, and the switch has to come in a certain point in
228 time that would be difficult to capture. (More see [section 6.4](#))

229

230

231 5.2.3.1 TT dialogue field

232

233 The TT field (top right of AAPS home screen) is the primary daily interface, and a dialogue field
234 opens when pressing on it:

HOME ACT INS COMBO

DIA7 90

Temporary target

Target ES

- 74 +

mg/dl

Target AC

- 139 +

mg/dl

Duration

- 165 +

min

iobTH%

- 70 +

-> 7.0 U

bGAccel%

- 100 +

-> 0.24

hiC

loC

niz

snk

ACTIVITY

wlk

ord

bik

mtb

h1

h2

Event time

10/31/23

21:30

CANCEL

OK

Initially, the form can be entirely empty re. TT inputs and just show the default iobTH and bGAccel parameters.

Two lines appear for target and get labeled ES, AC or HY depending on what was pushed at the bottom buttons. (HY => red frame; evtl. empty => grey)

TT entries can be made or overridden.

iobTH calculates from Target AC and other settings, shifts away from default set in /preferences (here 0.7) but could be overridden here. 0!..200 % is allowed

bGAccel ISF_weight can be modulated here, too. Note: it can change again if % profile is also changed. --0...200 % is allowed

CANCEL allows to start fresh (select one or two of the square buttons, ES, AC or HY)

OK needed to use the settings

Duration input is made in minutes. In the exceptional case that both, ES and AC targets are defined, the duration input is pr AC and framed blue. (This is because the preceding AC mode is automatically determined in length by the loop observing when iobTH is exceeded)

235

236

237 This looks complicated but only because it allows 4 different modes of use. Each user will primarily
238 use her/his preferred one.

ES

239 (1) Who is happy with the initially well tuned FCL and does not have huge variations in daily eating
240 and moving around, will **not use** the TT **at all**. FCL is possible without an intervention via the
241 TT button in your cockpit. Actually 4 of 8 modes (list at line 700) are not making use of TT.

242 (2) Super easy is also, to just input **any odd-TT** (odd-numbered temporary target) that will shut out
243 any SMBs for the set duration. *That can be a good idea when having a snack, for instance.*

244 Super quick access to stop SMBs is possible also via the loop icon ([section 5.2.1](#)).

245 Specifically, an **EatingSoon TT** can be activated here (*limited relevance see [section 2.5](#)*). It is
246 time-un-critical, can be manually set, or come up via an Automation.

247 *The cockpit enables you to set the iobTH differently (override) for the current meal.*

248 *(If the described feature is not yet included in your software version, change setting in*
249 */preferences or using an Automation)*

250 Temp. iobTH will always revert to default when the TT expires. If another TT immediately
251 follows, like in the example of the screen above, it will calculate, (then) show and use a new
252 temp. iobTH.

253 (3) The third way is to **use the input mask** (*if already ncluded in your software version* see picture
254 above) **to freely modulate the loop aggressiveness** for a declared number of minutes. Click
255 the bottom big square(s): Either HYPO, or ACTIVITY, or EATING SOON, or ACTIVITY and
256 EATING SOON (*example in the pictured screen above*). Make or override entries in the offered
257 fields. Press OK.

258 (4) The fourth way is to exclusively use one of the 4+4+2 little buttons seen in the bottom part of
259 the TT dialogue box (*if already included in your software version*). They provide a set of
260 settings (as will immediately show in all input fields above) that the user has set up in
261 Preferences/SMB/autoISF/FullLoop (refer to [section 6.3](#)), and can freely label there. *For*
262 *instance „hiC“ at high carb EatingSoon, „piz“ for Pizza/fatty meals, „grd“ for garden work,*
263 *„mtb“ for mountain biking ...*

264 **Capturing good settings for not-everyday situations in /preferences** (*if already included*)
265 **allows calling them up within 1 second**, from your cockpit on the AAPS home screen (...and
266 won't ruin the FCL experience at all , especially because in most cases it is not time-critical,
267 how long before the intended exercise the buttons are pressed).

268 The example picture given above is the most complicated (but also most useful) case, **when**
269 **exercise follows after a sizeable meal**. It is then that you need (a) aggressive FCL initial
270 performance at the meal, but, exactly when (!) a (for the intended sport already temp.lowered)
271 iobTH is exceeded, you need (b) to have SMBs automatically switched off and go into the

272 „milder“ mode, as defined for the exercise (with high instead of lowTT, that automatically
273 significantly reduces iobTH again, and insulin sensitivity(resistance) settings too).
274 Pressing exercise related buttons will automatically also light the **exercise button** on the main
275 screen yellow.
276 To summarize, the TT dialogue field offers easy but powerful ad-hoc modulation of loop
277 aggressiveness for FCL (*if already included*).
278
279 (tbd: Which of these UI elements will also be available in HCL autoISF ?)..

280

281 5.2.3.2 Exercise button (see more in [section 4](#).)

282

283 The exercise button automatically lights yellow when exercise related TTs are activated in the TT
284 dialogue box. 4 of 8 modes (list at [line 670](#)) are making use of the exercise button.

285 If pressing on the exercise button, a dialogue box appears

286 with info on exercise setting first (and opportunity to override), plus below the
287 activity monitor (experimental for auto-tracking of lighter movement during the day, and effects on
288 sensitivity that may have. See [section 4.5](#)).

289 So, first the exercis settings (as set under TT) are there to read. Example :

activity	TT	dura	%-profile	iobTH	bgAccel	%overall
mtb	171	180	70	3.0-U	0.24	70

Mode set to run for 134 more minutes

Mode starting after meal when iob > 4.4-U or n/a

290

291

292 The exercise (here mtb) is selected in the dialogue box of the neighboring TT field, and there auto-
293 filled with settings made in the set-up and tuning stage by the user under preferences (see
294 above,...). They are reported also under the exercise button here, and TT, duration, and % sens
295 (which also shows active on the %-profile field on the left side of the exercise button) can be temp.
296 changed there. iobTH, bgAccel_ISF and overall resulting sensitivity ratio is given in the other fields.
297 The **middle field** of the table, „% profile“ either picks up the % set under the %-profile button, or
298 an input can be made here, in the exercise button domain, which will:

- 299
- turn the neighboring %-profile button on yellow and show that inputted % on it, too
- 300
- be multiplied with the result from the exercise mode settings per se, and change the %
- 301 overall, accordingly.

302 So, if this middle field of above table (dialogue box of sports button) contains a figure other than
303 100, input field becomes yellow, and you are operating with a combination of traditional PLUS new
304 exercise mode (with all three top buttons of your FCL cockpit yellow). This maximally will soften
305 aggressiveness, for which you get an idea by the last calculated figure.

306

307 The mode is either running already (for another number of minutes, as probably also shown in the
308 yellow TT field anyways). Or it is scheduled to run, after insulination for a started meal reaches
309 iobTH (as in table). Or, no exercise is scheduled (both points red, no entries).

310
311 The lower part of the exercise dialogue box (not pictured above, but see in [section 6.5](#)) is
312 dedicated to the Activity Monitor

313 314 5.2.3.3 Profile button

315 The profile button can still be used to set a different profile, or profile%, for instance to adjust for
316 days with sickness (as you are used to from hybrid closed looping). 4 of 8 modes (list at [line 670](#))
317 are not making use of the profile button.

318 Any inputs made here will be used to modify profile_ISF on which all further changes are made on
319 (multiplied with).

320 The profile field remains grey if standard profile is applied.

321 It turns yellow, displaying a %figure relating to any altered loop overall aggressiveness:

- 322 • When no inputs (changes from 100% profile) are made here, but inputs in the TT field,
323 e.g. for exercise, automatically lead to different insulin sensitivity ratio (see [sections ...](#)
324 [and ...](#)), that ratio is shown here
- 325 • when% is changed by input in the profile button itself, it will be multiplied with with
326 profile_ISF and be used in place of profile_ISF by the algorithm.

327 However, for exercise (sports) you no longer must make an entry here, because
328 reasonable %reductions should be automatically provided, driven by your set TT (and half-basal
329 exercise target), see [section 6](#).

330 331 5.3 Recognizing your loop state in the AAPS home screen

332 333 5.3.1 Color scheme of top cockpit buttons tells kind of closed loop that is running

334 3 Buttons (%profile; exercise; TT) each in 2 states (yellow Y, or grey G) makes $2 \times 3 = 8$ possible
335 combinations:

336 GYY = dynamic exercise mode

337 YGY = not-dynamic „traditional“ exercise mode (if <100%) or hypo mode (if >100%)

338 GYG = basic closed loop with Activity Monitor running

339 GGG = basic closed loop (FCL or HCL) without any altered sensitivities etc

340 YGG = basic closed loop but with a „long wave“ sensitivity shift (e.g. sickness)

341 GGY = temp. target like e.g. EatingSDoonTT is set; or Hypo mode

342 YYG = closed loop with „long wave“ sensitivity adjustment and Activity Monitor running

343 YYY = dynamic exercise mode in time with additional „long-waved“ sensitivity shift

344

345 5.3.2 Information printed on the top buttons

346

347 The yellow TT field shows the currently valid TT (and further duration):

348 (profile) stands for the abbreviation you labeled your selected running profile



349

350 In the special case of settings for meal preceding sports, the field will look slightly differently:



351

...and ...

352 ... when iobTH is first time exceeded, this automatically switches to:



353

354

355 Likewise, if on the AAPS main screen just an **EatingSoonTT** is set (e.g.72), this is entered with the

356 desired duration. Afterwards, it automatically reverts to profile target and the display turns grey

357 again there with e.g. 90 on it (and no time limit).

358 Without sports context, the middle field remains grey.



359

360

361 Independently from setting a TT, the user can choose to set a **%profile in the left top field**, for an

362 independent number of minutes, e.g. 70% in this screen example: Also, or additionally, this will

363 influence the resulting ISF and sensitivity%

364



365

366

367 The % might change and turn yellow also in context of making TT inputs in the related dialogue
368 box (see chapter TT dialogue field, above). Still, the % (or the length of time the profile switch shall
369 be active) can be independently overridden in the top left field, if so desired.

370

371 If an **Automation** sets a %profile, and/or a TT (e.g. *automatic detection of meal start at condition*
372 *e.g. when delta >10*), this would automatically show in respective field(s) turning yellow and
373 showing the temp. setting. To show the set parameter comes from an Automation, „ **(A)** „ is added
374 in the end of button text.

375 Note that an Automation is *usually/often/always (?)* only permitted to temp. change default
376 profile settings, not other pre-existing temp. settings. This is for a good reason : Why should
377 a sometimes in the past thought-out Automation supersede your - just for the occasion
378 specified – temp.settings that you consciously activated for the day?

379 Advice: Try to stay away from Automations that also aim at temp. modifying
380 aggressiveness. For the reason just given in above note, they often will not kick in anyways.
381 Generally, it also is no good idea to double up sub-algorithms for tweaking loop behaviour.

382

383 Try to **keep things as simple and clear as possible**.

384

385 That said, a limited number of Automations can be of help in distinct scenarios (that differ in
386 purpose and in applicable time of day).

387 *A good one could be for night time, when your odd profile TT has SMBs shut off, but your*
388 *experience after pizza nights tells you that, under certain condition patterns (bg, iob), an*
389 *SMB or two should be „allowed in“. Another good example, if you go usually FCL without*
390 *any use of the TT button (which you could call a meal announcement of sorts), is to define*
391 *an Automation that, after detecting a meal start, automatically sets a low TT to get*
392 *maximally aggressive first SMBs.*

393

394 5.3.3 FCL related indicator fields in the AAPS home screen

395 In extra data fields of the AAPS main screen you can always see (not change) the key
396 „aggressiveness“ parameters your loop currently operates with (see also home screen
397 example below):


398 • how profile sensitivity (**ISF**) changes by the %profile input, by autoISF, and/or a set
399 exerciseTT.

400 • next to current available iob number is an indication of your **valid iobTH** (the iob above
401 which no more SMBs will be given)

402 • The AAPS home screen additionally shows, above the deltas, the current **acceleration**

403
404
405 5.3.4 Overall home screen:
406

Overall home screen:



<- Cockpit: yellow fields=>temp. modulated sens.
are additionally: acceleration-factor
<- violet <-> green circle for FCL <-> HCL
dotted if SMB temp off
<- % reduced insulin supply... because of resulting...
<- higher (weaker) ISF
<- buttons „bolus“ „carbs“ etc. eliminated
(auto-re-appearing when going violet -> green loop)

407
408
409