

## 4. Meals: Setting ISF\_weights in /Preferences V 2.1



**Please note that with autoISF 3.0 you are in an early-dev. environment,**

where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)

### Warning regarding importance of proper profile ISFs.

Starters on autoISF FCL who are coming from using HCL with **dynamicISF** must be aware of the following: It is absolutely essential to build your FCL on properly set **profile** ISFs (likely a circadian pattern over 24 hrs). It may not apply to you, but many dynamicISF users did never bother to determine their ISFs that would maximize their HCL performance, but employ dynamicISF so to speak for going „dynamically“ through a wide range of possible ISFs, until eventually hitting a sweet spot, and the whole thing works better than before with their profile ISF (often only one, e.g. coming from Autotune).

The following is important to understand, as it also leads straight into the core idea behind FCL with autoISF, too: It is a good idea to establish a well-running hybrid closed loop with set (non-dynamic) **ISF (as in profile each hour of the day)**. That ISF must be **aggressive enough** that it gets you down from a high around 200 mg/dl to target. That is roughly also the way you experimentally determined it (I hope).

Using that value also at lower bg, on the way up (after meal start), is very positive, as it is probably **stronger** than you would use, if you had just that (lower) bg to correct. autoISF will also do just that, but in a much more pronounced and elegant way.

On the way down from peak to glucose target, a somewhat too strong ISF will not hurt because much of the time your loop (well supplied with insulin before, „on the way up“) is zero temping or at least has only a small gap to correct from predicted bg to target bg.

You have no business to be much above 200 mg/dl where an even stronger ISF may or may not help. It sure does not help at an occlusion which is about the only reason to see super high values as an experienced looper.

Pegging ISF strength to bg level therefore does not make much sense for us. You will use the autoISF toolbox to get strongest ISF **at low** but beginning-to-rise bg,

Disclaimer: There are very much refined versions of dynamicISF that can have beneficial applications. But going to autoISF FCL, you absolutely must anchor on the proper profile\_ISF (which in times of illness etc. you can temp. change via profile switch, also when using autoISF in FCL).

## 39 Warning not to simply copy settings from others

40

41 When setting **your** parameters, **don't use any given numerical example**, but data from **your**  
42 *successful* Hybrid Closed Loop!

43

44 Most *examples given in this paper* are from an adult diabetic (Lyumjev, G6) whose insulin  
45 sensitivity can be characterized as follows: approximately 37 U TDD, thereof 13 U profile  
46 basal, at about 200g daily carbs from mainly lunch and dinner; no couch snacks or sweet  
47 drinks. The user also participates in multiple instances of daily moderate activity such as  
48 dog walking, biking and gardening. In Hybrid Closed Loop, a typical meal bolus was 8 U  
49 that was sometimes reduced such as when activity followed the meal.

50

51 After seeing some more inputs from a variety of users we might put together a profile helper  
52 for some rough orientation and plausibility cross-checking in [section 4.6](#)

53

## 54 Warning. Importance to starting from a well-performing Hybrid Closed Loop

55

56 A **satisfying performance in Hybrid Closed Loop** mode is a pre-requisite. Expect to reproduce  
57 about the same %TIR also in your FCL, but with less daily interaction, once established.

58 Note that this refers to prior use of „vanilla“ software, without fancy „dynamic add-ons“ (such as:  
59 Autotune determined factors, dynamicISF etc). that probably will not be compatible with autoISF  
60 use, and may have introduced bias into the profile settings you bring with you into FCL now.

61

62 To reach a satisfying performance you must start from a hybrid closed loop in which you did  
63 **master your meal management well** using the oref(1) algo SMB+UAM.

64 This is a pre-requisite **to be able to forget it** ... - because the initial tuning we now turn to  
65 demands that you analyze your prior best practice, in an attempt to find appropriate settings and  
66 „teach“ your FCL to come up with the necessary iob.

67 .

68 This is the main subject of this [section 4](#) (finding settings for automatic meal management) and  
69 [sections 5](#) and [6](#) (finding settings for highly automatic management also of other potential  
70 disturbances).

71

## 72 4.1 Getting started

73

74 Make sure you have studied the preceding [sections 1-3](#) on the general pre-requisites for FCL and  
75 on the workings of autoISF. Notably make sure you have set your default iobTH (refer to [section](#)  
76 [2.4](#) and if available [4.6](#))

77 In the early test phase, it is recommended to:

78 • run the system as dummy, not connected to your body (or, on own risk, connect only as  
79 long as you watch closely)

80 • in AAPS preferences, switch your autoISF FCL ( = **autoISF/"Enable adaptation of ISF to**  
81 **glucose behaviour"**) ON only during daytime hours of a meal, e.g. 11-18h, for fully  
82 automatic "full closed loop" management of lunches.

83 You can do this switching manually at 11 h and 18 h every day, or set up an  
84 Automation that does that (see [section 3.4](#) ).


85 • take typical but not extreme lunches. Omit sweet drinks, or drink only slowly

86 • do not use the Activity monitor (see [section 6.6](#)), unless it is already well calibrated.

87 In case you use an EatingSoonTT at meal start: Any active TT shuts activity monitor  
88 automatically off.

89 It is then essentially a matter of your UAM Full Closed Loop recognizing a meal start from the  
90 glucose trend, and ramping up iob.

91 When setting up your autoISF Full Closed Loop, you must set several ISF\_weight parameters in  
92 AAPS Preferences/OpenAPS SMB/autoISF settings. They relate to different stages of the typical  
93 glucose curve after starting a meal:



94 Make sure you studied the related flowcharts in [section 3](#). before you proceed.

96

97 **Warning:** Any bolus you „sneak in“ will severely distort the glucose curve. That can render  
98 your tuning of weights (see below) useless, and could make your loop act in unpredictable  
99 and potentially (see last bullet point below) dangerous ways:

- 100 • Your FCL ideally runs without an insulin button at the bottom of the AAPS home screen.
- 101 • Issuing a bolus should kick you out of the FCL mode, back into Hybrid Closed Loop. We are  
102 working on improving the User Interface (see e.g. [section 5.3.1](#)) that would facilitate and  
103 secure the bi-directional transitions.
- 104 • How proper autoISF settings would differ for your meal management, depending on no-  
105 bolus, bolus like in HCL, or very small pre-bolus (Meal Announcement), is not well  
106 investigated at this point.

107 Therefore it is best to remain consistent = use autoISF strictly for no-bolus FCL, and if you  
108 want to bolus for a meal, switch ISF adaptation to glucose behavior (autoISF) temporarily  
109 OFF.

- 110 • Maybe we are too cautious here, and in fact the autoISF adaptation to glucose behavior is  
111 tolerant enough of disturbances by **user boli**. Please report your findings in case you  
112 collect data of “mixed use” (FCL / Meal Announcement / HCL use with meal bolus). (A n=1  
113 finding, and guide how to evaluate, is reported here: [https://github.com/ga-](https://github.com/ga-zelle/autoISF/blob/A3.2.0.2_ai3.0/To%20prebolus%20or%20not%20to%20prebolus.pdf)  
114 [zelle/autoISF/blob/A3.2.0.2\\_ai3.0/To%20prebolus%20or%20not%20to%20prebolus.pdf](https://github.com/ga-zelle/autoISF/blob/A3.2.0.2_ai3.0/To%20prebolus%20or%20not%20to%20prebolus.pdf) ).

115 Once we have a body of data, including from those who moved from *HCL with*  
116 *autoISF* to FCL, we may need to re-define what the bi-directional transitions FCL <  
117 - > HCL in detail shall mean, and whether or not this has implications for needing  
118 different autoISF settings in /preferences for FCL and for HCL. .

119  
120 After doing the prep work as outlined in [section 2](#) you now get to calibrate your FCL to your **normal**  
121 **meal spectrum** by initially **setting and tuning the various \_ISF\_weights**, that dynamically  
122 change with bg curve characteristics as sketched in the chart on the previous page.

123  
124 Depending how satisfied you will be with the result, or which more extreme meals (smaller?  
125 faster/slower carbs? totally different fat/protein content?) you would like to manage, as well,  
126 you have a variety of options to deal with that, and this will be the topic in [sections 5](#) and,  
127 focused on exercise, in [section 6](#).

128  
129 In a nutshell, this will be about manual *or* (aided by Automations you would set up) semi-  
130 automatic (user triggered) *or* fully automatic **temporary modulation** of your FCL to deal

with *different* disturbances than presented by the meal spectrum you were calibrating for. This “nudging” will often involve:

- \* the **%profile button** (top left on your AAPS home screen). Note that the set % multiplies with both, the ISF resulting from autoISF and also with the default iobTH you have set, so both are nicely modulated in a linear way with the % temporarily chosen

- \* the **TT button** (top right on your AAPS home screen). Note that a lowered (relative to profile glucose target) TT signals lowered sensitivity (more insulin need), and an elevated TT (as often used with exercise) increases sensitivity and hence works in the direction of a lowered % profile to also reduce insulin given by the loop.

Moreover, the **exercise button** ((top center on your AAPS home screen) can be activated (turns yellow, then). This will **further boost** how your set TT elevates the resulting ISF, and sharply lowers iobTH, as often desired for sports. See [section 6.1](#)).

Taken together with a couple of more features ([section 5.2](#) and [6.3](#)), these functions make the AAPS home screen your **FCL cockpit**.

So, yes, FCL is about fully automatic cruising. However, you have a cockpit to check how everything is running, to sometimes “nudge” – and in exceptional situations also to temporarily intervene, or even take over.

Also, like a pilot, you need to learn a bit, so everything will fly well.

But: You should do some fair weather stuff first \*), which brings us back to our intended next step:

Researching your standard meal patterns, and finding settings for the various ..\_ISF\_weights.

\*) If, to keep the motivation up for your project, you are itching to see what fancy stuff can all be done, you might peek for instance into [case study 6.2](#). And if that looks like way too much, decide to be just a fair weather flyer for now - or, no offense taken, give up now before spending too much effort. [Section 13](#) is about (maybe) “easier” alternatives.

## 4.2 bgAccel\_ISF\_weight

When looping without carb inputs and without giving a bolus ourselves, the first crucial setting is to set the **bgAccel\_ISF\_weight** so that large SMBs are requested immediately when the loop detects an acceleration in your blood glucose (bg) that is starting to rise.

Ideally within about 20 minutes after acceleration detection, which would be the first up to 4 SMBs, as much iob should automatically be supplied as we would have given with our bolus in hybrid closed loop.

.

If you double the bgAccel\_ISF\_weight, the ISF strengthening is doubled.

169 Rule of thumb: Two of the first three SMBs each should be about  $\frac{1}{4}$  (max  $\frac{1}{3}$ ) the size of a  
170 previous meal bolus in your HCL „career“.

171       Going over  $\frac{1}{3}$  can be problematic if your diet contains occasional low carb (or only  
172 snacking), and generally of course if your CGM quality is sometimes unreliable, and might  
173 produce an artefact that could be mistaken for a meal start. Be vigilant about this topic!

174

175 For hands-off FCL, your settings have to fit the whole **range of your meals**. In extreme cases you  
176 will have to balance too high running iob with additional carbs (a late additional snack against  
177 going too low), and in the opposite case, you will have to reckon with temporarily exceeding the  
178 glucose target range and losses of the achieved %TIR for this day.

179

180 **If your meals vary very strongly**, there are avenues to ease your initial tuning job, or to optimize  
181 overall resulting loop performance:

182       • Automations allow you to differentiate. For instance it is possible to apply different  
183 iobTH\_percent and/or different bgAccel\_ISF\_weights for meals in different **time windows**  
184 or geo locations (details see [sections 3.4](#) and [5.1](#))

185       In case you use autoISF 3.0 on the iAPS platform for i-phones, you need to use a third  
186 party automation software (! [call for a case study 4.X](#) )

187       • You can pre-program **custom buttons for special** meal (or snack) **types**, with different  
188 underlying FCL settings (see “cockpit”, [section 5.2.1](#))

189       • You can **modulate FCL aggressiveness manually** making use of temporary switches  
190 of %profile and/or set glucose target ([section 5.2.3](#))

191       • [In an update, autoISF 3.x might provide the option to pre-program settings for 4 different](#)  
192 [meal type clusters, accessible from the TT button \(presented in \[section 5.2.3.1 \\(4\\)\]\(#\) and \[6.3\]\(#\)\).](#)

193

194 In search of appropriate settings, you must keep (real-time) track of the **SMB tab** when tuning. This  
195 can be impractical. You probably will end up making a lot of screenshots (quickly in the crucial  
196 minutes where the SMBs were given, or when you thought they should be given), for later analysis.

197

198 The superior method is to just copy **logfiles** about once a day from your phone/internal  
199 memory/AAPS/logs (all zip files there), and analyze them at your convenience later, using the  
200 **emulator** (see [section 10](#)). Some emulator-based analysis is also possible within AAPS on your  
201 phone ([section-11](#)).

202

203 Already when tuning the bgAccel\_ISF\_weight it can become evident that safety restrictions (as  
204 discussed in [section 2](#)) must be widened further:



- Especially if your *profile basal* rate is very small, the **smb\_delivery\_ratio** and/or the **smb\_max\_range\_extention** "must" often be increased further.
- Furthermore, the **smb\_delivery\_ratio** provides more leeway to increase the aggressiveness (e.g. 0.6 -> 0.72 results in another +20%).

In the end you should not set the limits too tight, so "nudging" aggressiveness by another 10 or 20% from your cockpit later will not bounce into your set limits.

In any case, it is worth the effort to tune the **bgAccel\_ISF\_weight** in such a way that high glucose increases are already nipped in the bud, so to speak.

This also facilitates the tuning task for the subsequent phases of the meal, because there is then largely zero-tempering, as well known from HCL-times after YOUR administered bolus. Also, the lower and shorter lasting the glucose peak, the lesser the hypo danger from the activity tail of SMBs given when glucose was „stuck“ high.

Default bgAccel\_ISF\_weight is set to zero in autoISF. **To start**, I would try 0.05 or **max 0.1**, and keep trying in max 0.05 steps. Soon move to 0.02 steps (which still means 10-20% change). From my (very limited) overview, many use around 0.2, but possibly higher if their hourly basal rate is 0.1U or lower. ([Consult section 4.6 when available](#)). Do not be tempted to rush this setting by using large jumps in adjustments.

Ideally, one should set the bgAccel\_ISF\_weight so that for meals that are in the **lower** range of the "fast **carb load**" of your cluster, the necessary insulin supply is already approximately provided with 3 SMBs. The glucose curve, at such meals, begins to flatten early in this SMB phase, so a deceleration follows very soon (-> section 4.4).

### 4.3 pp\_ISF\_weight

With **higher carb load** meals, or meals that come with a sweet drink, the acceleration phase will last longer, and BG will rise further, which will require a higher insulin supply.

Between acceleration and deceleration there is a more or less linear further increase of insulin need in these cases.

Our autoISF should now "fight" this with the help of the post-prandial ISF, set via **pp\_ISF\_weight**, after we have set a halfway suitable bgAccel\_ISF\_weight.

Select **pp\_ISF\_postprandial all day = ON**

In full closed loop mode, this parameter is preferred over deltaISF (.. and highly beneficial also for managing meals with gastroparesis)).

Tune your **pp\_ISF\_weight** after you have set a halfway suitable bgAccel\_ISF\_weight. You should check meals in the upper spectrum of your g carb, and carefully start tuning with a weight of 0.01.

Normally the SMBs triggered by bgAccel\_ISF\_weight and pp\_ISF\_weight should be sufficient to reach and slightly exceed the **iobTH** (see [section 2.4](#)) so all *the other* autoISF parameters are relatively unimportant for now.

A reason why this can work at all, also for quite a variety of meals, lies in the fact that there is an hourly carb absorption limit of about 30g/h (reference: Dana Lewis: <https://github.com/danamlewis/artificialpancreasbook/blob/master/8.-tips-and-tricks-for-real-life-with-an-aps.md#heres-the-detailed-explanation-of-what-we-learned> ). So while meals might wildly vary in composition and size, what is digested and needs insulin in the first <90 minutes (when FCL tries to catch up with insulin need and differs strongly from HCL, and bgAccel\_ISF and pp\_ISF play the leading role) will be relatively close (...for meals with similar *initial* glucose acceleration and rises, anyways, and *that* is information the loop does have).

Depending on the type of meal and "aggressiveness" of your bgAccel\_ISF\_weight and pp\_ISF\_weight tuning, the iob will already be so high that in the phase of decelerated glucose rise towards the peak (the "last part of the rise") that no insulinReq is seen by the loop.

Therefore the **bgBrake\_ISF\_weight** is often unimportant (-> [section 4.4](#))

**Warning:** Occasionally consult the SMB tab to see how your settings really work.

A setting that is actually set too aggressive might be masked. **Tuning only works if** the effects of the settings being tuned are **not** unintentionally **limited by other** (e.g. „safety“) **settings**.

Also, always look at two or three different meals before deciding whether a tuning "fits" („good enough“ for each of them)

[Case Study 4.1](#) (Pizza Meal) contains, towards the end, an example how you can go about tuning the \_weights for various \_ISF factors of autoISF.



#### 4.4 bgBrake\_ISF\_weight

At a **low carb** meal, or an attempt at doing a **weight reduction diet**, the glucose goes up only sluggishly and iobTH should not be reached at all.

Acceleration and the phase of strong glucose rise are quickly over in these cases, and there is mainly a decelerating bulge of insulin action that projects over the next few hours.

Now the importance of the **bgBrake\_ISF\_weight** comes in. In full closed loop, the bgBrake\_ISF\_weight is often only about half as large as the bgAccel\_ISF\_weight (but that would also depend on your personal diet pattern and eating/digestion speed). Also here, one should approach the tuning gradually, increasing the weight from small values.

[Case Study 4.2](#) shows a user example of a low carb meal managed in FCL by autoISF.

What is very helpful for us in any case is that the loop calculates the situation every 5 minutes, and corrects it.

However, if there is too much insulin in the system, the loop can only correct to a very limited extent, namely only to the extent that it can set basal to zero.

Therefore, the core problem is that the Full Closed Loop must build up iob very quickly, but not too much, in the initial phase of a meal.

But high BG values (out of range, >180 mg/dl) can not always be avoided..

**Note regarding acceleration happening in late part of dropping glucose:**

In version 2.2.8.2 there was a potential deficiency in situations where glucose was falling and the glucose acceleration was already positive. That meant a minimum glucose level can be extrapolated. If that happens to be less than target and expected in less than 15 minutes then there should be no strengthening of ISF as it would lower glucose even more. Therefore bgBrake\_ISF\_weight is used now instead of bgAccel\_ISF\_weight. But those situations were rare and less critical than might be expected at first sight. The reason is that in most cases the predictions ended up even below their threshold meaning SMB were disabled.

#### 4.5 High Glucose Values and dura\_ISF\_weight

With **large or high fat/protein** meals, a 2nd hill of glucose will form, or a long high plateau.

For such situations there is in autoISF the modulation of ISF depending on bg level or duration of plateau formation.

314 High bg values and a plateaus in bg values are tuned using the **dura\_ISF\_weight** and associated  
315 parameters. This feature is also very useful in Hybrid Closed Loop. It elegantly manages, fully  
316 automatically, temporary fatty acid resistance. Please refer to other papers for details (for instance,  
317 section „Late stage of meals“ of:

318 [https://www.facebook.com/download/649096606100188/MealMgt.Basics\\_09Dec21.pdf](https://www.facebook.com/download/649096606100188/MealMgt.Basics_09Dec21.pdf) ).

319

320 Since in Full Closed Loop we "turn up" our loop to give the maximum SMB size we can at the  
321 beginning of a rise, it is crucial to **resist the temptation to continue with a particularly strong**  
322 **ISF in the meal phase with the highest glucose values** .

323 This is a reason why in Full Loop we do not make much use of the **bg\_ISF** component of autoISF.  
324 Near glucose peak, zero-temping usually prevails anyway, so the settings we try might often not be  
325 used really by the loop.

326 This also means that too aggressive settings might not come into play most of the time. However,  
327 some other time they might come into play, and produce a hypo 1-2 hours later.

328

329 Therefore, **carefully study the SMB tab to see what the selected weights would do, if there**  
330 **was no zero-temping** at the time. Also, try a completely different meal to see how your settings  
331 work there.

332

333 The UAM Full Closed Loop doesn't get any information from you as to how many grams of  
334 carbs will be absorbed late. Not knowing when your steady-state max carb absorption  
335 phase (the earlier mentioned 30g/h), and even sometimes a brief episode of insulin  
336 resistance to fats, might end, the FCL will struggle to provide desired amounts of insulin,  
337 facing potential hypo danger later because of the DIA of the insulin in use.

338

339 Actually, the UAM Full Closed Loop is *not completely clueless* regarding how carb  
340 absorption will go on. It will work with a prediction of further carb absorption building on the  
341 **carb deviation** (=hypothesis of how much got absorbed in the past 5 minute segments),  
342 and phase out more carb decay in the course of the next 1 to max 3 hours. For more detail  
343 see

344 <https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand-determine-basal.html#understanding-the-basic-logic-written-version> (or study your  
345 SMB tab info).

346

347  
348 This UAM prediction about further carb absorption can be worse, but can also be better  
349 than a prediction based on the user's „e-Carb“ input in Hybrid Closed Loop.

350

351 In any case, and even when having perfect knowledge about how exactly the carbs fade  
352 out in the next hours, there would still be a principal problem for the loop: Heavy insulin  
353 „fire“ against highs will not work immediately (depending on the insulin's time-to-peak), and  
354 notably it comes with a significant hypo danger (from the „tail“ of insulin activity.)  
355 A big bolus, or also a series of boli, will rarely work exactly for several hours matching the  
356 absorption of carbs (from what, how much and and how fast the user ate).

357  
358 Once your BG sits high, neither you, nor a hybrid closed loop with all the carb info, nor your FCL  
359 can work wonders. Resist the temptation to elevate the **dura\_ISF\_weight** very high.

360  
361 The author is sceptical about using the **bg\_ISF** (at least be careful, use small weight, or shut-off).  
362 Highs will take time to resolve. Interestingly, an after-dinner walk can work wonders sometimes.

363  
364 As often, the best solution is to be pro-active:

365  
366 **The earlier large SMBs come** (driven by bgAccel\_ISF and pp\_ISF), the less high the  
367 overall increase in BG will be, and (provided you set a proper iobTH) **the lesser the risk**  
368 **will be for a hypo after the meal.**

369 Therefore, put most of your FCL tuning effort into determining suitable weights for  
370 bgAccel\_ and for pp\_ISF, and a suitable iobTH.

371  
372 Your FCL cockpit will give you access to modulate 2 of these 3 essential parameters (see  
373 [section 5.2.](#)), providing you an opportunity for more research on the fly, so to speak.

374  
375 The experience of the author is that it is possible to tune the above mentioned weights for very  
376 different meals in such a way that the glucose almost always remains acceptably in range.

377  
378 However, if you need **differentiated settings** for different meals or meal time clusters, you can

- 379 • manually modulate FCL aggressiveness via setting temp. %profile and/or TT (see [section](#)  
380 [5.2.3](#))
- 381 • or install and activate a user defined extra button in your cockpit for it (see [section 5.2.1](#))
- 382 • or pre-program 4 different clusters in /preferences, and call them up within a second from  
383 the TT button in your AAPS home screen (*only after implementation of an improved cockpit,*  
384 see [section 6.4.3](#))

385  
386

## 387 4.6 Profile helper

388

389 xls based tool is still under development / needs more user data / chapter will follow later