

7. Kids: Mastering additional challenges... (NN) V2.2

Please note that with autoISF 3.0 you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



- - - - - call for contributions, chapter and case studies - - - - -
-

Points from section 1.6 that need to be detailed here :

To establish and maintain a FCL for kids brings about some extra challenges if:

- Lyumjev is not available or well tolerated
- Hourly basal rate is very low, providing a poor basis for big SMBs => *Measures to get sufficient initial SMB sizes?*
- Diet is rich in sweet components. With the typical low blood volume of a small body, strong tendency towards very high bg spikes! => *Provisions for erratic hi carb snacking/ the occasional sweet drink? (see [case study 5.2](#) - is this applicable ?)*
- Going through marked changes of insulin sensitivity or of circadian pattern makes it difficult to keep the FCL appropriately tuned.

This problem is about the same in Hybrid Closed Looping. However, now you might expect miracles from the FCL. This is not going to happen. You still should try to set appropriate temp. changed profiles, that serve also as a basis for your autoISF FCL.

- Discipline is poor regarding keeping Bluetooth connectivity and infusion sites perfectly running (*=> alarms on parents phone, see section 1.3 Automation example ?*)
- Between kid and supervising parent it must be guaranteed, especially in the initial weeks, that an eye is kept on whether the FCL is working about as to be expected.

FCL => implications for remote parent monitoring ? (may have to come at the end, when all potential problems were discussed)

- Extra caution needed re. elevating the **SMB delivery ratio** ?
 - The fixed 0.5 value in AAPS was installed also with a consideration on user/follower (parent) set up and limiting potential problems from a bolus being initiated from both phones in parallel
 - As in FCL neither phone should initiate a manual bolus, this precaution seems not really needed (**dev – check**)

36
37 The **DIY Cockpit** (see quoted below, from [section 5.2.2.3](#)) allows to have buttons dedicated to
38 kindergarten hours etc that can in a time-uncritical manner activate special routines – **should be**
39 **very helpful to custom program Automations for kids!**

40 -----quoted section 5.2.2.3 -----

41 5.2.2.3 Triggered Automations: Grey extra DIY cockpit buttons for pre-programmed “responses”

42

43 Recognizing conditions for fully automatic handling by the loop may not be not possible, or come
44 too late for the loop to act on. Examples would be

- 45 • *exercise*: Minimum an hour before starting “the loop should know” to be able to lower iob
46 and elevate bg by the time exercise starts.
- 47 • *snacks*: High carb snacks, sweets, consuming ice cream or having a sweet drink comes
48 with the problem of even steeper glucose rises, but overall a lesser insulin need, compared
49 to major meals (for which we tuned our FCL according to [section 4](#)).

50

51 This not necessarily implies that snacks need different settings than a meal. After all, autoISF
52 was designed to act to all available data, especially to where the developing glucose curve is
53 headed. So, depending on your effort to set parameters for a broad variety of meals (notably:
54 how well you avoid to invariably bounce fast against your iobTH), you might be able to accom-
55 modate low carb, snack, and major meals with *one* set of settings.

56

57 In FCL autoISF, this is a bit more difficult than in HCL autoISF applications, because FCL
58 involves revving up iob supply (largely via big bgAccel_ISF-weights) often too much to be
59 balanced by just a snack getting absorbed.

60

61 For that reason, or just for increased comfort and safety, you might want to differentiate, and make
62 use of what follows for the *sweet snack* example.

63

64 **Tuning aggressiveness**

65 Key is that a sweet snack likely benefits from even more aggressive initial FCL
66 performance than the meals in your normal spectrum of diets require.

67 Therefore, you could set

- 68 • a higher **temp. profile%** and/or
- 69 • a temp.elevated **bgAccel_ISF-weight** (see screenshot of my Automation).
- 70 • a **low temp. target** (76 for instance; this additionally helps maximize the first SMBs
71 that will automatically be triggered at detection of acceleration)..

72

73 When first defining and testing this Automation, also check:

74 • that the safety limits as discussed in [section 2](#) will not block the intended elevated
75 aggressiveness

76 • SMBs will not get outrageously big and iobTH sometimes exceeded by too much

77 Note that “the last SMB” is allowed to overshoot the effective iobTH by 30%


78

79 Limiting iob

80 For “just a snack”, total insulin need will be lower than for a meal.

81 If you would just have your sweet drink, and your meal-oriented FCL would “attack”,
82 iob likely would become too high, and a glucose rollercoaster would start, with you
83 needing to consume more =>

84 If you just have a snack, or drink a glass of juice, you can lower the **iobTH_percent**
85 accordingly.




„User action“ is always ticked-on

This will, in the defined time space ..

..offer the “DIY cockpit” button..

..which I must press any time (~90...30 minutes) before my snack.

Note: Make sure that iobTH and bgAccel_ISF_weight revert to default afterwards



86

87

88 So, this can be a little **extra “project” when setting up your FCL.**

89 You need to research your snack habits (if any), and over time find out which settings in
90 the snack-related Automations work well.

91

92 **In everyday life** you **then just** must **press the related button in your cockpit** (which is
93 not time critical at all, except it should be clicked *latest* a couple of minutes after you took
94 the drink or snack).


95
96 If you consume more, and also eat something with your sweet drink, this will more
97 resemble a full meal... however, with unusual amounts of fast carbs.

98
99 Note: Pressing your snack button *a second time* would **not** help because the
100 lowered iobTH does not allow iob going high enough. So you are better off just
101 letting your normal FCL meal routine run, after your snack mode expired.

102
103 Other options when snacks keep extending would require a manual modulation
104 regarding %profile and/or bgAccel_ISF, but keeping the full default set
105 iobTH_percent, or even elevating it (refer to [section 5.2.3](#)). If that happens often,
106 define for yourself an extra User action Automation for a bigger snack (= another
107 grey DIY cockpit button).

108
109 **Caution:** If (as in autoISF 3.0) setting a different iobTH or bgAccel_ISF_weight can not be done
110 temporarily (i.e. with a duration attached) you **must** define a suitable **additional Automation, that**
111 **must be active in tandem, that restores the default**
112 **set iobTH or bgAccel-ISF_weight again.** Else, once
113 your Automation set in, it will *forever* shift this
114 important parameter setting!

115 If for instance you have several Automations that, in
116 combination with a set elevated TT also set a lower
117 iobTH: Don't be fooled, the duration only applies to
118 the TT. You need an extra Automation for all of them.
119 I picked out the highest of the altered iobTH values
120 that these Automations can set (45_percent), and
121 then I can automatically restore my default desired
122 60% via this one Automation (see screenshot - - >)



123
124 **Installing the DIY cockpit button**

125 In the related Automation, just keep the “User action” box clicked at all times, and define in
126 the Conditions when you want to see that button available for cockpit use (see screenshot
127 above) => you will see that button offered.

128

129 Besides snacks, also any **other recurring special situations can be addressed via a**
130 **DIY cockpit button, and receive different aggressiveness up to a suitable iobTH**
131 **level.**

132

133 Over time you can have a big number of User action Automations, and keep them
134 “shelved” rather invisibly (clicked in-active, top left box) in your long list of potential
135 Automations. Even when active, they only show in your cockpit (bottom grey field of your
136 AAPS home screen) in the time slot you assigned as potentially relevant.

137

138 -----end of quote from chapter 5.2.2.3 -----

139

140 *End with something like:*

141 There are a couple of parents and kids pioneering this area, see also [Case studies 7.x](#).

142 We highlighted areas that would require some minimum compliance. In the end it comes down to
143 compare the achieved ease in daily use and achieved %TIR to how it was when hybrid closed
144 looping.

145

146

147 **If users pre-bolus for kids**always, or sometimes? How(much)? =>

148 Make sure to give **details** ...

149 and possibly a **warning**, regarding implications for settings (I could imagine bgAccel_ISF-
150 weight tuning should be factor 5 if not 10 milder in that case = ... and not be good if sometimes you
151 pre-bolus a bit, and sometimes not ????)

152

153 *Re. pre-requisites, **CGM, for kids**: what works well ??*

154 *If not G6 is used:* (re CGMs, observe latest developments re 1 minute vs 5 minute readings, see
155 github/ga-zelle pls)

156

157 *Pls include implications of (from wiki aaps):*

158 **Dexcom G6 or Dexcom ONE with xDrip+ Companion Mode**[S](#)

159 The noise-level data is not shared with AAPS using this method. Therefore ‘Enable SMB
160 always’ and ‘Enable SMB after carbs’ are disabled.

161 **xDrip+ with FreeStyle Libre**[S](#)

None of the FreeStyle Libre systems (FSL1, FSL2, or FSL3) broadcast any information about the level of noise detected in the readings, and therefore 'Enable SMB always' 'Enable SMB after carbs' are disabled for all setups using the FreeStyle Libre.

In addition, many people have reported the FreeStyle Libre often produces noisy data. In xDrip+ there are a few options to help with this:

Smooth Sensor Noise. In xDrip+ Settings > xDrip+ Display Settings ensure that Smooth Sensor Noise is turned on. This attempts to apply smoothing to noisy data.

Smooth Sensor Noise (Ultrasensitive). If you are still seeing noisy data in xDrip+ you can apply more aggressive smoothing using the Smooth Sensor Noise (Ultrasensitive) setting. This will attempt to apply smoothing even on very low levels of detected noise. To do this, first enable [engineering mode](#) in xDrip+. Then navigate to Settings > xDrip+ Display Settings and turn on Smooth Sensor Noise (Ultrasensitive).