

3 1.1 Well-tuned hybrid closed loop

5 It is advisable to first establish a well-tuned hybrid closed loop before considering the transition to
6 FCL. There are two important reasons for this:

- 7 • The UAM full closed loop requires a highly personalized (individual) tuning of settings, so
8 the loop will give insulin mimicking YOUR successful hybrid closed loop mode.
- 9 • The UAM full closed loop comes with new parameters to be set and tuned. It would be
10 problematic to set and tune several new parameters before the basics were tuned "right".
11 Errors could easily be balanced with counter-errors. This can work in single scenarios, but
12 would create a highly unstable system, hard to re-calibrate better later.

14 1.2 Fast insulin (Lyumjev, Fiasp)

16 If the user does not bolus for meals, clearly a very fast insulin is needed so, upon realization of a
17 starting meal-related glucose rise, the loop has any chance to eventually keep glucose in range (by
18 common definition, under 180 mg/dl (10 mmol/l))

19 A modelling study (details see

20 <https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#fast-insulin-lyumjev-fiasp>)

21 can show in quantitative terms that **faster insulins**

- 22 • will result in significantly **lower** glucose **peaks** than slower insulins
- 23 • **tolerate** a couple of minutes **delayed** first meal bolus while not incurring unacceptable
24 height of peaks
- 25 • minimize the effect on glucose peak from **different** carb load (**meal sizes**).

26 In conclusion, do not attempt FCL with other insulin than Lyumjev or Fiasp, unless, maybe, if you
27 are on a very moderate to low carb diet. (See also [case study 1.2](#)).

29 1.3 Good Tolerance (of Lyumjev or of Fiasp) because occlusions threaten the
30 function of the full closed loop.

32 It is very important to have an eye on the time a **cannula (or pod)** is in use (many find **48 hrs** to be
33 the **limit**), and whether hard-to-explain glucose rises happen at ever increasing „fake“ iob (even
34 before a 48 hr routine replacement). (See [case study 1.1](#): You easily lose 25% TIR that day)

35 It is absolutely contra-indicated to attempt FCL coming from leaking pods and associated erratic
36 sensitivity swings that may or may not have been somewhat controlled and tolerable by
37 dynamicISF or other measures when you were Hybrid Closed Looping,

1.4 Excellent CGM

You do not give a meal-size-related bolus any longer. That leaves all insulination jobs to the algo! As glucose values are the very basis for this, please **inform yourself well about** how **your CGM** 1) principally performs 2) whether and how this may depend on data flow and intermediate apps you use 3) specifically, how and where any smoothing is done, and what this might imply for the ISF boosting method you will be using. See for instance here:

<https://androidaps.readthedocs.io/en/latest/Usage/Smoothing-Blood-Glucose-Data-in-xDrip.html>

In AAPS Preferences/OpenAPS SMB/autoISF/Glucose source data for parabola fit, you must select between 4 options relating to your CGM (1 or 5 minute values, raw or smoothened).

Around meals, a stable Bluetooth connectivity is absolutely essential, too, so CGM, loop, and pump can do their job without losing more valuable time. (See [case study 1.4](#)).

Then, but even more importantly in *all other* day and night *times*, the CGM should not produce any artefacts (jumpy values; see [case study 1.3](#)) that the loop could **misinterpret** as sign of a starting meal. Note that also calibrations could produce jumps.

But, in any case, a CGM with more scatter will make the loop lose more time, and lead to higher peaks and lower %TIR.

The best way currently is to use Dexcom G5 or **G6**, and to ensure via **overlapping** right and left arm sensor and transmitter utilization, that always good quality values can be used by the loop. Other ways are possible, but come with a lot of monitoring effort (via watch) and occasional time-outs for the loop.

Also FS **Libre 3** is useable. Observe info in the general section about autoISF (<https://github.com/ga-zelle/autoISF>) regarding how it should be implemented:

Hint for users with 1-minute CGM (Libre 3) (*status of 05 Dec. 2023*)

In AAPS 3.2.0.2 the glucose history only uses data at 5 minute intervals and ignores the ones in between. That has awkward consequences for the graph display and for the exponential smoothing method: Each minute a new reading is received the 5-minute pattern moves forward by 1 minute and a completely fresh and new subset of glucose values and times is used.

Alternatives of determining the glucose acceleration based on the 1 minute data are still under evaluation. Once a promising method can be found and validated then an interim release will be provided.

As a libre user you have 3 alternatives:

- AAPS gets values from xDrip+ where you smooth and reduce the 1-minute data to 5-minute data before sending to AAPS
- you should stick to the 5 minute mode in your Libre setup if possible
- AAPS gets values from Juggluco every minute but AAPS uses only the 5-minute subset for looping and for parabola fitting

1.5 Meal-related limitations

Setting up a full closed loop is relatively easy for people whose diet does not consist **mainly** of components with rapid high effect on blood glucose (more see <https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#meal-related-limitations>)

Meals do not have to be low on carb (provided you use a fast insulin for your FCL)

Fat or protein rich diets, or slow digestion/gastroparesis, make things easier rather than harder for the full closed loop because late carbs nicely cover for inevitable “tails” of late action from SMBs needed around peak time.

Erratic consumption of snacks with fast resorbing carbs can be a problem.

In autoISF you can reduce this problem to some extent via two or three keystrokes on your AAPS main screen. While certainly being a deviation from the FCL idea(I), this would be one of the exceptional situations where you better do a quick “nudging” step from your “FCL cockpit”. Details see in [section 5.2](#) /TT dialogue field / (4).

1.6 Lifestyle-related limitations

Technically stable system

Full closed looping requires a 24/7 technically stable system, especially regarding reliable **CGM** signals, but also **Bluetooth stability** with the **pump** (see [case study 1.4](#)), and avoiding (or at least early recognition of) occlusion. This could require attention to details like keeping all components

116 well charged and in close proximity; making cannula (or pod) changes always early enough to
117 lower the risk of occlusion (see [case study 1.1](#)); having always potentially needed parts with you.
118 **Depending on your system, your experience with it, but also on your acceptance and general**
119 **lifestyle, these aspects may or may not limit you.**

120

121 Preparing for exercise

122

123 To prepare for exercise (sports, heavy work), the normal protocol with a pump or hybrid closed loop
124 is to take actions that reduce insulin on board prior to exercise

125 With your full closed loop, the algorithm is tuned to detect meals and to give you insulin to counter
126 glucose rises automatically. Setting a high temp. target and lower %profile right away (effective al-
127 ready around meal start) could be a problem.

128 Unusual activity levels therefore likely require **disciplined preparation** (especially **if you want to**
129 **keep the need to snack during sports low**)

130 In autoISF you can reduce this problem to some extent via two or three keystrokes on your
131 AAPS main screen. While certainly being a deviation from the FCL idea(l), this would be
132 one of the exceptional situations where you better “flick a lever” from your “FCL cockpit” to
133 have temporarily adjusted settings for the planned exercise. Details see in [section 6.2-6.3](#)

134

135 Extra hurdles to establish FCL for kids

136

137 To establish and maintain a FCL for kids brings about some extra challenges if:

- 138 • Lyumjev is not available or well tolerated
- 139 • Hourly basal rate is very low, providing a poor basis for big SMBs
- 140 • Diet is rich in sweet components. With the typical low blood volume of a small body, strong
141 tendency towards very high bg spikes!
- 142 • Going through marked changes of insulin sensitivity or of circadian pattern makes it diffi-
143 cult to keep the FCL appropriately tuned.

144 This problem is about the same in Hybrid Closed Looping. However, now you might
145 expect miracles from the FCL. This is not going to happen. You still should try to set
146 appropriate temp. changed profiles, that serve also as a basis for your autoISF FCL.

- 147 • Discipline is poor regarding keeping Bluetooth connectivity and infusion sites perfectly run-
148 ning
- 149 • Between kid and supervising parent it must be guaranteed, especially in the initial weeks, that
150 an eye is kept on whether the FCL is working about as to be expected.

151 More see [section 7](#).

152

153

154 155 1.7 Time required for setting-up 156

157 Lastly, before enjoying a functioning full closed loop you need to have a period of a some weeks
158 with some free time and „free head“ for set-up –. Can you get, in the time you are willing to invest,
159 to a result that you consider good-enough is really the question. Depending on your „habits“, and
160 which – if any - compromises (like doing cannula/pod changes more often, never starting meals
161 when bg sits high ...) are you willing to make (and everyday able to stick to), for the ease of not
162 having to deal with assessing meals and bolussing for them?

163
164 While setting up your personal FCL using autoISF is a substantial project, there is no need
165 to implement it fully in one step. There is nothing wrong to go in your well running Hybrid
166 Closed Loop mostly, while switching to FCL only for dinners, for instance, or only for
167 weekend lunches, as a start. Once you found feasible settings, you can expand to other
168 meal times, and lastly towards figuring out your best strategies for challenges outside of
169 meal windows, as we shall discuss in [sections 5. and 6.](#)

170
171 There are alternatives to using autoISF for FCL, as well. See [section 13.](#) for more info.

172
173 Notably [13.1](#) FCL using AAPS Master and Automations could be a much easier and more
174 error-tolerant way of stepping into FCL. In a clinical study with 16 participants about 80%
175 TIR was achieved without much tuning effort (source: see in [section 13.1](#)).

176
177 To close the circle to where we had started ([section 1.1](#)): A very time consuming pre-requisite might
178 actually be to first sort out your Hybrid Closed Loop, so your profile parameters are set „right“, and
179 your data really can serve as a blueprint for what, now, you would like your loop to do in FCL mode
180 (- see warnings at start of [section 4](#)).

181 Note that if you had used dynamic parameters or special Automations („loops inside the loop“) this
182 might have balanced some principal errors, but leaves you now without a good starting point as
183 you must get rid of these over-patches...

184
185 You will see also success stories of loopers who just jump into using more powerful tools, in
186 kind of a trial and error mode, and frequently add the latest add-on, or self-constructed
187 patch (often in form of an Automation) to counter-balance problems.

188 So, yes, you can also continue in that spirit. Resulting solutions may be good-enough. But
189 they tend to be unstable and not well-understood. That is a poor basis for managing arising
190 problems (-> fine tuning), and for adjusting to special situations (-> which setting to
191 temporarily change). But it certainly is an alternative avenue for the impatient, less
192 analytically, and more adventurous inclined.

193 In any case, PLEASE always observe the safety settings/instructions coming with the DIY
194 dev- variant of software you select.
195
196 One key safety measure every AAPS user going towards FCL should have in place is to set an **iob**
197 **threshold** (size a bit below what you used as a bolus for bigger meals in HCL) above which no
198 more SMBs can be given by your FCL.
199 This is an integrated feature of autoISF (from 3.0 version onwards), but can easily be implemented
200 in any AAPS variant via an Automation (as discussed for instance here;
201 <https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#iob-threshold>).
202 *How is it done in iAPS ??... IFTTT ??*