

# 1. Pre-requisites for Full Closed Loop

V 2.1

Please note that with autoISF 3.0 you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



## 1.1 Well-tuned hybrid closed loop

It is advisable to first establish a well-tuned hybrid closed loop before considering the transition to FCL. There are two important reasons for this:

- The UAM full closed loop requires a highly personalized (individual) tuning of settings, so the loop will give insulin mimicking YOUR successful hybrid closed loop mode.
- The UAM full closed loop comes with new parameters to be set and tuned. It would be problematic to set and tune several new parameters before the basics were tuned “right”. Errors could easily be balanced with counter-errors. This can work in single scenarios, but would create a highly unstable system, hard to re-calibrate better later.

## 1.2 Fast insulin (Lyumjev, Fiasp, Apidra?)

If the user does not bolus for meals, clearly a very fast insulin is needed so, upon realization of a starting meal-related glucose rise, the loop has any chance to eventually keep glucose in range (by common definition, under 180 mg/dl (10 mmol/l))

A modelling study (details see

<https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#fast-insulin-lyumjev-fiasp>)

can show in quantitative terms that **faster insulins**

- will result in significantly **lower** glucose **peaks** than slower insulins
- **tolerate** a couple of minutes **delayed** first meal bolus while not incurring unacceptable height of peaks
- minimize the effect on glucose peak from **different** carb load (**meal sizes**).

In conclusion, do not attempt FCL with other insulin than Lyumjev® or Fiasp®, unless, maybe, if you are on a very moderate to low carb diet. (According to [case study 1.2](#), Apidra® might work, too, but Humalog® would not work well).

### 1.3 Reliable insulin delivery from the used pump/cannula/insulin system

#### Good Tolerance of Lyumjev (or Fiasp): Occlusions threaten the function of the full closed loop.

It is very important to have an eye on the time a **cannula (or pod)** is in use (many find **48 hrs** to be the **limit**), and whether hard-to-explain glucose rises happen at ever increasing „fake“ iob (even before a 48 hr routine replacement). (See [case study 1.1](#): You easily lose 25% TIR that day)

It is absolutely contra-indicated to attempt FCL coming from leaking pods and associated erratic sensitivity swings that may or may not have been somewhat controlled and tolerable by dynamicISF or other measures when you were Hybrid Closed Looping,


#### Stable pump connection

In FCL you absolutely rely on your pump delivering, without any further delay, the much needed insulin, after any meal start.

Hence it is absolutely essential not to “miss” any problems from a lost Bluetooth connection.

An Automation similar to the one pictured here → should help recognizing eventual problems.

See also [case study 1.4](#)



### 1.4 Excellent CGM

You do not give a meal size-related bolus any longer. That leaves all insulination jobs to the algorithm! As glucose values are the very basis for this, please **inform yourself well about** how **your CGM** 1) principally performs 2) whether and how this may depend on data flow and intermediate apps you use 3) specifically, how and where any smoothing is done, and what this might imply for the ISF boosting method you will be using See for instance here:

<https://androidaps.readthedocs.io/en/latest/Usage/Smoothing-Blood-Glucose-Data-in-xDrip.html>

In AAPS Preferences/OpenAPS SMB/autoISF/Glucose source data for parabola fit, you must select between 4 options relating to your CGM (1 or 5 minute values, raw or smoothened).

Around meals, a stable Bluetooth connectivity is absolutely essential, too, so CGM, loop, and pump can do their job without losing more valuable time (see [case study 1.4](#)).

Then, but even more importantly in *all other* day and night *times*, the CGM should not produce any artefacts (jumpy values; see [case study 1.3](#)) that the loop could **misinterpret** as sign of a starting meal. Note that also calibrations could produce jumps.

76 But, in any case, a CGM with more scatter will make the loop lose more time, and lead to higher  
77 peaks and lower %TIR.

78

79 The best way currently is to use Dexcom G5 or **G6**, and to ensure via **overlapping** right and left  
80 arm sensor and transmitter utilization, that always good quality values can be used by the loop.  
81 Other ways are possible, but come with a lot of monitoring effort (via watch) and occasional time-  
82 outs for the loop.

83

84 Also FS **Libre 3** is useable. Observe info in the general section about autoISF  
85 (<https://github.com/ga-zelle/autoISF>) regarding how it should be implemented:

86

87 *Hint for users with 1-minute CGM (Libre 3) (status of 05 Dec. 2023)*

88

89 In AAPS 3.2.0.2 the glucose history only uses data at 5 minute intervals and ignores the  
90 ones in between. That has awkward consequences for the graph display and for the  
91 exponential smoothing method: Each minute a new reading is received the 5-minute pattern  
92 moves forward by 1 minute and a completely fresh and new subset of glucose values and  
93 times is used.

94 *Alternatives of determining the glucose acceleration based on the 1 minute data are still un-  
95 der evaluation. Once a promising method can be found and validated then an interim re-  
96 lease will be provided.*

97

98 As a libre user you have 3 alternatives:

99

100 • AAPS gets values from xDrip+ where you smooth and reduce the 1-minute data to 5-mi-  
101 nute data before sending to AAPS

102

103 • you should stick to the 5 minute mode in your Libre setup if possible

104

105 • AAPS gets values from Juggluco every minute but AAPS uses only the 5-minute subset  
106 for looping and for parabola fitting

107

108

## 109 1.5 Meal-related limitations

110

111 Setting up a full closed loop is relatively easy for people whose diet does not consist **mainly** of  
112 components with rapid high effect on blood glucose (more see

113 <https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#meal-related-limitations> )

114

115 Meals do not have to be low on carb (provided you use a fast insulin for your FCL)

116 Fat or protein rich diets, or slow digestion/gastroparesis, make things easier rather than harder for  
117 the full closed loop because late carbs nicely cover for inevitable “tails” of late action from SMBs  
118 needed around peak time.

119

120 Erratic consumption of snacks with fast resorbing carbs can be a problem.

121 In autoISF you can reduce this problem to some extent via one or two keystrokes from your  
122 AAPS home screen. While certainly being a deviation from the FCL idea(l), this would be  
123 one of the exceptional situations where you better do a quick “nudging” step from your “FCL  
124 cockpit”. Details see in [section 5.2.1](#) and [5.3.3.1 \(4\)](#)

125

126

## 127 1.6 Lifestyle-related limitations

128

### 129 Technically stable system

130

131 Full closed looping requires a 24/7 technically stable system, especially regarding

- 132 • reliable CGM signals
- 133 • Bluetooth stability with the pump (see [case study 1.4](#))
- 134 • keeping your phone in sufficient proximity at all times
- 135 • avoiding (or at least early recognition of) occlusion.

136 This requires a habit (or, unlikely, permanent attention to details) like keeping all components well  
137 charged and in close proximity; making cannula (or pod) changes always early enough to lower the  
138 risk of occlusion (see [case study 1.1](#)); having always potentially needed parts with you.

139 **Depending on your system, your experience with it, but also on your acceptance and general**  
140 **lifestyle, these aspects may or may not limit you.**

141

### 142 Preparing for exercise

143

144 To prepare for exercise (sports, heavy work), the normal protocol with a pump or hybrid closed loop  
145 is to take actions that reduce insulin on board prior to exercise

146 With your full closed loop, the algorithm is tuned to detect meals and to give you insulin to counter  
147 glucose rises automatically. Setting a high temp. target and lower %profile right away (effective al-  
148 ready around meal start) could be a problem.

149 Unusual activity levels therefore likely require **disciplined preparation** (especially **if you want to**  
150 **keep the need to snack during sports low**)

In autoISF you can reduce this problem to some extent via two or three keystrokes on your AAPS home screen. While certainly being a deviation from the FCL idea(l), this would be one of the exceptional situations where you better “flick a lever” from your “FCL cockpit” to have temporarily adjusted settings for the planned exercise. Details see in [section 6.2-6.3](#)

### Extra hurdles to establish FCL for kids

To establish and maintain a FCL for kids brings about some extra challenges if:

- Lyumjev is not available or well tolerated
- Hourly basal rate is very low, providing a poor basis for big SMBs
- Diet is rich in sweet components. With the typical low blood volume of a small body, strong tendency towards very high bg spikes!
- Going through marked changes of insulin sensitivity or of circadian pattern makes it difficult to keep the FCL appropriately tuned.

This problem is about the same in Hybrid Closed Looping. However, now you might expect miracles from the FCL. This is not going to happen. You still should try to set appropriate temp. changed profiles, that serve also as a basis for your autoISF FCL.

- Discipline is poor regarding keeping Bluetooth connectivity and infusion sites perfectly running
- Between kid and supervising parent it must be guaranteed, especially in the initial weeks, that an eye is kept on whether the FCL is working about as to be expected.

More see [section 7](#).

## 1.7 Time required for setting-up

Lastly, before enjoying a functioning full closed loop you need to have a period of a some weeks with some free time and „free head“ for set-up –. Can you get, in the time you are willing to invest, to a result that you consider good-enough is really the question. Depending on your „habits“, and which – if any - compromises (like doing cannula/pod changes more often, never starting meals when bg sits high ... ) are you willing to make (and everyday able to stick to), for the ease of not having to deal with assessing meals and bolussing for them?

While setting up your personal FCL using autoISF is a substantial project, there is no need to implement it fully in one step. There is nothing wrong to go in your well running Hybrid Closed Loop mostly, while switching to FCL only for dinners, for instance, or only for weekend lunches, as a start. Once you found feasible settings, you can expand to other meal times, and lastly towards figuring out your best strategies for challenges outside of meal windows, as we shall discuss in [sections 5. and 6.](#)

There are alternatives to using autoISF for FCL, as well. See [section 13](#). for more info.

193 Notably [13.1](#) FCL using AAPS Master and Automations could be a much easier and more  
194 error-tolerant way of stepping into FCL. In a clinical study with 16 participants about 80%  
195 TIR was achieved without much tuning effort (source: see in [section 13.1](#)).  
196

197 To close the circle to where we had started ([section 1.1](#)): A very time consuming pre-requisite might  
198 actually be to first sort out your Hybrid Closed Loop, so your profile parameters are set „right“, and  
199 your “old” data really can serve as a blueprint for what, now, you would like *your loop* to do in FCL  
200 mode

201 Note that if you had used dynamic parameters or special Automations („loops inside the loop“) this  
202 might have balanced some principal errors, but leaves you now without a good starting point as  
203 you must get rid of these over-patches (see also warnings at start of [section 4](#))..  
204

205 You will see also success stories of loopers who just jump into using more powerful tools, in  
206 kind of a trial and error mode, and frequently add the latest add-on, or self-constructed  
207 patch (often in form of an Automation) to counter-balance problems.

208 So, yes, you can also continue in that spirit. Resulting solutions may be good-enough. But  
209 they tend to be unstable and not well-understood. That is a poor basis for managing arising  
210 problems (-> fine tuning), and for adjusting to special situations (-> which setting to  
211 temporarily change). But it certainly is an alternative avenue for the impatient, less  
212 analytically, and more adventurous inclined.

213 In any case, PLEASE always observe the safety settings/instructions coming with the DIY  
214 dev- variant of software you select.  
215

216 One key safety measure every AAPS user going towards FCL should have in place is to set an **job**  
217 **threshold** (jobTH; size a bit below what you used as a bolus for bigger meals in HCL) above which  
218 no more SMBs can be given by your FCL.

219 This is an integrated feature of autoISF, from 3.0 version onwards (see [section 2.4](#)).