

7. Advanced HCL: Meal Announcement (MA) via a pre-bolus...

V2.5next

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



- 7.1 Frequently encountered hurdles
- 7.2 Getting ready to advance
- 7.3 Reduced pre-bolus
- 7.4 Tuning autoISF in HCL
- 7.5 Dealing with disturbances/ins. sens/resistance
- 7.6 Exercise management
- 7.7 Remote control (small children)
- 7.8 Other methods w/ meal announcement (MA)

[Available related case studies:](#)

[Case study 7.1:](#)

[Case study 7.2:](#)

See also [Case study 13.3](#) from a user of Boost

----- incomplete first draft -----

Originally it was planned to provide an extra section on FCL for kids here. However, we came to realize that the approach is no different for kids than already laid out. It just seems some implementation hurdles are significantly higher for implementing a safe FCL for minors.

For that reason, we like to explore in this section 7. how a **hybrid closed loop without carb inputs, using autoISF**, might get you to a **solution that removes most of the everyday burden** associated with having to co-manage meals.

This could (not only for kids) also be an **intermediary step, from which to progress into FCL as soon as a currently missing pre-requisite resolves** for you in the future.

7.1 Frequently encountered hurdles

Deficit making FCL difficult or unsafe	Bridging solution in Advanced HCL
Lyumjev or Fiasp (also in 50% mix w.slower insulin) not tolerated/too many occlusions; poor discipline re. scheduled infusion site changes	Difficult. Experiment with different insertion (site, depth, angle, cannula material), injection speed, site exchange frequency. (all difficult w/ pod pumps). Low carb diet
Poor discipline regarding keeping 100% Bluetooth connectivity	Given bolus (and profile basal running) will reduce the problem significantly; install alarm on (parent) phone
Leaking pods	(still a “no go”)

Deficit making FCL difficult or unsafe	Bridging solution in Advanced HCL
Jumpy CGM	Use strong smoothing (HCL does not rely on early aggressive action upon first signs of rising bg)
Very low hourly basal	No problem as HCL does not require super boosted SMBs
Erratic patterns of sweet drinks and snacks	Much less of a problem as a bolus is given with it

29

30 To establish and maintain *any loop* for kids brings about some extra challenges if:

- 31 • Going through marked changes of insulin sensitivity or of circadian pattern makes it difficult
32 to keep the FCL appropriately tuned.

33 This problem is about the same in all loops (Autotune, self-learning systems
34 etc...??). You should try to set appropriate (temp.?) changed profiles, that serve
35 also as a basis for your autoISF loop.

- 36 • Between kid and supervising parent it must be guaranteed, especially in the initial weeks,
37 that an eye is kept on whether the “Advanced HCL (MA) loop” is working about as to be ex-
38 pected.

39 On CGM might expand a bit ~ *pre-requisites, CGM, for kids: what works well ??*

40 If not G6 is used: (re CGMs, observe latest developments re 1 minute vs 5 minute readings, see
41 github/ga-zelle pls)

42

43 *include implications of (from wiki aaps):*

44 **Dexcom G6 or Dexcom ONE with xDrip+ Companion Mode**

45 The noise-level data is not shared with AAPS using this method. Therefore ‘Enable SMB
46 always’ and ‘Enable SMB after carbs’ are disabled.

47 **xDrip+ with FreeStyle Libre**

48 None of the FreeStyle Libre systems (FSL1, FSL2, or FSL3) broadcast any information
49 about the level of noise detected in the readings, and therefore ‘Enable SMB always’
50 ‘Enable SMB after carbs’ are disabled for all setups using the FreeStyle Libre.

51 In addition, many people have reported the FreeStyle Libre often produces noisy data. In
52 xDrip+ there are a few options to help with this:

53 **Smooth Sensor Noise.** In xDrip+ Settings > xDrip+ Display Settings ensure that Smooth
54 Sensor Noise is turned on. This attempts to apply smoothing to noisy data.

55 **Smooth Sensor Noise (Ultrasensitive).** If you are still seeing noisy data in xDrip+ you can
56 apply more aggressive smoothing using the Smooth Sensor Noise (Ultrasensitive) setting.
57 This will attempt to apply smoothing even on very low levels of detected noise. To do this,
58 first enable [engineering mode](#) in xDrip+. Then navigate to Settings > xDrip+ Display Set-
59 tings and turn on Smooth Sensor Noise (Ultrasensitive).

60

61

62
63
64
65
66
67
68
69
70
71
72
73
74
75
76
77
78
79
80
81
82
83
84
85
86
87
88
89
90
91
92
93
94
95
96
97
98
99

7.2 Ready to advance?

Switch off dynamicISF, forget what Autotune tries to tell you, and make sure your profile parameters are set right. Refer to guidance given in the HCL repo here:

Optimize meal management, notably watch that your ISFs are set right to deal with rising bg once your given meal bolus loses power.

With properly set ISFs, you should be able to expand allowed SMBV sizes to 120 minutes worth of basal.

Next, introduce a method that allows your loop take care of temp. insulin resistance from fats. (In the past, dynamicISF might have helped you for that.)

Now you have two options:

- Temporary increase of %profile via an Automation at signs of post-meal fatty acid resistance. See: (readthedocs)
- Or step into using the AAPS dev variant with autoISF, but make exclusively use of the dura_ISF component there.

Make sure your HCL now works at satisfying performance.

In the next steps, you try to get same performance, but with only a very rough idea, what you will eat (and no carb inputs)

7.3 Reduced pre-bolussing

Note that by continuing to make (albeit reduced) boli at meal start we remain in Hybrid Closed Loop (HCL) territory. Do **not boost** your SMBs *nearly as much as* enabled if you already work with the AAPS dev variant, and advised in section 02 *for FCL!*

- [Add link to ga-zelle "To pre-bolus or not to pre-bolus?"](#)

7.4 Implications for autoISF tuning (HCL vs FCL)

- **warning**, regarding implications for settings `bgAccel_ISF-weight` tuning should be factor 5 if not 10 milder ??)

- 100 • FCL => implications for remote parent monitoring ?
- 101
- 102 • Extra caution needed re. elevating the **SMB delivery ratio** ?
- 103 ○ The fixed 0.5 value in AAPS was installed also with a consideration on user/follower
- 104 (parent) set up and limiting potential problems from a bolus being initiated from both
- 105 phones in parallel
- 106 ○ As in FCL neither phone should initiate a manual bolus, this precaution seems not
- 107 really needed (dev – check)
- 108
- 109 Discuss risk from a “loose” bolussing habit in MA
- 110 * shooting from the hip against highs ...
- 111 * Need to learn interpret iob, insulin activity curve, prediction curves etc ...
- 112 7.5 Dealing with special situations / insulin sensitivity / Disturbances
- 113 7.5.1 manually, via temp%profile and TT Tuning aggressiveness
- 114 Key is that a sweet snack likely benefits from even more aggressive initial FCL
- 115 performance than the meals in your normal spectrum of diets require.
- 116 Therefore, you could set
- 117 • a higher **temp. profile%** and/or
- 118 • a **low temp. target** (76 for instance; this additionally helps maximize the first SMBs
- 119 that will automatically be triggered at detection of acceleration)..

120
121 • 7.5.2 via Automations a temp.elevated **bgAccel_ISF-weight** (see screenshot of my
122 Automation).

123 7.5.3 via custom buttons

124 The **DIY Cockpit** (shorten quoted [section 5.2.2.3](#) that follows!) allows to have buttons dedicated to
125 kindergarten hours etc that can in a time-uncritical manner activate special routines – **should be**
126 **very helpful to custom program Automations for kids!**

127 Grey extra DIY cockpit buttons for pre-programmed “responses”

128

129 Recognizing conditions for fully automatic handling by the loop may not be not possible, or come
130 too late for the loop to act on. Examples would be

- 131 • *exercise*: Minimum an hour before starting “the loop should know” to be able to lower iob
132 and elevate bg by the time exercise starts.
- 133 • *snacks*: High carb snacks, sweets, consuming ice cream or having a sweet drink comes
134 with the problem of even steeper glucose rises, but overall a lesser insulin need, compared
135 to major meals (for which we tuned our FCL according to [section 4](#)).

136 You need to research your snack habits (if any), and over time find out which settings in
137 the snack-related Automations work well.

138

139 **In everyday life you then just must press the related button in your cockpit** (which is
140 not time critical at all, except it should be clicked *latest* a couple of minutes after you took
141 the drink or snack).

142

143 **Installing the DIY cockpit button**

144 In the related Automation, just keep the “User action” box clicked at all times, and define in
145 the Conditions when you want to see that button available for cockpit use (see screenshot
146 above) => you will see that button offered.

147

148 Besides snacks, also any **other recurring special situations can be addressed via a**
149 **DIY cockpit button, and receive different aggressiveness up to a suitable iobTH**
150 **level.**

151

152 Over time you can have a big number of User action Automations, and keep them
153 “shelved” rather invisibly (clicked in-active, top left box) in your long list of potential
154 Automations. Even when active, they only show in your cockpit (bottom grey field of your

155 AAPS home screen) in the time slot you assigned as potentially relevant. radically need to
156 shorten this quote from chapter 5.2.2.3 -----
157

158 *End with something like:*

159 There are a couple of parents and kids pioneering this area, see also [Case studies 7.x](#).

160 We highlighted areas that would require some minimum compliance. In the end it comes down to
161 compare the achieved ease in daily use and achieved %TIR to how it was when hybrid closed
162 looping.

163

164

165