

## 6. Temporary Modulation for Exercise and lighter (In-)Activity V 2.3

Please note that with autoISF 3.0 you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them.

This is not a medical product, refer to disclaimer in [section 0](#)



Exercise management in autoISF builds on the „historic“ exercise mode of OpenAPS, and integrates the iobTH aspect for full closed looping.

### Preliminary remarks

This section is **no easy read** because it attempts to describe *all options* to deal with various types of exercise. Fortunately,

- you might need none of them
- you can set any of them up at your leisure, one at a time, for any of your occasional or regular exercise events
- then pick one or two of the described options, how to go about it.

This is a **toolbox** and, in green, it sketches even *not yet developed tools*, so, unless you are deeply development interested, you can just skip over these **green** passages until a new product version is announced that may include aspects relevant to you.

Looking at related case studies may be easier to digest. –

Please also report *your* experience by supplying a [case study](#).

Staying in contact with the related discord/github community should help greatly to find suitable ways to manage your type(s) of exercise.

## 6.1 Dynamic iobTH and sensitivity ratio in „exercise mode“

iobTH is a threshold you can set, above which AAPS will no longer deliver additional SMBs.

For exercise, we like to limit how high iob can go, therefore automatic “dynamic” reduction of your default iobTH (= iobMAX x iobTH%) is a benefit, notably as you can individually tune it.

Note: When transitioning to autoISF 3.0 from a previous version, de-activate (but keep for a while) the Automations you had for iobTH in previous autoISF versions. autoISF 3.0 totally changes how iobTH is accessed and modulated. (This can affect your automatic meal management, too).

In autoISF 3.0 and later, a default setting for iobTH is made in AAPS preferences, defined there as fraction (e.g. 0.6) of your set maxIOB:

/OpenAPS\_SMB/autoISF\_settings/Full\_Loop\_settings: iob\_threshold\_percent,

=> default iobTH = iobMAX x iob\_threshold\_percent

### 6.1.1 Manual (direct) iobTH modulation

„Manual“ routes to temporarily change iobTH would be

- changing the setting for the new parameter „iob\_threshold\_percent „
- or changing the setting for iobMAX

in /Preferences.

This is not a preferred route for temporary adjustment, because it would not automatically revert to default, after use.

A future improved FCL cockpit (-> [section 5.3](#)) may also give direct access to

- [override iobTH temporarily, at any point of time.](#)

## 6.1.2 Automations for temporary iobTH modulation

You can define Automations that sets a different iobTH under pre-defined conditions, or for a defined period of minutes to hours.

If your Automation has the User action box clicked, you get a grey button into your AAPS home screen from which you can activate that changed iobTH manually.

Note that this is the iobTH you tell the loop to use in place of the default iobTH



- it will still be modulated further if %profile and TT are set (see below).



Watch out for a potential stumbling block, because many Automations only work under the condition that no TT is already running.

## 6.1.3 Dynamic iobTH: Fully automated iobTH modulation via setting a temp. glucose target

Note that in AAPS preferences, you need to set high TT raises sensitivity = TRUE.

Then, setting a temporary glucose target (TT), modulates iobTH the same way as it modulate sensitivity (ISF). sensitivity).

When also the **exercise button** is ON (lit yellow), iobTH gets reduced particularly strong, and ISF weakened (as desired for exercise).

That effect is the stronger (**ISF gets the weaker, iobTH the lower**), **the lower you set the half-basal exercise target** for your exercise mode in AAPS/preferences/OpenAPS SMB:

The following table shows, for a profile target of 100 mg/dl, the effects of your set:

- half\_basal\_exercise\_target (set in AAPS/preferences/OpenAPS SMB)

Choose a low number if you later want a high dynamic range of sensitivity modulation

- and your current exercise TT.

Set your TT (on the day you do the respective exercise) with an eye on how you wish sensitivity auto-adjusted. Higher TT = lesser insulin delivered

Note that:

- $\text{temp. basal} = \text{profile basal} * \text{sens.ratio}$

*Example: At a half-basal\_exercise\_target of 120, setting a TT of 120 gives only half (0.5) of profile basal (hence the name of the parameter)*

- $\text{temp.ISF} = \text{profile ISF} / \text{sens.ratio}$

- $\text{temp.iobTH} = \text{set iobTH} * \text{sens.ratio}$

The following table gives some examples for resulting sensitivity ratios.

Half basal ex.target	180	150	120
TT	sens.ratio	sens.ratio	sens.ratio
100	1	1	1
120	0,8	0,71	0,5
140	0,67	0,56	0,33
160	0,57	0,45	0,25
180	0,50	0,38	0,20

The exact calculation for any combination of profile target, set TT, and half-basal\_exercise\_target is given in [section 3.3](#)

You do not really have to deal with these details, though. Just sit back in your cockpit, and watch the effects of various inputs on iobTH and %sens on your AAPS home screen

Try to determine good settings for the kinds of exercise that you frequently engage in.

Later, just press the grey DIY button, [or the TT button, and make an exercise-related entry there \(sections 6.2 or 6.3; case study 6.2\)](#). This will automatically switch the exercise button to yellow (ON), and it lowers basal and iobTH as (in your experience) suitable.

Note that

(1) setting a TT often shuts out other Automations. Choose the duration wisely (and also the sequence, in which all your Automations are listed).

(2) (assuming, you use the even/odd differentiation for SMB on/off:) Consciously decide whether you set an even or an odd numbered TT.

- Pick **odd**, if you do not want SMBs during exercise. (Despite you softened ISF, SMBs still might „attack“ a sports snack too strongly).

However, odd cannot be set too early, when your meal digestion still requires SMBs. Likewise, you might want the option for a few automatically delivered SMBs against unforeseen spikes (e.g. from excitement) also later. An **Automation** that switches from odd to even for a couple of minutes might sneak in a desired SMB or two .

However, you are probably out of luck because an already set odd (or any) TT would preclude such Automation from kicking in. Then you need to develop additional ideas, another detour, like to first define an Automation that briefly shuts your odd TT down.

- Working with an **even** TT can sometimes be preferable, notably of course if your exercise is one that can get you totally excited, with glucose spikes. While this mode generally does allow SMBs, the loop softens the ISF (by the sens.factor like in the table given above), and will temp. shut SMBs down, when **iobTH** (which also got lowered by the sens. factor) is exceeded.

Whether odd or even TT is better depends on the kinds of exercise you are doing, and probably depends on the protein and fat load of your meal and snacks, as well.

(3) Timing can be critical as to when you do this exercise announcement, especially relative to a preceding hi-carb meal. Then you want the reduced iobTH in place latest after you received the first SMB. See [section 6.4](#) and [case study 6.2](#)

You always can see the valid iobTH your loop is working with in your AAPS home screen, next to the current iob status.

You can use any of the above discussed methods, or also the one that now follows in [section 6.2](#), to further tweak iobTH temporarily, should you see a need.

## 6.2 Temporary % profile switch

A complementary measure you can take from the AAPS home screen is to set a **reduced temp.% profile sensitivity**.

This setting would **multiply** with the results in above table and further reduce basal and iobTH (whenever exercise button AND profile button both are yellow).

Temp. reduction of basal will proportionally also reduce the *max. allowed* size of SMBs (which is two hours worth of basal x SMB\_range\_extention, see [section 2.1](#))

Note that the **time windows** for doing this profile switch (which was the main ingredient of going into exercise in hybrid closed loop) can differ from your TT-related exercise settings. Using all available tools then allows a nearly surgical approach to what you want to achieve for and during your favorite exercise(s).

- Often the %profile modulation is used for several hours if not days to accommodate “long waved” sensitivity swings (See e.g. in [case study 6.2](#)).
- Instead, or even additionally, the percentage might be modified for just a couple of minutes, or for one special snack or meal duration, to “nudge” the proportionally modulated aggressiveness of the FCL (see [section 5.2.3](#)).

You can prepare yourself for anything you see coming up, or potentially coming up, in your daily life, so, from the comfort of your cockpit you get ready for it within just a second or two, doing a few „clicks“.

## 6.3 Managing exercise via Cockpit inputs

### 6.3.1 Basic Settings for Exercise

Coming from FCL with no TT set (both top fields, TT and exercise, are grey), you best prepare for an intended exercise by **pressing the TT field** of your AAPS main screen (your looping cockpit; presented in [section 5.2](#)).

There, you can **freely select** TT and duration.

Alternatively, you can press on **one of 4 offered exercise presets**. (Note: This, and many other – in this green color - described cockpit features are yet to be developed)


When you do either one, the exercise button in the top middle of your AAPS main screen will turn yellow: (It also does turn yellow, or remains yellow, whenever you **make a new selection or input in these fields (or when you** just press on the exercise button, when a TT is set.)

### 6.3.2 „Dynamic“ exercise mode off = traditional AAPS exercise mode (YGY)

When the dynamic exercise mode is off, you still have the instruments for *exercise management just as you always had it in the past* = a combination of manually softened aggressiveness via setting a temp. %profile change, and orienting corrections towards an elevated TT.

By selecting an odd numbered TT you now have the *additional option* to shut SMBs temporarily off, too.

*If improved cockpit is launched, the top part of the dialogue box looks about like this when the exercise field is grey:*



% profile can be changed:

- either here => neighboring %profile button turns yellow too (with the % info on it); or
- under the %profile button; or
- it had already been changed using AAPS / Action / Profile switch

In all 3 cases, you see the number < 100 or >100 in the middle of above table, on a yellow colored field, too.

In this “YGY” mode, the % temp. set profile is the effectively applied “%overall” sensitivity

TT and duration can be entered or changed (= traditional mode to set exercise targets).

If there is a desire to try, for the **remaining duration**, a different iobTH or bgAccel\_ISF-weight, this can be overridden in the table; field turns yellow, and the algorithm uses temp. iobTH and/or temp bgAccel\_ISF\_weight as modified in the exercise button (and reports this also in the SMB tab).

*In the dialogue box pictured above, 70% profile was set for 3 hours, and the default iobTH of 60% \* 10 U was cut by 50% down to 3.0 U.*

The effective iobTH shows also in the AAPS home screen, next to the actual iob (e.g. „1.2 U < 3.0 U“)

The remaining duration shows below the table (in the example: 134 minutes and counting down).

TT and % profile will also show on the yellow labels of the neighboring %profile (left top of AAPS home screen) and TT (right side), respectively.



The middle (exercise) field remains grey because the automatic sensitivity tuning (that would use TT and half-basal exercise target) are off.

### 6.3.3 Dynamic exercise mode ON (GYG or YYY)

By pressing the yellow exercise button on the AAPS home screen, *you have the option* to switch the **dynamic exercise mode ON**, in which case the middle field/exercise button of your AAPS main screen will go from grey to yellow.

In a version update you could do your setting for the upcoming exercise under the **dialogue box of the TT button**

Then, when you look into the exercise button in the middle of your FCL cockpit the dynamic exercise mode will automatically be „ON“, and all entries made:

activity	TT	dur	%sens	iobTH	bgAccel	%overall
mtb	171	180	100	4.0-11	0.16	67

● Mode set to run for 134 more minutes

● Mode starting after meal when iob > iobTH

The data for the kind of exercise (*here mtb; could also be n/a or ?*) are coming from prior selections made in the dialogue box of the neighboring TT field. There, as well as in this window here, the resulting iobTH and bgAccel\_ISF\_weight are shown. Also the overall aggressiveness (% overall insulin sensitivity factor) is calculated.

The **middle field** of the table in this dialogue box, **% profile** either picks up the % set under the %profile button, or an input can be made here, in the exercise button domain, which will:

- turn the neighboring %profile button on yellow and show that inputted % on it, too
- be multiplied with the result from the exercise mode settings per se, and change the % overall, accordingly.

So, if this middle field of above table (dialogue box of exercise button) contains a figure other than 100, the input field becomes yellow, and you are operating with a combination of traditional PLUS new exercise mode (with all three top buttons of your FCL cockpit yellow). This maximally will soften aggressiveness, for which you get an idea by the last calculated figure.

247

The mode is either running already (*for another 134 of the total 180 minute in the picture*) as also the label on the neighboring yellow TT field will show 171 (134, and counting down),

Or (*see at the red dot in picture above*), it is scheduled to run, after insulination for a started meal surpasses iobTH (*as in table*).

Note that, when the TT expires or is changed, your overriding input (if you made any) is automatically erased, forgotten.

254

#### 6.3.4 Dynamic exercise mode ON plus %profile change (YYY)

256

The **middle field** of the table in the dynamic exercise mode dialogue box (see above), **% profile**“ either picks up the % set under the %profile button, or an input can be made here, in the exercise button domain, which will:

- turn the neighboring %profile button on yellow and show that inputted % on it, too
- be multiplied with the result from the exercise mode settings per se, and change the % overall, accordingly.

So, if this middle field of above table (dialogue box of exercise button) contains a figure other than 100, input field becomes yellow, and you are operating with a combination of traditional plus new exercise mode (with all three top buttons of your FCL cockpit yellow). This maximally will soften aggressiveness, for which you get an idea by the last calculated figure.

267

It is advisable to find good settings within the dynamic exercise mode and NOT use profile switches on top – unless the profile switch is meant, also outside of the temporary exercise context, to provide for other, „longer waved“, health or hormonal situations.

Also, that middle field offers easy access for temporarily tweaking the aggressiveness without immediately changing core settings like the half-basal-exercise target etc.

## 6.4 Option to pre-set for 4 kinds of exercise or meals (for 1 button operation)

### 6.4.1 iob\_threshold\_percent

In AAPS preferences/OpenAPS SMB/autoISF settings / Full Loop Settings, the default iob\_threshold\_percent used for the normal meal spectrum is defined.

In an updated later autoISF version you might be able to differentiate there for up to 4 meal clusters (see next section)

### 6.4.2 Pre-settings for (up to) 4 kinds of exercise:

In AAPS preferences/OpenAPS SMB/autoISF settings / Full Loop Settings: follows next input fields for pre-settings you can define for (up to) 4 kinds of exercise:

The following table gives an example of settings you may find well-suited for 4 of your favourite exercises

#1-4	give name (max 3 characters)	duration for TT ( min)	TT (AC) (mg/dl)	% profile	iobTH	bgAcce:weight	Approx % ins reduct.
1	wlk	60	111	100			
2	grd	120	131	90			
3	bik	300	151	90			
4	mtb	180	171	70			

Input fields (during tuning phase to determine good settings) are only the columns 2-5.  
The last 3 columns will be calculated from TT and %profile inputs, using also the half-basal exercise target and the default weight setting. In this setting.  
The last is only an approximation to get a feel for a reasonable setting of the other parameters. Here in preferences they should never be overridden, but TT or % profile should be adjusted to reach desired result when tuning for FCL.  
Likewise, you find tables to make pre-settings for meals and for hypo treatments:

#### 6.4.3 Pre-settings for (up to) 4 kinds of meals:

In AAPS preferences/OpenAPS SMB/autoISF settings / **Full Loop Settings**: follows next:  
Input fields for pre-settings you can define for (up to) 4 kinds of meals. For instance:

TT# 1-4	give name (3 letters)	TT (Eating Soon) (mg/dl)e	Duration for TT (min)	iobTH (0---130% and < iobMAX)	bgAcce factor 200...0%	
1	hiC	72	120	110	110	
2	loC	74	180	67	67	
3	piz	76	300	100	100	
4	snk	78	60	100	50	

Input fields (during tuning phase to determine good settings) are all columns  
Difference in TT is fairly unimportant (unless you do not give a name and memorize the set TT number instead, for which meal type it codes.  
Logic why not having a % profile column here: %profile switch should be set extra, potentially for another time period (e.g. „reserved“ for periods of exercise, or for entire days of altered insulin sensitivity, for instance due to illness, fasting, extensive sports week.)

310

311 6.4.4 Pre-settings for (up to) 4 kinds of Hypo treatment:

312

313 In AAPS preferences/OpenAPS SMB/autoISF settings / **Full Loop Settings**: follows next:

314 Input fields for pre-settings you can define for (up to) 4 kinds of HYPO treatment. Example:

TT (ES) (mg/dl)	give name (3 letters)	TT (AC) (mg/dl)	Duration for TT (AC) (min)	bgTH (mg/dl)	
1	Hy1	131	55	none	
2	Hy2	131	55	200	

315

316 Input fields (during tuning phase to determine good settings) are all columns, 2-5.

317 Choosing an odd-numbered TT is recommended as it can shut-out SMBs (with the

318 appropriate setting in preferences/Open APS SMB/autoISF settings/smb\_delivery

319 settings/"enable alternative activation...".

320

321 Those of us who tend to over-treat hypos may prefer to set Hy2 (unless for night snacks->

322 Hy1): Reverting to standard loop aggressiveness with SMBs after/if a certain bg level

323 („threshold“, similar to our iobTH for meals) is surpassed, and we want our loop to react

324 again with SMBs before the set duration expires.

325

326

## 6.5 Mastering Exercise after a Meal

In Hybrid Closed Loop, we gave less insulin at meals (a reduced bolus) before exercise.

Since we now get our meal insulin automatically from the loop, we would have to at least somehow tell it that exercise follows this time.

Simply setting an exercise profile *before* the meal would make our full closed loop too weak in the "treatment" of the first glucose rise. What we want is, to get our (already, compared to HCL, delayed) meal insulin delivered as fast as possible by SMBs. It just should be capped at the desired iob reduction.

### 6.5.1 Manual mode requires 2 user interventions

What we can do, is (1) **reduce** the **iobTH** (e.g. by one third).

- In the example we were using, this would mean to reduce by 2 U to  $iobTH^* = 4U$ .
- Do that estimate for your data, and think back how you did bolus reduction in hybrid closed loop before same exercise.
- Likewise, you can use your profile ISF, e.g. 30 mg/dl/U and „translate“ by how much (  $2U * 30 \text{ mg/dl/U} = 60 \text{ mg/dl}$ ) this „pulls you away from going into a hypo“.
- Using your IC (e.g. 8g/U) you can also translate the iobTH reduction (2 U) into a „snack equivalent“ (  $2U * 8 \text{ g/U} = 16 \text{ g}$ ) that you „replace“ by thinking ahead and „budgeting“ for some exercise with your iobTH modulation.

In this senario, our loop delivers SMB insulin as fast as always, only that when the last SMB has passed the iobTH, the loop only has elevated %TBR to work with, meaning it cannot raise iob by much any longer. This provides an elevated glucose level on which we enter exercise, and saves us hypo danger or snack need (as calculated in above examples).

After this reduced iobTH is reached, step (2) must follow = an increased exercise **bg target** is set (see [section 6.2](#)).

The problem with this approach is that it requires **two** user interventions, first **setting the lower iobTH**, later (**and this *in a time-critical manner***, after iobTH is exceeded), to **input a exercise TT** or activate a related setting.

To eliminate this problem, the following refined solutions are suggested:

#### 6.5.2 DIY cockpit: Using pre-set meal / exercise settings from a User action Automation

The „DIY cockpit“ user interface allows a *one-step* setting for meal + exercise that can be selected in time-uncritical fashion, any time before the meal starts.

Summary from detailed example given in [case study 6.2](#):

A sequence of 3 Automations must be set up, of which only the first one must be manually triggered, in just one time-uncritical key stroke from the AAPS home screen.

The others come on automatically when the respective Conditions are met.

Automation #1 provides, for a meal that precedes exercise, the full loop aggressiveness, but makes sure that this aggressiveness stops immediately after a (reduced) iobTH is exceeded. The reduced iobTH ensures that not too much insulin is on board for exercise after the meal. Also it provides an elevated bg level at (re-)start of exercise.

In this Automation, the box “User action” should be permanently ticked. This will automatically provide a grey button on the bottom of the AAPS home screen (“DIY cockpit”) that can be freely named (= headline of Automation #1).

When the reduced iobTH is exceeded, two things need to be provided :

- (1) a milder running FCL (reduced exercise %profile, after the meal rise had been managed based on 100% profile boosted further by bgAccel\_ISF driven full loop aggressiveness) => Automation #2 sets e.g. 70% profile and ends TT

(2) setting an exercise TT (not possible with Automation #2. But *after* it ended the TT, an Automation #3 can immediately follow and set the desired exercise TT=125 (which implies the exercise mode

Note that Automations 2 and 3 are fully automatic, no User action is involved. See [case study 6.2](#) for an example

Should during the exercise a need arise to modulate the loop aggressiveness (jobTH, effective ISF) that can be done within 1-2 seconds also right from the AAPS home screen („FCL cockpit“) by setting a higher or lower temp. %profile, and/or by setting a higher or lower temp. exerciseTT.

To make the loop temporarily act a bit more aggressive, switching the exercise button OFF (from yellow to grey) could also be considered

[Defining User action - Automations to build your FCL cockpit](#)

If you want to develop your **DIY User Interface**, make sure you define suitable settings that reflect **your** personal insulin sensitivity and data patterns.

As mentioned in other places, Automations can be tricky as to whether they actually will ever work, because the loop goes through the exact **sequence of all your active Automations**, and might be switched into a direction that no longer is compatible with the conditions that must be a given, for the Automation you think that should kick in.

To have a clean AAPS home screen (and also to prevent unnecessary accidental activation), define reasonable time windows for each of your shelved special routines, or keep them entirely dormant (de-activated) in the list of all your Automations, and activate them only for the day when you think you might need them



### 6.5.3 Improved cockpit: Using pre-set meal / exercise combination from TT dialogue box

The improved „FCL cockpit“ User Interface (when available) also allows a one-step setting for meal + exercise that can be selected in time-uncritical fashion, any time before the meal starts.

It manages the meal with an appropriately reduced iobTH, and is programmed to automatically activate the exercise settings when iobTH is exceeded:

If in addition to meal, one of the 4 pre-programmed exercises is also selected from the bottom of the TT dialogue box, (for example, in case of biking after a hi carb lunch, hiC + bik, see [section 5.3.3.1](#).) then meal gets superseded /overridden with condition „duration = until when iobTH is first time exceeded“. Plus, that is the other important point, the activity-related reduced iobTH is taken over for the meal, too.

All this happens from the AAPS home screen and associated dialogue box from the TT field there.

Actual valid settings can at any time point be seen in the AAPS home screen (see [section 5.3.3.1](#) on extra data fields, above).

When your FCL is in this meal + exercise mode, you first see at the TT field ([section 5.3.3.1](#)) of your AAPS main screen:



... and when iobTH is the first time exceeded, this automatically switches to:



438

439 That 1-step setting can either be freely done according to [section 6.2](#)

440 Or you can just press one of your frequent meal [and](#) frequent exercise „codes“, as described  
441 in [sections 6.3.2 and 6.3.3](#)

442 *Example: For mountain biking after pizza lunch* press two buttons, *piz* and *mtb*, in the  
443 dialogue box of your AAPS home screen's TT field. That's all (...after, one time, you figured  
444 out what settings suit that scenario, and you put it into /preferences, see [sections 6.3.2 and](#)  
445 [6.3.3](#)).

446

447 6.5.4 Laissez-faire alternative

448

449 You could also just use an exercise setting and accept a reduced loop aggressiveness  
450 already before meal start. You would go a bit higher in your glucose peak. As, in principle, a  
451 higher glucose level is desirable for starting exercise, this can be a viable route, too.

452

453 This depends on your meal's carb load also, and should be viable if you do the often  
454 recommended protein-rich meal before exercise.

455

456 Note that making the exercise setting after meal start is problematic in case the first  
457 SMBs already exceeded the iob you see as limit for starting your exercise (which is  
458 not the limit for the meal *per se*).

459

460

461

## 6.6 Activity Monitor

An optional feature for times without serious exercise, but still suspected **effects on insulin sensitivity (max +20% to minus 30%)** is the **activity monitor**.

It can be generally activated under /preferences/OpenAPS SMB/Activity modifies sensitivity)

If the user

- has scaling factors set there (in preferences/OpenAPS SMB/Activity modifies sensitivity)
- has **no TT running**
- (and, regarding nighttime: did not opt for „ignore\_inactivity\_overnight“)


then AAPS automatically modulates for sensitivity changes **based on movement intensity** for the last minutes to 1 hour time frame.

**Personalized tuning of the two scaling factors** is necessary in your FCL set-up phase. For details see [section 3.4](#).

The Activity Monitor can also be used (overridden/ used for tuning the scaling factors) from a dialogue box (*if already launched*) coming up from the exercise button (top middle of AAPS home screen).

**Note that Activity Monitor only works if no exercise (or other) TT is active** (which would influence insulin sensitivity ratio much stronger than the tweeking done by the Activity Monitor, for slighter everyday effects).

487 In this dialogue box (connected *in a future update* with the exercise button), the two scaling  
 488 parameters (set as default by the user during initial set-up in preferences) are displayed, and  
 489 can be temp. over written. (These settings will expire and revert to default as set in  
 490 /preferences, whenever the Activity Monitor closes (goes auto-off, or is pushed off)).



491

492

493 The resulting sensitivity effect is the roughly expected effect of requiring >100% insulin if  
 494 moving around a bit (activity), or needing a lesser %age when being very stationary.

495 It is displayed in the right side column of the dialogue box (*if already launched*) to give the  
 496 user a feeling for the expected effects from her/his „weight“ inputs.

497

498 The exact impact is calculated by the loop and shown on top of the autoISF results in the  
 499 SMB tab (every 5 minutes).

500