

11. Emulator on your AAPS smartphone

V 2.8

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



11.1 Installing the emulator on your smartphone

- 11.1.1 Download QPython 3L
- 11.1.2 .py files in phone internal memory
- 11.1.3 .config file
- 11.1.4 .vdf files
- 11.1.5 Customization of output table

11.2 Checking loop decisions on the smartphone

- 11.2.1 Principal purpose
- 11.2.2 Generating **results table**
- 11.2.3 Analyzing results

11.3 Options available on i-Phone (Trio or iAPS)

11.4 Real-time checking a „what-if“ question using **speech synthesis**

[Available related case studies:](#)

Case study 11.1: none available yet

Skip what is in green writing:

= Drafted fragments or not implemented ideas.
Please contribute, or wait for update with the missing info

The emulator **on the PC** was presented in [section 10](#) as

- a very good tool for making your **initial** tuning for a meal spectrum (= when weighing different effects over the entire course of time after each meal, and for a variety of your meals).

A very useful *additional* tool, is the emulator running **on the AAPS loop smartphone**


- for a quick look how autoISF triggers SMBs after starting a meal ([section 11.2](#)).
- and especially for real-time checking „what-if“ you implemented a specific change idea (developed on your PC, or after analyzing many SMB tabs) ([section 11.4](#)).

Regarding **i-Phone** options, see section 11.3

The emulator for the AAPS phone is described here: <https://github.com/ga-zelle/APS-what-if>

There you find the files needed to download on phone or PC, and the primary instructions.

Join <https://discord.gg/n3tD5eXExC> for seeking (and giving) help with the Emulator set-up or use, and to share experience.



Sketch of Running the Emulator on AAPS Phone

Github/ga-zelle /
APS what-if

11.1 Installing the emulator on your Android smartphone

Note that iOS based autoISF variants cannot use the Emulator on an i-Phone. However, some tabular outputs of ISF-factors are available (see [section 11.3](#) to be provided by iAPS user)

11.1.1 Installing QPython 3L

On your smartphone, go to Google Playstore and download the QPython 3L app. Put the app icon next to your other looping related app icons on the main screen of your smartphone.



With many QPython 3L versions, phone and Android OS versions etc around, you might run into problems and may need to consult detailed installation instructions from [Github/ga-zelle / APS what-if](#), or seek advice via <https://discord.gg/n3tD5eXExC>

By long pressing on the QPython 3L app icon, go into “app info” and make the settings like for all your other looping related apps, so they do not get killed all the time by energy savings routines.

11.1.2 Copy .py files from your PC into your phone's QPython/skript3 folder:

Connect your phone to your PC for USB data transfer

67 Locate your .py files on your PC (in Emulator file).

68

PC: - your path to the Emulator files may differ -

Name	Änderungsdatum	Typ	Größe
__pycache__	24.05.2024 10:05	Dateiordner	
5minute_emulator_std.config	24.05.2024 22:06	CONFIG-Datei	1 KB
determine_basal.py	17.05.2024 17:14	Python File	149 KB
emulator_batch.py	25.05.2024 10:59	Python File	21 KB
emulator_core.py	25.05.2024 10:59	Python File	164 KB
emulator_GUI.py	07.05.2024 20:04	Python File	42 KB
find_string_batch.py	07.05.2024 19:31	Python File	6 KB
find_string_core.py	07.05.2024 19:31	Python File	8 KB

69

70 **Copy** all Python related files **except** emulator:GUI.py from your PC over the internal memory / QPython

71 / **Scripts3** of your phone:

Name	Typ
bluetooth_chat.py	Python File
determine_basal.py	Python File
emulator_batch.py	Python File
emulator_core.py	Python File
emulator_GUI.py	Python File
find_string_batch.py	Python File
find_string_core.py	Python File
hello_world.py	Python File
pip_console.py	Python File
say_chat.py	Python File
say_time.py	Python File

72

73

74

75

76

77

78


79

11.1.3 Put configuration file and noChange.vdf into the phone's AAPS logfiles folder

With your phone connected to your PC for data transfer, retrieve **5minute_emulator.std.config** (or 1minute_emulator.std.config if you use Libre3/1minute) on your PC, (picture, left), and put a **copy** into the phone into the logfiles (not the QPython!) folder (picture, right):

- Internal memory/AAPS/logs/info.nightscout.androidaps

This config file contains the programming as to which hours of day there can be speech messages. Default 07-23 h ("your" time, not "Z"scale). How to change, see in [section 11.4.3](#)



11.1.4 Put noChange.vdf into the phone's AAPS logfiles folder

Retrieve the **noChange.vdf** on PC in the parallel Studies file

- C:\...\Documents\ Looping\ Logfiles_Emulator\ Emulator_Studies

...and put it on the phone, also into the folder (picture, right):

- Internal memory/AAPS/logs/info.nightscout.androidaps

Later, in [section 11.4.1](#), you will add also (yourChange).vdf files into the same folder. One is already included, line under noChange.vdf, in picture above, right.

107 11.1.5 Customization of output table

108

109 The table should contain the most relevant information that can be displayed reasonably on most
110 smartphone screens

111 Consult [https://github.com/ga-zelle/APS-what-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)

112 [if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)
113 [phone.pdf](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf) in case you see a need to customize .config files and output tables for you.

114

115

116

117

118

119 11.2 Inspecting loop decisions on your smartphone

120

121 11.2.1 Principal purpose

122

123 The result table from the emulator on your smartphone allows you much easier insight
124 than the SMB tab can offer into current and recent determinants of given SMB sizes (e.g.
125 which of the 4 autoISF categories contributed).

126 So, if you would have benefitted from a bigger (or smaller) SMB at times where, say,
127 bgAccel_SF was the dominant factor, you would increase (or decrease) the associated
128 weight. Before actually making such changes, look in other lines of the table to estimate
129 how this would influence decisions in other time points.

130 Testing your hypothesis on the PC ([section 10.3](#)) would allow multi-day multi-(kind-of-)-meal
131 judgement on feasibility of your tuning idea.

132 Before firmly deciding on a settings change, it is advised to run on your phone a what-if
133 emulation ([section 11.3](#)) using the less aggressive settings for your active loop run, and the
134 more aggressive defined in (yourChanges).vdf.

135

136

137 11.2.2 Generating a results table (for last 75 minutes' loop decisions)

138

139 Note: 1 minute **Libre3 users** will get data only on a significantly shorter time span.

140 15 minutes is really too short to analyze what is going on, and Libre3 users probably should

141 mostly use the “what-if” part, see [section 11.4](#). That part is principally not impacted (except,

142 cutting one 5 minute change into 5 very tiny changes, often will make it hard to see and

143 interpret “what-if” effects).

144

145 Display setting on your phone should be set for automatic switch between landscape and

146 upright viewing (depending how you hold your phone).

147 Reducing selected font size will *not* help to get more table info on your phone screen, or to

148 avoid broken/double lines. Go to [11.1.4](#) if, in the following, you are not happy with lay-out.

149

150 1). On your main phone screen, press the QPython 3L app icon:

151

152 The first black screen then popping up asks to make a **language** selection

153 In case you don't get an alpha-numerical input field (with <- enter button), touch the

154 upcoming black screen

155

156 All black screens have a keypad at the bottom:

157

158 2). Press your selected number and <- (for enter).


159 3). Then proceed in the same screen

160 to make your capital-”N” entry

161 and finally “ <- “ for Next (see illustration)

162

163 In the following the instructions from the developer how to navigate through the screens:




After some Android12 update the previously used GUI dialogues no longer worked and I had to create a keyboard based version similar to the typical telephone dialogues “for option X dial 2”. This dialog system consists of two parts:

- The top part has numerical keys for each option that can be selected. One of those options has **(default)** at the end of the line which indicates that this option is the current selection. If you enter a different number the dialog screen is redrawn and that indicator moved to your new selection. Once your intended selection is OK you focus on the ...
- bottom part, which has letters as keys for the action to be selected. Again, the **(default)** indicator highlights what would be done next if you just press enter without any digit or any letter. Those actions typically are Next, Test and Exit.

The first dialogue is used to select the **language** for the speech synthesis.

Select „Test“ to listen to a sample speech synthesis.


Select „Next“ to proceed.



16:58 39%
 ← No. 2 NEW CTRL
 List of Languages:
 1-Dieses Smartphone spricht Deutsch (default)
 2-This smartphone speaks English
 When done, list of actions:
 N-Next (default)
 T-Test
 E-Exit
 Enter key for option or action: 1

4-5). Repeat steps 2) and 3) with the next screen:


The next dialogue is for selecting your **variant definition file**. All vdf-files found in the logfile folder will be listed.




16:59 39%
 ← No. 2 NEW CTRL
 List of vdf-files:
 1-CR50pct.vdf (default)
 2-acce_pm0p03.vdf
 3-aimax4.vdf
 4-noChange_ai227.vdf
 5-noChange_aisf226.vdf
 6-bgAcce_dura_quadrat.vdf
 7-bgAcce_dura_linear.vdf
 8-TT_boost_50pct.vdf
 When done, list of actions:
 N-Next (default)
 E-Exit
 Enter key for option or action:

164 6-7). Repeat steps 2) and 3) with the next screen:

The last dialogue is used to select your favourite **configuration file** with the content discussed in the preceeding section. After you activate the selection you get an informational message of how many columns the selected tabular output will occupy.



17:01 40%
 ← No. 2 NEW CTRL
 List of config-files:
 1-vary_settings_normal.config
 2-vary_settings_gazelle.config (default)
 When done, list of actions:
 N-Next (default)
 E-Exit
 Enter key for option or action: N



17:01 40%
 ← No. 2 NEW CTRL
 List of config-files:
 1-vary_settings_normal.config
 2-vary_settings_gazelle.config (default)
 When done, list of actions:
 N-Next (default)
 E-Exit
 Enter key for option or action:
 Total width of output table is 72
 Press Enter to continue

8) Before activating the “N(ext)” selection now is a convenient time to rotate the phone to landscape to prepare for the multi column result table display.

165

166



167

168

169 9).Now a table comes up detailing the loop decisions on SMB size for the last (15 * 5 =) 75
 170 minutes.

171 In case you forgot to do it in step 4), turn the phone now 90 degrees for landscape format;
 172 in this case, give it 5 minutes for the format to straighten out (after a new value has arrived)

173 In case you want other, or less, columns, see [section 11.1.4](#).

01:39 84 58%

← N

1

2

3

4

5

6

7

8

9

III

UTC

eff.

acce

bg

pp

dura

-----ISFs-----

insulin Req

---SMB---

time

bg

IOB

iobTH

ISF

ISF

ISF

orig

prof

emul

orig

emul

orig

emul

23:23Z

97

0.93

4.0

-0.28

0.65

1

1

375

150

375.0

0

0

0

0

23:24Z

94

0.9

4.0

-0.44

0.61

1

1

375

150

375.0

0

0

0

0

23:25Z

92

0.89

4.0

-0.21

0.58

1

1

375

150

375.0

0

0

0

0

23:26Z

92

0.85

4.0

0.25

0.58

1

1

375

150

375.0

0

0

0

0

23:27Z

91

0.84

4.0

0.38

0.7

1

1

375

150

375.0

0

0

0

0

23:28Z

89

0.82

4.0

0.4

0.67

1

1

372.8

150

372.8

0

0

0

0

23:29Z

88

0.79

4.0

0.43

0.65

1

1

348.7

150

348.7

0

0

0

0

23:31Z

86

0.77

4.0

0.54

0.61

1

1

275.3

150

275.3

0

0

0

0

23:32Z

85

0.74

4.0

0.63

0.6

1

1

248.2

150

248.2

0

0

0

0

23:33Z

90

0.72

4.0

0.63

0.68

1

1

239.8

150

239.8

0

0

0

0

23:34Z

91

0.69

4.0

3.53

0.71

1

1

59.9

150

59.9

0

0

0

0

23:35Z

87

0.67

4.0

0.78

0.64

1

1

234.1

150

234.1

0

0

0

0

23:37Z

85

0.65

4.0

1.08

0.61

1

1

229.7

150

229.7

0

0

0

0

23:38Z

86

0.63

4.0

0.79

0.62

1

1

241.5

150

241.5

0

0

0

0

23:38Z

88

0.62

4.0

1.2

0.65

1

1

192.1

150

192.1

0

0

0

0

Waiting 60sec for next loop at 01:40; Variant "noChange"

174

175 The line at the bottom of the table says the time (hh:mm) when the next bg result and loop
 176 decision are expected. Also it shows the (yourChange).vdf file investigated

177

178 In sequence of the time(Z!), the table consists default of the columns with info on (1) bg (2)
 179 IOB and eff.iobTH (3-6) the adaptation factors on ISF suggested by the 4 autoISF
 180 categories (7) resulting ISF that was used, profile_ISF, and emul ("what-if" ISF) (8)
 181 insulinRequired calculated by your running loop, and "what-if" result (9) same for resulting
 182 SMB

183

184 The columns marked “**orig**”.in the table of results shows the ISF actually used to
185 determine insulinReq and SMB size in the actual run.

186

187 In the columns marked “**emul**” you find the calculated emulation results that is calculated
188 every 5 minutes

189 • If a (yourChange).vdf was clicked (in step 2), see also bottom right of the table), the emul
190 columns show the result, how the investigated changes *would have* changed SMB sizes (strictly
191 always looking at the one decision, in that line of the table).

192 • If only a noChange.vdf was run, the emul columns contain the same results as the orig.
193 columns.

194

195 **Caution when interpreting the values in the acce, bg, pp and dura ISF columns (3-6)**

196

197 The factors given there are always those for the emul run.

198 • So when using (yourChange).vdf, you see in your phone's table exclusively the ISFs that
199 *would result-if* (yourChange) were made.

200 We get to “what-if” projects later, in [section 11.4](#).

201 • Only in the noChange scenario, the values there would be the *orig ones, corresponding to*
202 *what could be seen also in the SMB tab* at the times

203 You easy recognize whether you run the noChange: it says so at the bottom of the table.
204

205 orig.ISF is called **sens** in the SMB tab, a couple of lines below “end autoISF”. It is the ISF
206 that replaces, for that one decision you are looking at, the profile ISF (called profile.sens in
207 the SMB tab),

208 We get back to this topic at end of [section 11.4.2](#) “Understanding how the ISF is emulated by
209 (yourChange) and how SMB or TBR would differ”

210 “

211

212

213 11.2.3 Analyzing the results

UTC	bg	IOB	iobTH	eff. ISF	acce ISF	bg ISF	pp ISF	dura ISF	-----ISFs-----	insulin Req	---SMB---
time	bg								orig prof emul	orig emul	orig emu
13:53:42Z	116	2.19	6.0	1.54	1.01	1	1	26.1	40 26.1	0 0	0 0
13:58:42Z	120	2.0	6.0	2.08	1.01	1.12	1.04	19.2	40 19.2	0.16 0.16	0.1 0.1
14:03:57Z	117	1.97	6.0	1.37	1.01	1	1.06	28.5	39 28.5	0 0	0 0
14:08:43Z	113	1.81	6.0	1.18	1.0	1	1.08	33	39 33.0	0 0	0 0
14:13:41Z	115	1.64	6.0	1.11	1.01	1	1.1	35.2	39 35.2	0 0	0 0
14:18:42Z	125	1.47	6.0	1.86	1.02	1.3	1	20.9	39 20.9	2.3 2.3	1.7 1.7
14:23:46Z	131	3.19	6.0	1.44	1.03	1.18	1	27	39 27.0	0.26 0.26	0.1 0.1

214 In above example (table), all given SMBs were driven by bgAccel_ISF, when glucose rose.
 215 The biggest SMB in the time we are looking at (actual local time = Z + 2 hours, so we are
 216 looking at late small rises, like 3 hours after lunch) was 1.7 U = 0.74 SMB delivery rate *
 217 2.3 U insulinRequ.

218 The insulinRequ. Is proportionally amplified by the effective ISF, called sens in your SMB
 219 tab, or ISF“orig” in this table. The amplification of 39 profile_ISF / 20.9 effective_ISF = 1.86
 220 happens to be the dominating bgAccel_ISF amplification factor

221
 222 Note **not** always just the biggest ISF factor “wins”. Consult flowchart in [LINK](#), and
 223 occasionally read the real-time explanations in your SMB tab as to which other factors are
 224 contributing to the amplification resulting from profile_ISF to effective_ISF (“orig”, or sens).
 225 For instance, the Activity Monitor, or a set %temp. profile, or TT, could contribute, or also
 226 the question whether glucose already decreases.

227
 228 Depending on your selected “safety” settings, you might occasionally bump into
 229 restrictions. **Tuning up factors that make the system bounce into restrictions is a**
 230 **completely useless, and potentially even dangerous, exercise!**

231
 232 Therefore, **for your initial set-up** of parameters ([section 2](#) and [section 4](#) of the FCL e-
 233 book) it is **advisable to not do this analysis** on the smartphone, but **on the PC**, where
 234 you can inspect the complete info on each loop decision (see [section 10](#))..

235
 236 If for some reason you cannot do this on the PC, you must frequently take screenshots (in
 237 very many decisive 5 minute segments), and analyze **more complete data**, than the table
 238 on our smartphone could offer, **in the SMB tab**)

239

240 In your initial set-up of your FCL, you probably must “loop” a couple of times back into
241 [section 2](#) to adjust the safety settings made for max. SMB sizes.
242
243

11.3 Options available in Trio or iAPS


The emulator will not work in the iOS world.

However, iAPS / Open iAPS offer in their autoISF variants also access to a tabular as well as graphical representation of **autoISF contributions to resulting SMB sizes**

Note that on the iPhone, so far, the what-if emulation and speech synthesis (see next [section 11.4](#) for AAPS) are currently not possible. Also not the PC analysis as in section 10.

(insert example of such table and graph)

... is this an example from i-APS ?



How to get table and graphics (double click on statistics button ?)

Screenshots and description must be provided by iAPS / Trio user please

270

271 11.4 Real-time checking a „what-if“ question using speech synthesis

272

273 The emulator on your **Android** smartphone can help clarify "what if..." you implemented a
274 considered change.

275

276 11.4.1 Put a (yourChange).vdf into the phone's AAPS logfiles folder

277

278 In running the emulator on the phone, you can define in the .vdf file of the emulator, which
279 setting(s) you would like to be differently aggressive than in the active AAPS.

280

281 How to write .vdf files was already explained in the section "Emulator on PC". See in end of
282 [section 10.2.1](#). You could also pull a vdf file example from the developer's Github that you
283 could customize further: Access see [section 3.8 /5](#).

284



285 Produce or retrieve **(yourChange).vdf**. on PC from one of your studies files.

286

- 287 • C:\....\Documents\ Looping\ Logfiles_Emulator\ Emulator_Studies\ Study_n

288 PC: your path to the Emulator files may differ -

PC: C/ Documents / Logfiles_Emulator > Emulator_Studies > Studi



289

290 ... and put a **copy** of that .vdf into the **smartphone folder with the AAPS logs** :

- 291 • Phone: Internal memory/AAPS/logs/info.nightscout.androidaps

292

293

294 Switch between python scripts running at the same time

295

296 In case you have more than one (yourChange).vdf to investigate, you look at the data of
297 your currently running loop (last 15*5 minutes) by just switching between the related vdf
298 files used for emulation.

299 Details see [section 3.8, 9](#)) – or try to go direct via [https://github.com/ga-zelle/APS-what-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)
300 [if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)
301 [phone.pdf](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf) and there p.5, under above sub-headline “.Stop the emulator, or switch...”

302

303 [Stop the emulator](#)

304 (see in paper as above)

305

306 [11.4.2 Inspect emulated results](#)

307

308 Now, whenever you run QPython 3L emulation (following the steps as described in [section 11.2.2](#))
309 you additionally get the emul. results filled in in the resulting table, showing in which time points
310 your change would lead to adapted insulinReq, and SMB size data.

311

23:57

←

No. 1

NEW

CTRL

UTC	eff.	acce	bg	pp	dura	-----ISFs-----			insulin Req		---SMB---			
time	bg	IOB	iobTH	ISF	ISF	ISF	ISF	orig	prof	emul	orig	emul	orig	emu
20:41:16Z	112	1.93	6.0	1	1.0	1	1	40.9	41	40.9	0	0	0	
20:46:16Z	111	1.77	6.0	1.14	1.0	1	1.03	36	41	36.0	0	0	0	

etc.

21:11:17Z	108	1.06	6.0	1.12	1.0	1	1.08	38.3	43	38.3	0	0	0	
21:16:18Z	109	0.94	6.0	1.1	1.0	1.03	1.07	39.1	43	39.1	0.1	0.1	0	
21:21:20Z	109	0.87	6.36	0.95	1.0	1	1.08	41.9	43	39.5	-0.1	-0.08	0	
21:26:20Z	108	0.76	6.36	1	1.0	1	1	43	43	40.6	0	0	0	

312

313 etc

314 In the marked 21;21 Z line, dura_ISF was the dominant factor. The (yourChange).vdf would apply a
315 factor of 1.2 and lead to further strengthening the ISF: profile 43 -> orig.(noChange) 41.9 -> emul
316 (yourChange) 39.5

317 In this case, late after a meal, and bg=109 mg/dl, the loop saw in the orig. (noChange) case 0.1 U
318 insulin too much; and as the (yourChange) emul case asks for typically more insulin (all weights in
319 the (yourChange).vdf are > 1), now only 0.08 U are seen as too much (a 20% difference).

320

321 The table on your phone is too reduced to show each emulated ISF component. If you need to see
322 more details on how (yourChange).vdf would make a change in a point of time that you like to
323 analyze deeper: Note down interesting Z times, and later look it up in more detail in the Emulator
324 on the PC:

325

below are "orig."ISFs:

UTC		eff.	acce	bg	pp	dura	-----ISFs-----			insulin Req		---SMB---		
time	bg	IOB	iobTH	ISF	ISF	ISF	ISF	orig	prof	emul	orig	emul	orig	emu
21:21:20Z	109	0.87	6.36	0.95	1.0	1	1.08	41.9	43	39.5	-0.1	-0.08	0	

Understanding ISF orig.(the sensitivity used to determine insulinRequ)

To fully understand how acce, bg, pp and dura_ISF determined the orig.ISF, we do an **emulator** run **on the PC**, and look up the **.txt** result:

@ 21:21 Z / Script Debug -----

Activity monitor disabled inactivity detection: sleeping hours; Autosens ratio: 1;

Basal unchanged: 0.45; ISF unchanged: 43 CR: 9

start autoISF 3.0.1

Loop allows APS power level; SMB enabled due to enableSMB_always

acce_ISF adaptation is 0.95

bg_ISF adaptation is 1

pp_ISF adaptation is 1

dura_ISF adaptation is 1.08 because ISF 43 did not do it for 30 m

strongest autoISF factor 1.08 weakened to 1.03 as bg decelerates already

final ISF factor is 1.03

end autoISF

profile.sens: 43 sens: 41.9

Note: 43 / 1.026 = 41.9

Understanding how the ISF is emulated by (yourChange) and how SMB or TBR would differ

While the result for ISF_emul (and for the SMB size) is given in the table on your phone, the details behind it, and also for finer effects in %TBR, can only be inspected from the .csv table from the (yourChange),vdf **emulator** run on the **PC**:

360 The following table is an extraction of the most relevant data from the (yourChange).csv:

	B	C	E	F	L	Q	R	S	AA	AB	AC	AD	AE	AF	AG	AH	AI	AP	AQ
1			bg	bg		final	dura		acce	bg	pp	delta	dura	final					
2	UTC	AAPS	accel	brake		ISF	min-	dura	ISF	ISF	ISF	ISF	ISF	ISF	ISF	ISF	ISF	TBR	TBR
3	time	time			iob	orig	utes	avg.	emul	emul	emul	emul	emul	emul	orig	prof	emul	orig	emul
29	21:16:18	23:16	109		0,94	1,1	25	108,2	1,12	1	1,03	1	1,08	1,12	39,1	43	38,4	0,65	0,65
30	21:21:20	23:21	109	109	0,87	1,03	30	108,3	0,9	1	1	1	1,1	1,04	41,9	43	41	0	0,46
31	21:26:20	23:26	108		0,76	1	5	108,5	1	1	1	1	1	1	43	43	40,6	0	0
32	C=B+n/24 @ n=+2 h time diff.										If there were bigger orig->emul effects, we would show SMB columns, too								

361

362 It shows that the final (noChange) ISF factor of 1.03 (box Q30) becomes 1.04 (box AF30) in the
 363 (yourChange) emulated case. Given that, so late after dinner (box C30), we sit near the 90 mg/dl
 364 target (box E30), the difference the changes in the emulated case are minor. Still zero SMB, but
 365 the minor changes reflect in elevating the temp. basal rate for fine adjustment (box AP30 -> AQ30).

366

367 Note that the orig. values for the four autoISF components (acce, bg, pp, dura "ISF orig")
 368 are **not** contained in this csv table (only, in Q30, the final amplification factor for final_ISF in
 369 the noChange scenario).

370 However, you can fetch the autoISF details for 21:20 Z in the actual "noChange" run:

- 371 • from the txt result file (as shown above); or
- 372 • you also see them on the phone:

373

orig	orig	orig	orig	ISF
acce	bg	pp	dura	-----
ISF	ISF	ISF	ISF	orig
0.95	1.0	1	1.08	41.9

374

375 ...while the emul_ISFs come from .csv results @ (yourChange).vdf run on the PC:

AA	AB	AC	AE	AI
acce	bg	pp	dura	
ISF	ISF	ISF	ISF	ISF
emul	emul	emul	emul	emul
1,12	1	1,03	1,08	38,4
0,9	1	1	1,1	41,2

376

377

378 Apology: The above example was not well chosen to see relevant effects. The author is struggling to
 379 put this chapter first time together, quickly for the V.3.0.1 launch, and just picked from his phone what
 380 was available at the moment, I might patch it over with a better example in a later update, or I (or
 381 maybe you?) provide an adjunct case study.

382

383

384

385 Principal limitation of any of your emulations

386

387 Note that always the first biggest change regarding insulinRequ and SMB size in emul vs. orig. is
388 the most relevant. This is because:

- 389 ○ Doing that change would change, by the same amount, the iob basis for the next
390 following loop decision
- 391 ○ Doing that change would certainly change the course of the bg curve about half an hour
392 afterwards (to be precise: for the duration of DIA, for the extra insulin), so then it is
393 anybody's guess, for instance how difficult of a job presents itself to the loop to "attack"
394 with duraISF.

395

396 This is one of the reasons why setting up your FCL ([section 4](#)) should be a iterative process,
397 seeking solutions (with mainly bgAccel_ISF) for the first rise (and a range of different meals) first,
398 AND NOT concurrently already "tuning" the dura_ISF.

399 Both, bg peak height, and pattern of insulin activity from the present iob (that is only sluggishly
400 further adjustable driven by dura_ISF and evtl. by bg_ISF) depend on how the first bg rise stage was
401 managed by autoISF (your bgAccel_ISF_weight setting, and others).

402

403 Another conclusion you might draw, is to limit use of the emulator *on the PC* largely to analyzing
404 the (in FCL extremely important) first rise, to seek bgAccel_ and pp_ISF_weights.

405

406 Then to test these, and many other like ideas for changed settings "in real life" emulated *on the*
407 *smartphone*: A synthetic voice will announce to you how (if) each actual decision would differ. See
408 next section!

409

410

411 11.4.3 Real time speech synthesized treatment suggestions

412 Unfortunately this great feature is **not available on i-Phone**. Look in (updates of) [section 11.3](#) for
413 eventual alternatives.

414

415 At time points when the (yourChange) setting would result in smaller or greater difference in SMB
416 insulin delivery (compared to the real "noChange" run), you can get a real-time notification **via**
417 **speech synthesis**, and you can assess the situation in real-time yourself.

418 If for instance a suggested extra, or bigger, SMB makes sense, **you can add this portion**
419 **manually*and observe**, for this meal, whether this bolus was OK and you should switch in
420 direction of the different setting you were investigating (which would automatically give you that
421 extra in the future).

422 Likewise, you might choose **not** to intervene, but regret it an hour or so later, seeing the further
423 development **without** implementing the supposed improvement.

**In Full Closed Loop, you don't need any buttons at the bottom of the AAPS main screen. But for such test phases it is practical to re-install the insulin button at the bottom of the AAPS main screen (Preferences/Overview/Buttons/Insulin -> ON).*

After a couple of days, you will get a feel for whether you want to incorporate your investigated change (or a gradual step towards it) into your active AAPS settings.

Warning: Your settings must always work for a variety of meals. Do not put too much effort into optimizing one situation! (See [case study 8.2](#)).


Activating and silencing emulator suggestions

For silencing the suggestions from voice synthesis you have the following options:

- (1) Change **principal settings** what shall be announced (e.g. only if bigger SMB size is suggested, or also warnings about carbs eventually needed?), and in **which hours of day**, to make any announcements via speech synthesis .

These are set in the (1 or) 5minute_emulator_std.config file:

Go on **PC** into the config. file (see illustration). Open it with Notepad++ and edit the hours there for when you would want (no) announcements regarding: extra carb need (line 1), extra bolus need (line 2), or less bolus (line 3):



Save the changes, and copy the file also into your **phone** at Internal memory/AAPS/logs /info.nightscout.androidaps (see [11.1.3](#)) over the 5minute_emulator_std.config

- (2) Turn off **phone volume** (silence media + switch on do not disturb)

451 Of course, this also shuts off many other potential alerts that you might not want to shut off.

452

453 (3) **“Kill” (and later resume) the “what-if” emulation.** This could be done by de-
454 selecting the (yourChange).vdf in step . However, this stops (or interrupts, until you
455 re-start) the entire emulation and you will have no tabular data later for the silenced time.

456

457 (4) **De-activate Qpython 3L app** (temporarily?): Press on Qpython 3L app icon, force
458 close or remove necessary permission – re-activate (or need to re-start emulation then ??)
459 when you want to hear again .. However, this stops (or interrupts, until you re-
460 start) the entire emulation and you will have no tabular data later for the silenced time

461

462 (5) [□] Probably better than the 2 afore mentioned options would be to **run (for the**
463 **intended silencing period) exclusively the noChange.vdf,**
464 Then you will not get any speech outputs (because you are NOT investigating a what-if
465 question, in that case). But you get all data (the un-interrupted noChange actual run) and
466 later on the PC still can investigate any “what-if” scenarios

467 How to change the .vdf reference during a run, see [section 3.8, 9](#) – or go direct to the
468 instructions in the Github repo at: [https://github.com/ga-zelle/APS-what-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)
469 [if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf)
470 [phone.pdf](https://github.com/ga-zelle/APS-what-if/blob/A3.2.0.4_ai3.0.1/Documentation%20in%20English/How-to-run-the-emulator-on-the-phone.pdf) and there p.5, under above sub-headline “.Stop the emulator, or switch...”

471

472 Note that this area (silencing) has not been researched much, and good answers are likely
473 lifestyle and phone specific.

474

475

476 Please share your experiences with the emulator in Discord / Full-Closed-Looping / HOW
477 TO / _emulate-aaps, at: <https://discord.gg/n3tD5eXExC>

478

479

480