

Please note that with autoISF 3.0 you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



9.1 How to get back into Hybrid Closed Loop

You can go back to bolussing for meals and making carb inputs again at any time just by

- going into AAPS Preferences/OpenAPS SMB/autoISF and switch „Enable ISF adaptation by glucose behaviour“ OFF.

You might need to re-install your insulin button via AAPS Preferences/Overview/Buttons

- pressing on the **violet Full Closed Loop circle** and select the green Hybrid Closed Loop circle (easier, if that user interface element is already included).

This will automatically bring back your buttons “Insulin, Calculator...” you always had at the bottom of your AAPS HCL main screen

Be aware that now it is again up to you to bolus for meals

Note that even if you had an established Hybrid Closed Loop *with autoISF*, the HCL you switch back to would be normal OpenAPS SMB (without autoISF ISF modulation).

That is the safest bet in absence of sufficient data whether tuning the `_weights` in autoISF would have to differ between FCL and HCL mode. According to a n=1 observation they might work with the same setting (see https://github.com/ga-zelle/autoISF/blob/A3.2.0.2_ai3.0/To%20prebolus%20or%20not%20to%20prebolus.pdf)

but more data are needed, notably also the variance with size of pre-bolus, and kind of meal.

It can also be wise, especially in your initial months, to do FCL only for certain meal time slots, and use an **Automation** that shuts down

- either `bgAccel_ISF` utilization,
- or all of autoISF `_ISF` modulations („Enable ISF adaptation by glucose behaviour“) for the rest of the 24 hour period.

In that case, the loop button will automatically adjust its color violet <-> green to show which state your loop operates under (if that user interface is already included)..

For instance, it is perfectly possible, without any extra steps involved, to do FCL only for dinners, while breakfast and lunch are done in hybrid closed loop as you are used to.

9.2 Are the pre-conditions for FCL still given?

- * Is the basic profile still correct?
- * Has the CGM quality deteriorated?
- * etc (see [section 1](#) pre-requisites)

9.3 Glucose goes too high

- Meals are not recognized asap
 - Check regarding Bluetooth (in)stability
 - Experiment with an aperitif, soup a couple of minutes before meal start
- First SMB(s) seem a bit delayed
 - Check whether SMB got blocked by the 30% rule (refer to [section 1.3](#)), and what the underlying cause may be
 - Check whether an odd bg target or TT (maybe in context with an Automation) interfered
 - Check pump connection (BT, and physical)
 - Check stability of regular CGM values (notably: was phone in proximity at meal start?)
- SMBs are too weak
 - Check acceleration detection (e.g. CGM, BT or smoothing related)
 - Check (real-time) in SMB tab what ("safety"?) setting limits allowed SMB size
 - Check whether your autoISF_max is set too low in AAPS preferences
 - Check (real-time) in SMB tab whether bgAccel_ISF_weight or pp_ISF_weight should be set higher
 - Check potential interference from sensitivity modulations (e.g. exercise mode active)

64 ○ Check in preferences whether, after a FCL pause, you “forgot” to re-activate “Enable
65 ISF adaptation by glucose behavior” , or whether an Automation could have temp.
66 deactivated it, and hence you had received SMBs only from the basicoref(1)
67 SMB+UAM algo, without boost by autoISF

68 ○ In case you are (e.g. due to an on-coming infection) more insulin resistant currently,
69 consider – as in Hybrid Closed Loop – a temp. profile switch to >>100% and do not
70 shift other settings (that would be wrong when your general sensitivity bounced
71 back to your normal)

- 72 ● Default iobTH (or: modulated iobTH) might be (or: go) too low, and therefore cuts SMBs too
73 early

74 Check whether your FCL *really* operates with the iobTH you think it uses. Observe the
75 Caution notes (e.g. in [section 5.1.4](#)) about needing re-sets to default, after an Automation
76 had lowered the effective iobTH.

- 77 ● An important observation by pilot users was, that how your glucose and iob curves
78 approach meal start matters a lot regarding how you peak from carbs:
79 Going down (e.g. towards a set EatingSoonTT), building some iob, and curving already
80 towards strong positive acceleration seems very helpful to keep peaks low.

- 81 ● When all trouble shooting ideas are exhausted, notably all measures to pull powerful SMBs
82 earlier are exhausted, and the tail of insulin activity already pushes you close to a hypo:
83 Then you simply must **find your personal balance** between

84 ○ Accepting *sometimes* trending higher than you would like to, for not going low.

85 ○ Change diet (probably to something with lower amounts of carbs, and higher
86 amount of protein and fibre).

- 87 ● A pretty benign way to improve bg development after meals could be to just take a
88 walk (notably if starting before/ when glucose seems “stuck” high).

- 89 ● Some users resort to using a small pre-bolus in their “FCL” (maybe just for to-them-
90 known troublesome types of meals). However, this interferes with how glucose curve
91 and hence detection of rises and triggered SMBs behave. It is therefore not easy to
92 implement with convincing overall benefit.(See discussion in [section](#)).

93

94

95

96

97 9.4 Glucose goes too low

98

99 • Meals are falsely recognized

100 ○ Check whether this occurs outside of usual meal times and can be averted by e.g.
101 setting an odd profile target for that time of day.

102 ○ Check whether your bg_Accel_ISF driven first SMBs are too big and lead to too
103 much iob when the bg rise turns out just a temporary bump.

104 ○ Try to exclude the problematic situation via an Automation that, for a couple of
105 minutes only, sets an odd TT and thus prevents a SMB.

106 ○ To prevent snacks from triggering SMBs as for a meal, use the FCL cockpit for an
107 appropriate temp. setting (low iobTH ;or odd TT for SMBs off; or pre-set “snk”
108 button).

109 • SMBs deliver overall too much insulin

110 ○ Check whether you operate with a too high iobTH.

111 If an Automation might have elevated the effective iobTH, check whether your FCL
112 *really* operates with the default iobTH again afterwards. Observe the Caution notes
113 (e.g. in [section 5.1.4](#)) about needing re-sets.

114 ○ Check (real-time) in SMB tab whether SMB range extention or autoISF_MAX should
115 be set smaller

116 ○ Check (real-time) in SMB tab which of the autoISF ..._weight should be dialled in
117 smaller . Often it will be a too strong dura_ISF. However, that one inheritantly gets
118 stronger with higher and longer lasting highs. Therefore the best remedy is to first
119 try to be more aggressive before, in the glucose rise phase and limit height and
120 duration of the high, then tweak the dura_ISF_weight (downwards).

121 ○ SMB delivery ratio probably can be set smaller. Note in this case, it works across
122 the bord for all SMBs (all time slots),

123 ○ In case you are (e.g. due to a preceding sports day) more insulin sensitive currently,
124 consider – as in Hybrid Closed Loop – a temp. profile switch to >>100% and do not
125 shift other settings (that would be wrong when your general sensitivity bounced
126 back to your normal)

127 • Problems with insulin “tail” after meals

128 ○ See 2 bullet points higher up: dura_ISF tuned too strong?

- You may need to take a snack (seeing hypo prediction) or glucose tablets (if already in hypo zone). But note that the grams of carbs required that the loop might tell you at some point are very likely exaggerated as the loop has no info *) on your carb intake (while you may be able to guess how much more, incl. from fats and proteins) is still waiting to be absorbed. (*)The loop makes assumptions based on past minutes carb deviations, see reference given in [section 4.5](#))
- A valuable information would be whether the problem originates mostly in the bg rise phase already. Then setting a lower iobTH might be an easy remedy.
- If the need for additional carbs happens frequently, note down how many grams were needed (not counting what you eventually took too much which required extra insulin again). Then use your profile IC value to estimate how much insulin less the SMBs should have delivered, and go with this info into your tuning (regarding the % profile in the Automations, or maybe also your set iobTH). This may relate to the SMBs given when glucose was high, or also extend regarding the SMBs during the glucose rise.
- When all trouble shooting ideas are exhausted, you simply must **find your personal balance** between
 - Accepting *sometimes* trending too low and needing a <10 g snack. (There are worse things in life, and if you are weight conscious, eat these grams less at the meal itself. Over time you should learn at which type of meals – probably those *low* in fibre, fat, and protein - this scenario arises)
 - Accept on average a bit higher glucose peaks, for not going low.
 - Change diet (probably to something with lower amounts of carbs, and higher amount of protein and fibre).

154

155 9.5 Glucose goes too high and too low

156

157 Frequent roller coasters point to serious problems with your set-up.

- 158 • **Try not do master everything at once.** Is your eating and general life style just so
159 extreme, and your expectation into the capabilities of the system too high / your
160 understanding of its limits too low? Then you could consider giving the FCL another
161 serious try for periods that are less varied (e.g. just one of your daily meal times, only
162 weekends, not days with Xtreme sport ...). Make it work there, then gradually expand.
163 Sections 5. and 6. describe a mind boggling number of ways to deal with special
164 “disturbances”. Just occasionally try one that interests you. Stay connected with others that
165 are in the same boat.

166 Even when basic pre-conditions (see 9.2) seemed given, and you “tried already everything”:

- 167 • Was your autoISF FCL built based on true and experimentally proven ISFs?
- 168 • Did you follow the sequence of tuning steps (sections 2, then 4 ; bgAccel_ISF-weight first?)
- 169 • How often did you consult SMB tab or emulator, to gain an understanding what is/was
170 happening?

171 With the multitude of inter-acting parameters and settings (that already after a short time would be
172 burdened with counter-balancing errors) it is extremely difficult to untangle and correct this. Best
173 idea then might be

- 174 • an entire new start. (There is emphasis all over this paper that with autoISF FCL you are in
175 a development project. So, taking some steps back and starting over from there should be
176 acceptable.)
- 177 • going back into your prior Hybrid Closed Loop (potentially with some benefits from autoISF
178 also there)
- 179 • or resorting to another method as e.g. mentioned in section 13
- 180 • or switching to a simpler loop system as fully supported by your doctor, and wait for
181 improvements the industry will provide over time, too.

182

183

184 9.6 Staying out of Trouble...

185

186 In closing this chapter we like to remind everyone that interfering with a closed loop should be kept
187 at a minimum.

188 Also, „optimizing“ parameters for just *one* specific meal type or other experience is a **flawed „fine-**
189 **tuning“ concept** when you are in FCL, and can easy backfire (see [case report 8.2](#) as an
190 example). What we want is settings that get us „good-enough“ through (nearly) all scenarios in our
191 personal everyday lifes.

192

193 This is a good time to throw in a reminder, where the real world of T1Ds stands, and the consensus
194 in the medical community, regarding desirable %TIR (or HbA1c, as only a minority would have TIR
195 data).

196

197 Weigh for yourself what it is that you try to achieve. For instance, 80%TIR was reached in a
198 study even with a much simpler FCL (AAPS with Automations, Fiasp, no autoISF), and no meal
199 announcements whatsoever.

200

- 201 • **Do the basics right, keep it simple.** Resist the temptation to embark always on the latest
202 craze without knowing how extra features might topple your carefully set balance.
- 203 • Learn to use FCL in some times, and not in others that you believe may be too challenging,
204 or you have already a bad experience with (and no time, interest, skill, to resolve it for now.
205 That is fine, too.).
- 206 • Stay in touch with the community of developers and other users
- 207 • **Relax and enjoy** as/when/while good-enough. „Just eat!“

208