

1. Pre-requisites for Full Closed Loop

V 2.4

Please note that with autoISF 3.0 you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



- 1.1 Well tuned hybrid closed loop
- 1.2 Fast insulin
- 1.3 Reliable insulin delivery from pump and cannula
- 1.4 Excellent CGM
- 1.5 Meal-related limitations
- 1.6 Lifestyle-related limitations
- 1.7 Time required for setting-up

Available related case studies:

- Case study 1.1: Occlusion
- Case study 1.2: Comparing insulins for FCL
- Case study 1.3: Jumpy CGM
- Case study 1.4: Lost pump connection

1.1 Well-tuned hybrid closed loop

It is advisable to first establish a well-tuned hybrid closed loop before considering the transition to FCL. There are two important reasons for this:

- The UAM full closed loop requires a highly personalized (individual) tuning of settings, so the loop will give insulin mimicking YOUR successful hybrid closed loop mode.
- The UAM full closed loop comes with new parameters to be set and tuned. It would be problematic to set and tune several new parameters before the basics were tuned “right”. Errors could easily be balanced with counter-errors. This can work in single scenarios, but would create a highly unstable system, hard to re-calibrate better later.

1.2 Fast insulin (Lyumjev, Fiasp, Apidra?)

If the user does not bolus for meals, clearly a very fast insulin is needed so, upon realization of a starting meal-related glucose rise, the loop has any chance to eventually keep glucose in range (by common definition, under 180 mg/dl (10 mmol/l))

A modelling study (details see

<https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#fast-insulin-lyumjev-fiasp>)

can show in quantitative terms that **faster insulins**

- will result in significantly **lower** glucose **peaks** than slower insulins
- **tolerate** a couple of minutes **delayed** first meal bolus while not incurring unacceptable height of peaks
- minimize the effect on glucose peak from **different** carb load (**meal sizes**).

39 In conclusion, do not attempt FCL with other insulin than Lyumjev® or Fiasp®, unless, maybe, if
40 you are on a very moderate to low carb diet. (According to [case study 1.2](#), Apidra® might work,
41 too, but Humalog® would not work well).

42 43 1.3 Reliable insulin delivery from the used pump/cannula/insulin system

44
45 **Good Tolerance of Lyumjev (or Fiasp): Occlusions threaten the function of the full closed loop.**

46
47 It is very important to have an eye on the time a **cannula (or pod)** is in use (many find **48 hrs** to be
48 the **limit**), and whether hard-to-explain glucose rises happen at ever increasing „fake“ iob (even
49 before a 48 hr routine replacement). (See [case study 1.1](#): You easily lose 25% TIR that day)

50 It is absolutely contra-indicated to attempt FCL coming from leaking pods and associated erratic
51 sensitivity swings that may or may not have been somewhat controlled and tolerable by
52 dynamicISF or other measures when you were Hybrid Closed Looping,

53

54 **Stable pump connection**

55 In FCL you absolutely rely on your pump delivering,
56 without any further delay, the much needed insulin,
57 after any meal start.

58 Hence it is absolutely essential not to “miss” any
59 problems from a lost Bluetooth connection.

60 An Automation similar to the one pictured here →
61 should help recognizing eventual problems.

62

63 See also [case study 1.4](#)

64

65 1.4 Excellent CGM

66


67 You do not give a meal size-related bolus any longer. That leaves all insulination jobs to the
68 algorithm! As glucose values are the very basis for this, please **inform yourself well about** how
69 **your CGM** 1) principally performs 2) whether and how this may depend on data flow and
70 intermediate apps you use 3) specifically, how and where any smoothing is done, and what this
71 might imply for the ISF boosting method you will be using See for instance here:

72 <https://androidaps.readthedocs.io/en/latest/Usage/Smoothing-Blood-Glucose-Data-in-xDrip.html>

73

74 Around meals, a stable Bluetooth connectivity is absolutely essential, too, so CGM, loop, and
75 pump can do their job without losing more valuable time (see [case study 1.4](#)).

76



77 Then, but even more importantly in *all other* day and night *times*, the CGM should not produce any
78 artefacts (jumpy values; see [case study 1.3](#)) that the loop could **misinterpret** as sign of a starting
79 meal. Note that also calibrations could produce jumps.

80

81 autoISF has also a couple of in-built checks on the quality of the recent CGM values. Hence, a
82 CGM with more scatter will make the loop lose more time, and lead to higher peaks and
83 lower %TIR.

84 So, if you are unhappy with a slow reaction of your loop it could be because the loop is unhappy
85 with your CGM. Consult the detail info given (at the time) in your SMB tab, or look it up later in the
86 logfiles (using the emulator eventually).

87

88 The best proven way to stay out of trouble currently is to use Dexcom G5 or **G6**, and to ensure via
89 **overlapping** right and left arm sensor and transmitter utilization for always good quality values that
90 can be used by the Full Closed Loop.

91 Other ways (making use also of day-1 sensor values, G7, Libre2/3...) are possible, but come with
92 a lot of monitoring effort (best via watch) and occasional time-outs for the FCL.

93

94 One safety feature in autoISF is a **blockage of SMB delivery whenever delta bg** (within the last
95 two 5 minute values) is **higher than 30% of that bg**. *So from 74 mg/dl, a jump to 97 (+23) or more*
96 *would not receive SMB "response", or from 100 mg/dl to 131 mg/dl (+31) would neither.*

97

98 Check in your (HCL or FCL) data whether at meals or sweet drinks with rapid absorbing
99 carbs you could run into the problem that jumps are "too high" and much needed insulin will
100 be blocked (only come via very much smaller portions: *e.g. 400%TBR @ 0.6 U/h => 0.2 U in*
101 *5 minutes, instead of one ~3 U SMB. The difference of 2.8 U missed translates @ ISF~ 40*
102 *mg/dl/U into up to + 112 mg/dl higher bg peak! It will not become quite that bad, because*
103 *the loop will catch up to the insulinRequired with it's next couple of decisions).*

104

105 Instead searching in old data, you can also just have an eye on instances where you think a
106 first SMB was due, but blocked. Confirm that by looking in the SMB tab, and think about a
107 solution that would not require changing the 30% safety limit in the code. *For instance, not*
108 *drinking so much juice rapidly around meal start could be a likely "behavioral" correction to*
109 *get rid of the problem.*

110

111 This blockage (no SMBs) would likely last only 5 minutes (and go probably unnoticed - **in a**
112 **suggested improved User Interface it would be indicated for about 5 minutes via a dotted violet full**
113 **closed loop logo/button on the AAPS home screen**). However, not only would you lose 5 valueable
114 minutes to get your iob substantially elevated; all following deltas are likely much smaller, and, as a

115 consequence, you will miss some of the boost sought from bgAccel_ISF if the >30% delta was in
116 fact (largely) due to carb absorption.

117

118 This example also underscores that the CGM in use cannot be allowed random scatter that leaves
119 no reasonable room for safe detection of (smaller and) bigger “truly carb related” deltas

120

121 If or when (like: first half day of a new sensor) you are not sure about sufficient CGM
122 performance you might develop for yourself an Automation with User action ticked (along
123 the lines as used for other purposes in [section 5.2.2.3](#)). It would “ask you” before giving a
124 SMB whether you really want it delivered. That way you can a) have a look on your glucose
125 curve b) on the delta and acc data underneath the TT field of your AAPS home screen c)
126 think about what sense a SMB now makes with respect to your last meal, and the carbs to
127 be still absorbed. Ultimately, you could also d) consult some of the detailed info given
128 (every 5 minutes) in your SMB tab.

129 Such User action Automations need not be active at all times, but if you have it for
130 your first half day of a new G6 sensor for instance, you could activate that
131 Automation from your list of Automations; after the values have settled in, you can
132 deactivate (“shelve”) it again.

133 For a brief period, and if you are tech savvy, another way to deal with uncertainty about
134 CGM would be to employ the emulator method as presented in [section 11](#): Run a “too
135 mildly” tuned FCL, and in parallel run a “what-if” with your more aggressive settings that you
136 really would like to use once you are certain about your CGM.

137 However, I found it easiest to lay a solid groundwork by using 1 Anubis, and 2 overlapping
138 G6 to get rid of most problems that I saw (and keep seeing, on the worse sensor of the two
139 running for some days often in parallel) in my data.

140

141 With a sensible iobTH defined, and your standard alarms for going towards a hypo not silenced,
142 the worst consequence from any automatically “over-treated” glucose jump should be that you
143 need an unplanned snack for the balance of “missing” carbs.

144

145 Also FS Libre 3 is useable in the autoISF FCL context. Observe info in the general section about
146 autoISF (<https://github.com/ga-zelle/autoISF>) regarding established ways to use with autoISF.

147 As of Dec.2023, there is still development work going around what to best make of the 1
148 minute values.

149

150

151 1.5 Meal-related limitations

152

153 Setting up a full closed loop is relatively easy for people whose diet does not consist **mainly** of
154 components with rapid high effect on blood glucose (more see
155 <https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html#meal-related-limitations>)

156

157 Meals do not have to be low on carb (provided you use a fast insulin for your FCL)

158 Fat or protein rich diets, or slow digestion/gastroparesis, make things easier rather than harder for
159 the full closed loop because late carbs nicely cover for inevitable “tails” of late action from SMBs
160 needed around peak time.

161

162 Erratic consumption of snacks with fast resorbing carbs can be a problem.

163 In autoISF you can reduce this problem to some extent via one or two keystrokes from your
164 AAPS home screen. While certainly being a deviation from the FCL idea(l), this would be
165 one of the exceptional situations where you better do a quick “nudging” step from your “FCL
166 cockpit”. Details see in [section 5.2.1](#) and [5.3.3.1 \(4\)](#)

167

168

169 1.6 Lifestyle-related limitations

170

171 Technically stable system

172

173 Full closed looping requires a 24/7 technically stable system, especially regarding

- 174
- reliable CGM signals
 - Bluetooth stability with the pump (see [case study 1.4](#))
 - keeping your phone in sufficient proximity at all times
 - avoiding (or at least early recognition of) occlusion.
- 175
- 176
- 177

178 This requires a habit (or, unlikely, permanent attention to details) like keeping all components well
179 charged and in close proximity; making cannula (or pod) changes always early enough to lower the
180 risk of occlusion (see [case study 1.1](#)); having always potentially needed parts with you.

181 **Depending on your system, your experience with it, but also on your acceptance and general**
182 **lifestyle, these aspects may or may not limit you.**

183

184

185

186

187 Preparing for exercise

188

189 To prepare for exercise (sports, heavy work), the normal protocol with a pump or hybrid closed loop
190 is to take actions that reduce insulin on board prior to exercise

191 With your full closed loop, the algorithm is tuned to detect meals and to give you insulin to counter
192 glucose rises automatically. Setting a high temp. target and lower %profile right away (effective al-
193 ready around meal start) could be a problem.

194 Unusual activity levels therefore likely require **disciplined preparation** (especially **if you want to**
195 **keep the need to snack during sports low**)

196 In autoISF you can reduce this problem to some extent via two or three keystrokes on your
197 AAPS home screen. While certainly being a deviation from the FCL idea(l), this would be
198 one of the exceptional situations where you better “flick a lever” from your “FCL cockpit” to
199

200 Extra hurdles to establish FCL for kids

201

202 To establish and maintain a FCL for kids brings about some extra challenges if:

- 203 • Lyumjev is not available or well tolerated
- 204 • Hourly basal rate is very low, providing a poor basis for big SMBs
- 205 • Diet is rich in sweet components. With the typical low blood volume of a small body, strong
206 tendency towards very high bg spikes!
- 207 • Going through marked changes of insulin sensitivity or of circadian pattern makes it difficult
208 to keep the FCL appropriately tuned.

209 This problem is about the same in Hybrid Closed Looping. However, now you might
210 expect miracles from the FCL. This is not going to happen. You still should try to set
211 appropriate temp. changed profiles, that serve also as a basis for your autoISF FCL.

- 212 • Discipline is poor regarding keeping Bluetooth connectivity and infusion sites perfectly run-
213 ning
- 214 • Between kid and supervising parent it must be guaranteed, especially in the initial weeks,
215 that an eye is kept on whether the FCL is working about as to be expected.

216 More see [section 7](#).

217

218

219

220

221

222 1.7 Time required for setting-up

223

224 Lastly, before enjoying a functioning full closed loop you need to have a period of a some weeks
225 with some free time and „free head“ for set-up –. Can you get, in the time you are willing to invest,
226 to a result that you consider good-enough is really the question. Depending on your „habits“, and
227 which – if any - compromises (like doing cannula/pod changes more often, never starting meals
228 when bg sits high ...) are you willing to make (and everyday able to stick to), for the ease of not
229 having to deal with assessing meals and bolussing for them?

230

231 While setting up your personal FCL using autoISF is a substantial project, there is no need
232 to implement it fully in one step. There is nothing wrong to go in your well running Hybrid
233 Closed Loop mostly, while switching to FCL only for dinners, for instance, or only for
234 weekend lunches, as a start. Once you found feasible settings, you can expand to other
235 meal times, and lastly towards figuring out your best strategies for challenges outside of
236 meal windows, as we shall discuss in [sections 5. and 6.](#)

237

238 There are alternatives to using autoISF for FCL, as well. See [section 13.](#) for more info.

239

240 Notably [13.1](#) FCL using AAPS Master and Automations could be a much easier and more
241 error-tolerant way of stepping into FCL. In a clinical study with 16 participants about 80%
242 TIR was achieved without much tuning effort (source: see in [section 13.1](#)).

243

244 To close the circle to where we had started ([section 1.1](#)): A very time consuming pre-requisite might
245 actually be to first sort out your Hybrid Closed Loop, so your profile parameters are set „right“, and
246 your “old” data really can serve as a blueprint for what, now, you would like *your loop* to do in FCL
247 mode

248 Note that if you had used dynamic parameters or special Automations („loops inside the loop“) this
249 might have balanced some principal errors, but leaves you now without a good starting point as
250 you must get rid of these over-patches (see also warnings at start of [section 4](#))..

251

252 You will see also success stories of loopers who just jump into using more powerful tools, in
253 kind of a trial and error mode, and frequently add the latest add-on, or self-constructed
254 patch (often in form of an Automation) to counter-balance problems.

255 So, yes, you can also continue in that spirit. Resulting solutions may be good-enough. But
256 they tend to be unstable and not well-understood. That is a poor basis for managing arising
257 problems (-> fine tuning), and for adjusting to special situations (-> which setting to
258 temporarily change). But it certainly is an alternative avenue for the impatient, less
259 analytically, and more adventurous inclined.

260 In any case, PLEASE always observe the safety settings/instructions coming with the DIY
261 dev- variant of software you select.
262
263 One key safety measure every AAPS user going towards FCL should have in place is to set an **iob**
264 **threshold** (iobTH; size a bit below what you used as a bolus for bigger meals in HCL) above which
265 no more SMBs can be given by your FCL.
266 This is an integrated feature of autoISF, from 3.0 version onwards (see [section 2.4](#)).