

6. Temporary modulation for exercise and lighter (in-)activity V 3.7

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them.

This is not a medical product, refer to disclaimer in [section 0](#)



6.1 Dynamic iobTH and sensitivity ratio

6.1.1 Manual (direct) iobTH modulation

6.1.2 Automations for iobTH modulation

6.1.3 Dynamic iobTH

6.2 Temp. % profile switch

6.3 DIY cockpit based on User action Automations

6.4 Improved FCL cockpit

6.4.1 Manual (direct) iobTH modulation

6.4.2 pre-sets for 4 kinds of exercise

6.4.3 optional meal pre-sets

6.4.4 optional hypo management pre-sets

6.5 Mastering the exercise after meal challenge

6.5.1 Manual mode

6.5.2 DIY cockpit button for User action Automation

6.5.3 Laissez-faire mode

6.5.4 Using pre-sets in improved FCL cockpit

6.6 Activity monitor based on step-counter

[Available related case studies:](#)

Case study 6.2: Biking day with hi carb lunch;
..... DIY cockpit

Skip what is in green writing:

= Drafted fragments or not implemented ideas.

Preliminary remarks

This section is **no easy read** because it attempts to describe *all options* to deal with *various types* of exercise.

- Fortunately, you might need *none* of them,
- you can set *any* of them up at your leisure, later, one at a time, for any of your occasional or regular exercise events,
- then pick one or two of the described options, how to go about it.

- As long you were not able yet to define better ways, you should always be able to manage bg dropping during sports with **extra snacks** (keep those at hand).
- Also, and especially in your first weeks of FCL (with focus on meal management), you can always **exit into Hybrid** Closed Loop or even Open Loop, and manage certain exercise the way you always did before.

Staying in contact with the related discord/github community should help greatly to find suitable ways to manage *your* type(s) of exercise.

- Please report *your* experience by supplying a case study.

Looking at case studies that relate to *your* kinds of exercise might be easier to digest than working your way through *all the options laid out* in this section.

For detail insights into how to manage **various kinds of exercise**, listen-in this reference by looping pioneer and sportswoman Dana Lewis: https://bit.ly/DC1_631 (starts around minute 05:30), or read up in: <https://diyps.org/.../how-to-exercise-when-exercise-is...>

6.1 Dynamic iobTH and sensitivity ratio in „exercise mode“

iobTH is a iob threshold you can set, above which AAPS will no longer deliver additional SMBs.

(This overrides the SMB management via even/odd bg target differentiation).

(Regarding by how much “*the last SMB*” may shoot over iobTH, see [section 2.4](#)).

For exercise, we like to limit how high iob can go, therefore automatic “dynamic” reduction of your set iobTH (= iobMAX x iobTH%) is a benefit, notably as you can individually tune it.

In autoISF 3.0 and later, a setting for iobTH is made in AAPS preferences, defined there as fraction (e.g. 0.6) of your set maxIOB:

/OpenAPS_SMB/autoISF_settings/Full_Loop_settings: iob_threshold_percent,

=> default iobTH = iobMAX x **iob_threshold_percent**

63 So, while iobTH could also be modulated via iobMAX, we mostly adapt the
64 iob_threshold_percent, to do that.

65

66 In the following, 3 principal avenues to temporary adjusting iobTH to your exercise
67 requirements are described: **Manual** intervention ([6.1.1](#)), making use of individually defined
68 Automations ([6.1.2](#)), and relying on the **automatic** dynamic adjustments coming with
69 autoISF ([6.1.3](#)).

70 The author experimented with all of them, but rarely needs manual intervention
71 because the *automatic* “**dynamic**” adjustments work pretty well, after some individual
72 tuning (see e.g. case study 6.2).

73 In any case, it is good to educate yourself about manual tweaking options, should the
74 need arise.

75

76 6.1.1 Manual (direct) iobTH modulation

77

78 „Manual“ routes to directly change iobTH would be

- 79 • changing the setting for the new parameter „iob_threshold_percent „
 - 80 • or changing the setting for iobMAX
- 81 in /Preferences.

82 This is not a preferred route for temporary adjustment, because it is not easy accessible with
83 just a button stroke, and it would not automatically revert to your prior setting, after use.

84 A future improved FCL cockpit (-> [section 5.3](#)) might give direct access to

- 85 • **override iobTH temporarily, at any point of time.**

86 A bridging solution that can achieve nearly the same is: to construct your own “DIY
87 cockpit” button to change iobTH% from the AAPS main screen, see next section.

88

89 6.1.2 Automations for temporary iobTH modulation

90

91 You can define Automations that set a different iobTH% **under pre-defined conditions**

In a variation of this idea (if your Automation has the User Action box ticked), you get a grey button into your AAPS home screen, from which you can activate that changed iobTH manually (“DIY cockpit”, as was already presented in [section 5.2.2.3](#)).

Note that this is the iobTH you tell the loop to use **in place of** the previously set iobTH:


- it will still be modulated further if %profile and TT are set (see below)
- it will overwrite the iobTH% you had set in /preferences!

Caution: A different iobTH% or bgAccel_ISF_weight can *not* be set *temporarily* in Automations (i.e. a duration cannot be attached). You **must** define a suitable **additional Automation that** must be active *in tandem*, and **restores the prior set iobTH%** or bgAccel-ISF_weight later **again**. Else, once your Automation sets in, it will *forever* shift this important parameter setting!

If for instance you have several Automations that, in combination with a set elevated TT also set a lower iobTH: Don’t be fooled, the duration only applies to the TT. You need an extra Automation for all of them.

Example: My Automation that restores my prior set profile iobTH%:
I picked out the *highest* of the *lowered* iobTH values that *any of my* Automations can set (45 percent was the highest “of the low ones” in my case), and then I can automatically restore to my *prior* 60% via this one:

Caution: Watch out for another potential stumbling block, because many Automations only work under the condition that no TT is already running.



As temp. changing iobTH is quite tricky to automate, it is the author’s preferred route to only *indirectly* modify it – see next section.

6.1.3 Dynamic iobTH: iobTH modulation via setting a temp. glucose target (TT)

In AAPS/Preferences, set “**High TT raises sensitivity** = TRUE”. Then, setting an **elevated** temporary glucose target (TT), decreases iobTH by the same factor as it increases sensitivity (as it “softens” ISF). Both measures decrease insulin that the loop will give.

123 Likewise, In AAPS/Preferences, set “**Low TT lowers sensitivity** = TRUE”. Then, setting a
124 **low** temporary glucose target (e.g. a EatingSoonTT of 74 mg/dl), elevates iobTH by the
125 same factor as it also sharpens (lowers) the ISF. The loop will give more insulin.

126 6.1.3.1 How does automatic sensitivity and iobTH adaptation work in the exercise mode?

127 **When. additionally. the exercise button is ON** (lit yellow), **iobTH gets reduced**
128 **particularly strong, and ISF is particularly weakened** (as desired for exercise).

129 That effect is the stronger (**ISF gets the weaker, iobTH the lower**), **the lower you set the**
130 **half-basal exercise target** for your exercise mode in AAPS/preferences/OpenAPS SMB:

131 The following table shows, for a profile target of 100 mg/dl, how the set ...

- 132 • half_basal_exercise_target you set in AAPS/preferences/OpenAPS SMB...
133 Choose a low number if you later want a high dynamic range of sensitivity modulation
134 Lower half-basal exercise target = lesser insulin delivered
- 135 • ...and current exercise TT (that you set on the day you do the respective exercise,
136 with an eye on how you wish sensitivity auto-adjusted)...
137 Higher TT = lesser insulin delivered

138 ... determine the effective sensitivity ratio:

Half basal ex.target	180	150	120
TT	sens.ratio	sens.ratio	sens.ratio
100 = profile target	1	1	1
120	0,8	0,71	0,5
140	0,67	0,56	0,33
160	0,57	0,45	0,25
180	0,50	0,38	0,20

139 The exact calculation for any combination of profile target, set TT, and half-
140 basal_exercise_target is given in the autoISF Quick guide (see [section 3.3](#)).

141 Note that:

- 142 • **temp. basal = profile basal * sens.ratio**

143 *Example: At a half-basal_exercise_target of 120, setting a TT of 120 gives only half (0.5) of profile*
144 *basal (hence the name of the parameter)*

- 145 • **temp.ISF = profile ISF / sens.ratio**
- 146 • **temp.iobTH = set iobTH * sens.ratio**


147 Whereas in “vanilla” AAPS the sens ratio is simply coming from you (manually), or from
 148 Autosens (automatically) setting a temporary profile sensitivity other than 100% (and in the
 149 special case of dynamicISF with additional effects on ISF), here, in autoISF, we have strong,
 150 non-linear, and user scaleable effects on the sens.ratio.

151 6.1.3.2 How you recognize the real-time iobTH, and “aggressiveness” status of your FCL loop in 152 general

153 Rather than bothering with the math, you can just look into your **SMB tab** where your
 154 selected temporary settings put your iobTH, and the modified ISF (called **sens**):

155 The valid “effective iobTH” can also be
 156 seen in the SMB tab, see example (for a
 157 **TT=144** and exercise button clicked);

Same, with (via top left button in AAPS home screen) **additionally 70% profile** applied:

158 

159 Result

160 Script debug : d:

161 Activity monitor disabled:

162 tempTarget

163 Sensitivity ratio set to 0.53 based on temp target of 144

164 Adjusting basal from 0.75 to 0.398;

165 ISF from 39 to 73.6 =39 / 0.53

166 CR: 7

167 -----

168 start autoISF 3.0.1

169 -----

170 User setting iobTH=60% modulated to 32% or 3.18U

171 due to profile % and/or exercise/activity mode

172 SMB enabled; current target 144 is even number

173 Loop allows medium power

174 acce_ISF adaptation is 0.71

175 bg_ISF adaptation is 0.65

176 final ISF factor is 0.34 including exercise mode impact

177 -----


178 end autoISF

179 -----

180 currenttemp: 0 lastTempAge: 0

181 m tempModulus: 21 m

182 profile.sens: 39 sens: 113.2 =39 / 0.34

158 

159 Result

160 Script debug : d:

161 Activity monitor disabled:

162 tempTarget

163 Sensitivity ratio set to 0.53 based on temp target of 144

164 Adjusting basal from 0.52 to 0.276;

165 ISF from 57 to 107.8 =57 / 0.53

166 CR: 10.714285714285

167 -----

168 start autoISF 3.0.1

169 -----

170 User setting iobTH=60% modulated to 22% or 2.23U

171 due to profile % and/or exercise/activity mode

172 SMB enabled; current target 144 is even number

173 Loop allows medium power

174 Parabolic fit extrapolates a maximum of 107 in about 4 minutes

175 acce_ISF adaptation is 0.95

176 bg_ISF adaptation is 0.65

177 final ISF factor is 0.34 including exercise mode impact

178 -----

179 end autoISF

180 -----

181 profile.sens 39 sens: 162 = 39 / (0.34 * 70%)

173 These examples show that, **by just setting an exercise TT and a typical exercise**
174 **profile%** (two super easy “interventions” via the top buttons on our AAPS main screen,
175 turning yellow in response as an easy “reminder” we are in a special mode), **the iobTH will**
176 **be automatically very sharply reduced** (in our example, to about half just by the TT; and
177 further to about 1/3, by the % setting)

178 From the figures given (left side), the user’s iobTH calculates to 6.0 U (= (60%/32%)*3,18U)
179 which is 60% of iobMAX of 10.0 U. So, normally, autoISF FCL could give SMBs up to
180 anywhere between 6.0 U and 7.2 U (=6.0 +20%; see [section 2.4](#) at step 2.4: at bg>100, iob
181 can run max. 20% over with “last” SMB). For doing exercise, this window shrinks now to 3.18
182 – 3.82 U (left) or, even to 2.23 – 2.77 U (right).

183 In conclusion, these easy-to-make settings (TT, %profile) automatically provide the same
184 thing like would have been done in Hybrid Closed Loop, where a meal bolus of about 7 U
185 would get cut down to 4 U or even to 3 U, depending on type of exercise.

186 If you concurrently use QPython 3L and the emulator on your Android phone (see [section 11](#))
187 you need not look into the SMB tab, but could see more details (~ for the last hour, plus all
188 contributing ...ISF_categories from autoISF), in tabular form, on your phone.

189 For i-Phone autoISF users, double clicking on “Statistics” also provides similar information
190 (see [section 11.3](#)).

191

192 6.1.3.3 Customization and-tuning

193 Try to determine good settings for the kinds of exercise that you engage in:

194 Set your **half-basal exercise target** in /preferences that suits all of them...

- 195 • ... in tandem with reasonable TTs that you like to use later, for each of your
196 intended specific exercises
- 197 • Iterate through this a couple of times (whenever you happen to do *that*
198 exercise).

199 ⇒ **Remember** (“code” for yourself), **which TT stands for which exercise type, so**
200 **that just by setting *that* TT everything (ISF, iobTH) *automatically* will provide**
201 **the lower loop aggressiveness that you need *for that specific type of***
202 **exercise.**

203

204

205 When setting a TT please watch out for unintended implications and side-effects:

206 (1) Setting a TT often **shuts out other** Automations.

207 Therefore, choose the **duration** wisely (and also the **sequence**, in which all your
208 Automations are listed).

209 (2) You always must consciously decide whether you set an **even or an odd** numbered bg
210 **target** (TT or profile target). (This is assuming you use, as you should, the even/odd bg
211 target differentiation for SMB on/off).

212 • Pick **odd**, if you do not want SMBs during exercise. (Despite your softened ISF, and
213 lowered iobTH, SMBs still might „attack“ a sports snack too strongly).

214 ○ However, odd cannot be set too early, when your meal digestion still
215 requires SMBs.

216 ○ Likewise, you might want the option for a few automatically delivered SMBs
217 against unforeseen spikes (e.g. from excitement) also later.

218 In that case, an **Automation** might **sneak in a desired SMB** or two via
219 switching from odd to even, just for a couple of minutes, and under a well
220 thought-out set of conditions (that you must find in **your** data patterns,
221 when you do that kind of exercise that you try to find good settings for),

222 However, you are probably out of luck because an already set odd (or
223 any) TT would preclude such Automation from kicking in. Then you
224 need to **develop additional** ideas, another detour, like to first define an
225 **Automation that briefly shuts your odd TT down**.

226 ○ So, defining everything so you really can be happy with oddTT being your
227 primary way is a quite tricky project you should not under-estimate.

228 • Working with an **even** TT can be preferable, notably of course if your exercise is one that
229 can get you totally excited, with glucose spikes.

230 ○ While this mode generally does allow SMBs, the loop softens the ISF (by the
231 sens.factor like in the table given above), and will temp. shut SMBs down, when
232 **iobTH** (which also got lowered by the sens. factor) is exceeded.

233

234 Whether odd or even TT is better depends on the kinds of exercise you are doing, and
235 probably depends on the protein and fat load of your meal and snacks, as well.

236 (3) **Timing** can be **critical** as to when you do your exercise announcement, especially
237 relative to a preceding hi-carb meal. Then you want the reduced iobTH in place latest after
238 you received the first SMB. See [section 6.4](#) and [case study 6.2](#)

239 (4) Once you are familiar with the **dynamic range of your iobTH**,

- 240 • after you made your settings, notably set your half-basal exercise target
- 241 • knowing the range of TTs and %profile adaptations that you intend to use
- 242 before/during/after your types of exercise

243 please confirm or re-consider *your iobTH_percent setting in /Preferences*, [section 2.4](#).

244 (5) You always can **look the effective iobTH up in the SMB tab** (see screenshots given 3
245 pages earlier).

246 *In future releases you might see the valid iobTH that your loop is working with also in*
247 *your AAPS home screen, next to the current iob status.*

248

249 6.1.4 Tweaking iobTH

250

251 You can use any of the above discussed methods, or also the one that now follows in [section](#)
252 [6.2](#), to *further tweak* iobTH temporarily, should you see a need.

253

254 Also outside of exercise, setting an **even elevated TT plus** pressing the **exercise button**
255 gives easy access **to significantly reduce aggressiveness of your autoISF loop** via a
256 resulting lowered iobTH and, concurrently, elevated effective ISF.

257 This could be used for instance for 45-60 minutes **at low/medium carb snacks**, as an
258 alternative to shutting SMBs **entirely** off via an **odd** TT.

259

260 *When exercise follows a meal*, it might be smart to use the just discussed tweaking methods
261 right after you felt the sting from the first big Lyumjev SMB.

262 However, we will look at smarter and safer ways for this “exercise after meal” scenario in
263 [section 6.5](#) and in [case study 6.2](#)

264

265

266 6.2 Temporary % profile switch

267

268 A complementary measure you can take from the AAPS home screen is to set a **reduced**
269 **temp.% profile** sensitivity.

270 This setting would **multiply** with the results in above table and further **reduce basal and**
271 **iobTH**_(whenever exercise button AND profile button both are yellow).

272 An example was already given with the 2nd screenshot, 4 pages earlier

273

274 Note that temp. reduction of basal will proportionally also **reduce the max. allowed size of**
275 **SMBs** (which is two hours worth of basal x SMB_range_extention, see [section 2.1](#))

276 This is desirable at exercise, albeit seldom needed, because milder ISFs and radically lowered
277 iobTH anyways will reduce iob delivered (see at lines 141-146).

278

279 The **time windows** for doing a profile switch *can differ* from the time window (duration) of
280 your TT-related exercise settings. **Using all available tools then allows a nearly surgical**
281 **approach to what you want to achieve for and during your favorite exercise(s).**

282 • Often the %profile modulation is used for several hours if not days to accommodate “long
283 waved” sensitivity changes (See e.g. in [case study 6.2](#)).

284 • Instead, or even additionally, the percentage might be modified for just a couple of
285 minutes, or for one special snack or meal duration, to “nudge” the proportionally
286 modulated aggressiveness of the FCL (see [section 5.2.3](#)).

287

288 You can prepare yourself for anything you see coming up, or potentially coming up, in your
289 daily life, so, from the comfort of your cockpit ([section 6.3](#); [section 6.5.2](#)) you get ready for it
290 within just a second or two, doing a few „clicks“.

291

292

293 6.3 Managing exercise via Cockpit inputs

294

295 You may want to skip reading the *green texts* which are about suggested further features in
296 future software updates

297

298 6.3.1 Basic Settings for Exercise

299

300 Coming from FCL with no TT set (the top fields TT and exercise are grey), you best prepare
301 for an intended exercise by **pressing the TT field** of your AAPS main screen (your looping
302 cockpit; presented in [section 5.2](#)).

303 There, you can **freely select** TT and duration.

304 *Alternatively, you can press on **one of 4 offered exercise presets**. (Note: This, and many
305 other – in this green color - described cockpit features are currently not developed)*

306

307 In case you want maximal “softening” effect on basal, ISF, and iobTH, activate also the
308 **exercise button** in the top middle of your AAPS main screen. It then should turn yellow, to
309 indicate you are in the dynamic exercise mode (YY \underline{Y} or G \underline{Y} Y).

310

311 6.3.2 „Dynamic“ exercise mode off = traditional AAPS exercise mode (YGY)

312

313 When the dynamic exercise mode is off (YGY), you still have the instruments for *exercise*
314 *management just as you always had it in the past* = a combination of manually softened
315 aggressiveness via setting a temp. %profile change, and orienting corrections towards an
316 elevated TT.

317

318 By selecting an odd numbered TT you now have the *additional option* to shut SMBs
319 temporarily off, too.

320 Skip the next 7 pages, and **continue with [section 6.5](#)**

321 *The following **green texts** (and embedded screenshots) describe features **suggested for
further development** which the author meanwhile sees not really needed.*

322 If improved cockpit is launched, the top part of the dialogue box looks about like this when
 323 the exercise field is grey:

activity	TT	dura	%profile	iobTH	bgAccel	effect.sens.
n/a	171	180	70	3.0 U	0.24	70 %

Mode set to run for 134 more minutes

324
 325 % profile can be changed:

- 326 • either here => neighboring %profile button turns yellow too (with the % info on it); or
- 327 • under the %profile button; or
- 328 • it had already been changed using AAPS / Action / Profile switch

329 In all 3 cases, you see the number < 100 or >100 in the middle of above table, on a yellow
 330 colored field, too.

331 In this “YGY” mode, the % temp. set profile is the applied “effective sensitivity” (% ratio)

332

333 TT and duration can be entered or changed (= traditional mode to set exercise targets).

334

335 If there is a desire to try, for the **remaining duration**, a different iobTH or bgAccel_ISF-
 336 weight, this can be overridden in the table; field turns yellow, and the algorithm uses temp.
 337 iobTH and/or temp bgAccel_ISF_weight as modified in the exercise button (and reports this
 338 also in the SMB tab).

339 In the dialogue box pictured above, 70% profile was set for 3 hours, and the default
 340 iobTH of 60% * 10 U was cut by 50% down to 3.0 U.

341 The remaining duration shows below the table (in the example: 134 minutes and counting
 342 down).

343

344 The effective iobTH is given in the SMB tab. In later versions, the effective iobTH should
 345 show also in the AAPS home screen, next to the actual iob (e.g. „1.2 U < 3.0 U“)

346 TT and % profile will also show on the yellow labels of the neighboring %profile (left top of
 347 AAPS home screen) and TT (right side), respectively.

348 The middle (exercise) field remains grey because the automatic sensitivity tuning (that would
349 use TT and half-basal exercise target) are off.

350

351 6.3.3 Dynamic exercise mode ON (GY or YY)

352

353 By pressing the yellow exercise button on the AAPS home screen, *you have the*
354 *option* to switch the **dynamic exercise mode ON**, in which case the middle
355 field/exercise button of your AAPS main screen will go from grey to yellow.

356

357 In a version update you could do your setting for the upcoming exercise under the **dialogue**
358 **box of the TT button**

359 Then, when you look into the exercise button in the middle of your FCL cockpit the dynamic
360 exercise mode will automatically be „ON“, and all entries made:

Dynamic exercise mode **ON**

activity	TT	dur	%sens	iobTH	bgAccel	effect.sens.
mtb	171	180	100	4.0-11	0.16	67 %

• Mode set to run for 134 more minutes

• Mode starting after meal when iob > iobTH

361

362 The data for the kind of exercise (*here mtb; could also be n/a or ?*) are coming from prior
363 selections made in the dialogue box of the neighboring TT field. There, as well as in this
364 window here, the resulting iobTH and bgAccel_ISF_weight are shown. Also the overall
365 aggressiveness (% overall insulin sensitivity factor) is calculated.

366 The **middle field** of the table in this dialogue box, **% profile** either picks up the % set under
367 the %profile button, or an input can be made here, in the exercise button domain, which will:

368 • turn the neighboring %profile button on yellow and show that inputted % on it, too

369 • be multiplied with the result from the exercise mode settings per se, and change the
370 effective sensitivity %, accordingly.

371 So, if this middle field of above table (dialogue box of exercise button) contains a figure other
372 than 100, the input field becomes yellow, and you are operating with a combination of
373 traditional PLUS new exercise mode (with all three top buttons of your FCL cockpit yellow).
374 This maximally will soften aggressiveness, for which you get an idea by the last calculated
375 figure.

376 The mode is either running already (*for another 134 of the total 180 minute in the picture*) as
377 also the label on the neighboring yellow TT field will show *171 (134, and counting down)*,

378 Or (*see at the red dot in picture above*), it is scheduled to run, after insulination for a started
379 meal surpasses iobTH (*as in table*).

380 Note that, when the TT expires or is changed, your overriding input (if you made any)
381 is automatically erased, forgotten.

382

383 6.3.4 Dynamic exercise mode ON plus %profile change (YYY)

384

385 The **middle field** of the table in the dynamic exercise mode dialogue box (see above), %
386 **profile**“ either picks up the % set under the %profile button, or an input can be made here, in
387 the exercise button domain, which will:

388 • turn the neighboring %profile button on yellow and show that inputted % on it, too

389 • be multiplied with the result from the exercise mode settings per se, and change
390 the % overall, accordingly.

391 So, if this middle field of above table (dialogue box of exercise button) contains a figure other
392 than 100, input field becomes yellow, and you are operating with a combination of traditional
393 plus new exercise mode (with **all three top buttons** of your FCL cockpit **yellow**). This
394 **maximally will soften aggressiveness**, for which you get an idea by the last calculated
395 figure.

396 Maximal softening happens only **when >100% profile**. A <100% profile is not used for
397 exercise typically, and would counter-act the softening from the elevated TT.

398 It is advisable to find good settings primarily *within* the dynamic exercise mode. Use profile
399 switches only complementary as needed. notably to provide for other, „longer waved“, health
400 or hormonal situations.

401

402 Note that you can make use of the **exercise button** at any time for easy tweaking..

403 • yellow -> off/grey = increase

404 • grey -> on/yellow = decrease (only possible if a TT> 100 mg/dl is set)

405 ..the aggressiveness, without changing (and changing back) settings.

406

407 6.4 Option to pre-set for 4 kinds of exercise or meals (for 1 button operation)

408

409 You may want to skip the entire section 6.4 because *green texts* describe features *suggested*
410 *for further development* which the author meanwhile sees not really needed

411

412 6.4.1 iob_threshold_percent

413

414 In AAPS preferences/OpenAPS SMB/autoISF settings / Full Loop Settings, the default
415 iob_threshold_percent used for the normal meal spectrum is defined.

416 In an updated later autoISF version you might be able to differentiate there for up to 4 meal
417 clusters (see next section)

418

419 6.4.2 Pre-settings for (up to) 4 kinds of exercise:

420

421 In AAPS preferences/OpenAPS SMB/autoISF settings / Full Loop Settings: follows next input
422 fields for pre-settings you can define for (up to) 4 kinds of exercise:

423 The following table gives an example of settings you may find well-suited for 4 of your
424 favourite exercises

#1-4	give name (max 3 characters)	duration for TT (min)	TT (AC) (mg/dl)	% profile	iobTH	bgAcce:weight	Approx % ins reduct.
1	wlk	60	111	100			
2	grd	120	131	90			
3	bik	300	151	90			
4	mtb	180	171	70			

425 Input fields (during tuning phase to determine good settings) are only the columns 2-5.

426 The last 3 columns will be calculated from TT and %profile inputs, using also the half-basal
427 exercise target and the default weight setting. In this setting.

428 The last is only an approximation to get a feel for a reasonable setting of the other
429 parameters.

430 Here in preferences they should never be overridden, but TT or % profile should be adjusted
431 to reach desired result when tuning for FCL.

432 Likewise, you find tables to make pre-settings for meals and for hypo treatments:

433

434 6.4.3 Pre-settings for (up to) 4 kinds of meals:

435

436 In AAPS preferences/OpenAPS SMB/autoISF settings / **Full Loop Settings**: follows next:

437 Input fields for pre-settings you can define for (up to) 4 kinds of meals. For instance:

TT# 1-4	give name (3 letters)	TT (Eating Soon) (mg/dl)e	Duration for TT (min)	iobTH (0---130% and < iobMAX)	bgAcce factor 200...0%	
1	hiC	72	120	110	110	
2	loC	74	180	67	67	
3	piz	76	300	100	100	
4	snk	78	60	100	50	

438

439 Input fields (during tuning phase to determine good settings) are all columns

440 Difference in TT is fairly unimportant (unless you do not give a name and memorize the set
441 TT number instead, for which meal type it codes.

442 Logic why not having a % profile column here: %profile switch should be set extra,
443 potentially for another time period (e.g. „reserved“ for periods of exercise, or for entire
444 days of altered insulin sensitivity, for instance due to illness, fasting, extensive sports
445 week.)

446

447 6.4.4 Pre-settings for (up to) 4 kinds of Hypo treatment:

448

449 In AAPS preferences/OpenAPS SMB/autoISF settings / **Full Loop Settings**: follows next:

450 Input fields for pre-settings you can define for (up to) 4 kinds of HYPO treatment. Example:

TT (ES) (mg/dl)	give name (3 letters)	TT (AC) (mg/dl)	Duration for TT (AC) (min)	bgTH (mg/dl)	
1	Hy1	131	55	none	
2	Hy2	131	55	200	

451

452 Input fields (during tuning phase to determine good settings) are all columns, 2-5.

453 Choosing an odd-numbered TT is recommended as it can shut-out SMBs (with the
454 appropriate setting in preferences/Open APS SMB/autoISF settings/smb_delivery
455 settings/“enable alternative activation...“.

456

457 Those of us who tend to over-treat hypos may prefer to set Hy2 (unless for night snacks->
458 Hy1): Reverting to standard loop aggressiveness with SMBs after/if a certain bg level
459 („threshold“, similar to our iobTH for meals) is surpassed, and we want our loop to react
460 again with SMBs before the set duration expires.

461 6.5 Mastering Exercise after a Meal

462

463 In Hybrid Closed Loop, we gave less insulin at meals (a reduced bolus) before exercise.

464 Since we now get our meal insulin automatically from the loop, we would have to at least
465 somehow tell it that exercise follows this time.

466 Simply setting an exercise profile *before* the meal would make our full closed loop too weak
467 in the "treatment" of the first glucose rise. **What we want is, to get our** (already, compared
468 to HCL, delayed) **meal insulin delivered as fast as possible by SMBs. It just should be**
469 **capped at the desired iob reduction.**

470

471 6.5.1 Manual mode requires 2 user interventions

472

473 What we can do, is (1) **reduce** the **iobTH** (via the `_%` setting, *e.g. by one third*).

474 • *In the example we were using, this would mean to reduce by 2 U to $iobTH^* = 4U$.*

475 • Do that estimate for your data, and think back how you did bolus reduction in hybrid
476 closed loop before same exercise.

477 • Likewise, you can use your profile ISF, *e.g. 30 mg/dl/U* and „translate“ by how much
478 ($2U * 30 \text{ mg/dl/U} = 60 \text{ mg/dl}$) this „pulls you away from going into a hypo“.

479 • Using your IC (*e.g. 8g/U*) you can also translate the iobTH reduction (2 U) into a
480 „snack equivalent“ ($2U * 8 \text{ g/U} = 16 \text{ g}$) that you „replace“ by thinking ahead and
481 „budgeting“ for some exercise with your iobTH modulation.

482 In this senario, our loop delivers SMB insulin as fast as always, only that when the last SMB
483 has passed the iobTH, the loop only has elevated %TBR to work with, meaning it cannot
484 raise iob by much any longer. This provides an elevated glucose level on which we enter
485 exercise, and saves us hypo danger or snack need (as calculated in above examples).

486

487 After this reduced iobTH is reached (or up to 30% exceeded by the last SMB, up to 20% @
488 even $TT > 100 \text{ mg/dl}$), step (2) must follow = an increased exercise **bg target** is set (see
489 [section 6.2](#)).

490

491 The problem with this approach is that it requires **two** user interventions, first **setting the**
492 **lower iobTH%**, and later (**and this *in a time-critical manner***, after iobTH is exceeded), to
493 **input an exercise TT**, or to activate a related setting.

494 To eliminate this problem, the following refined solutions are suggested:

495

496 6.5.2 DIY cockpit: Using pre-set meal / exercise settings from a User action Automation

497

498 The „DIY cockpit“ user interface allows a *one-step* setting for meal + exercise that can be
499 selected in time-uncritical fashion, any time before the meal starts (or even shortly thereafter,
500 within the “grace period” after which we hope to already see the first SMB triggered).

501 A detailed example is given in [case study 6.2](#):

502 A sequence of 3 Automations must be set up, of which only the first one must be manually
503 triggered, in just one time-uncritical key stroke from the AAPS home screen.

504 The others are activated automatically, when the respective Conditions are met.

505

506 Automation #1 provides, for a meal that precedes exercise, the full loop aggressiveness, but
507 makes sure that this aggressiveness stops immediately after a (reduced) iobTH is exceeded.
508 The reduced iobTH ensures that not too much insulin is on board for exercise after the meal.
509 Also it provides an elevated bg level at (re-)start of exercise.

510

511 In this Automation, the box “User action” should be permanently ticked. This will
512 automatically provide a **grey button on the bottom of the AAPS home screen** (“DIY
513 cockpit”) that can be freely named (= headline of Automation #1).

514 For exercise that is not done frequently, I choose to get rid of that cockpit button by disabling
515 the Automation fully, in my list of Automations... until the evening before e.g. a bike tour, when
516 I will want to have my cockpit give me the optional button again.(See [case study 6.2](#))

517

518 As soon as the (reduced) iobTH is exceeded, two things need to be provided :

519 (1) a milder running FCL (reduced exercise %profile, after the meal rise had been
520 managed based on 100% profile boosted further by bgAccel_ISF driven full loop
521 aggressiveness) => Automation #2 sets e.g. 70% profile and ends TT

522 (2) setting an exercise TT (not possible with Automation #2. But *after* it terminated the
523 TT, an Automation #3 can immediately follow, and set the desired exercise TT=125
524 (which implies the exercise mode

525 Note that Automations 2 and 3 are fully automatic, no User action is involved. See [case](#)
526 [study 6.2](#) for an example

527

528 Should, during the exercise, a need arise to modulate the loop aggressiveness (jobTH,
529 effective ISF), this can be done within 1-2 seconds, also right from the AAPS home screen
530 („FCL cockpit“), by setting a higher or lower temp. %profile, and/or by setting a higher or
531 lower temp. exerciseTT.

532 To make the loop temporarily act a bit more aggressive, switching the exercise button OFF
533 (from yellow to grey) could also be considered

534

535 **Defining User action - Automations to build your FCL cockpit**

536

537 If you want to develop *your* **DIY User Interface**, make sure you define suitable settings that
538 reflect ***your*** personal insulin sensitivity and data patterns.

539

540 **Caution:** As mentioned in other places, Automations can be tricky as to whether they
541 actually will ever work, because the loop goes through the exact **sequence of all your**
542 **active Automations**, and might be switched into a direction that no longer is compatible with
543 the conditions that must be a given, for the Automation you think that should kick in.

544

545 To have a clean AAPS home screen (and also to prevent unnecessary/accidental activation
546 by kids playing around with offered buttons), define reasonable time windows for each of
547 your shelved special routines, or keep them entirely dormant (de-activated) in the list of all
548 your Automations, and activate them only for/on the day when you think you might need
549 them

550

551

552

553 6.5.3 Laissez-faire alternative

554

555 You could make your life easier: **Just use** (as in Hybrid Closed Loop) **an exercise setting**
556 and accept a resulting reduced loop aggressiveness **already before meal start**. You would
557 go a bit higher in your glucose peak. As, in principle, a higher glucose level is desirable for
558 starting exercise, this can be a viable route, too, **especially if you do a** (often
559 recommended) ***protein-rich meal before exercise***.

560 Logic: For a high carb containing meal, we wanted in the preceding sections a strong initial
561 FCL response, but only up to a certain (lowered) iob. (This resembles the reduced user bolus
562 in HCL around exercise.).

563 However, the more our meal shifts to high protein/low carb, and the more we are accepting of
564 bg going a bit higher before we start or resume our exercise, the better we can tolerate a
565 reduced FCL aggressiveness also at meal start, i.e. the entire day.

566 Continue with [section 6.6](#)

567 *Skip next [section 6.5.4](#) – The following **green texts** (and embedded screenshots) describe*
568 *features that were suggested for further development, and may actually not be needed.*

569 6.5.4 Improved cockpit: Using pre-set meal / exercise combination from TT dialogue box

570 The improved „FCL cockpit“ User Interface (**when available**) also allows a one-step setting
571 for meal + exercise that can be selected in time-uncritical fashion, any time before the meal
572 starts.

573 It manages the meal with an appropriately reduced iobTH, and is programmed to
574 automatically activate the exercise settings when iobTH is exceeded:

575

576 **If in addition to meal, one of the 4 pre-programmed exercises is also selected from the**
577 **bottom of the TT dialogue box, (for example, in case of biking after a hi carb lunch, hiC + bik,**
578 **see [section 5.3.3.1](#).) then meal gets superseded /overridden with condition „duration = until**
579 **when iobTH is first time exceeded“. Plus, that is the other important point, the activity-related**
580 **reduced iobTH is taken over for the meal, too.**

581

582 **All this happens from the AAPS home screen and associated dialogue box from the TT field**
583 **there.**

584 Actual valid settings can at any time point be seen in the AAPS home screen (see [section](#)
585 [5.3.3.1](#) on extra data fields, above).

586 When your FCL is in this meal + exercise mode, you first see at the TT field ([section 5.3.3.1](#))
587 of your AAPS main screen:



... and when iobTH is the first time exceeded, this automatically switches to:



588

589

590 That 1-step setting can either be freely done according to [section 6.2](#)

591 Or you can just press one of your frequent meal and frequent exercise „codes“, as described
592 in [sections 6.3.2 and 6.3.3](#)

593 *Example: For mountain biking after pizza lunch* press two buttons, *piz* and *mtb*, in the
594 dialogue box of your AAPS home screen's TT field. That's all (...after, one time, you figured
595 out what settings suit that scenario, and you put it into /preferences, see [sections 6.3.2 and](#)
596 [6.3.3](#)).

597

598

599

6.6 Activity Monitor

An optional feature for times without serious exercise, but still suspected **effects on insulin sensitivity** is the **activity monitor**.

It can be generally activated under /preferences/OpenAPS SMB/Activity modifies sensitivity)

If the user

- has scaling factors set there (in preferences/OpenAPS SMB/Activity modifies sensitivity)

- has **no TT running**

- (and, regarding nighttime: did not opt for „ignore_inactivity_overnight“)

then AAPS automatically modulates for sensitivity changes **based on movement intensity**

for the last minutes to 1 hour time frame.

Loopers whose insulin sensitivity is affected by erratic patterns of total in-activity (desk, couch) and moving around (walking, light house and garden work) can greatly benefit from this feature. Responding in a ~ 60 minutes (rather than in 8 – 24 hours) time frame makes the Activity Monitor far superior to Autosens.

Personalized tuning of the **two scaling factors** is necessary *in your FCL set-up phase*.

For details see p.9 of the autoISF Quick guide ([section 3.4](#)).

Later, *in your running FCL*, this **will automatically adjust insulin delivery** (basal, ISF, and iobTH; see 1st screen of AAPS SMB tab!) to suit activity state of the past minutes (up to 1 hour). See also [section 5.1.5](#).

Effects from the Activity Monitor are hard-limited to go **maximally**

- **plus 20%** insulin at detected resistance (in-sensitivity to insulin) **at in-activity**
- **minus 30%** insulin at detected increased sensitivity to insulin due to **activity**.

626 **Note that Activity Monitor only works when no exercise (or other) TT is active:**

- 627 • Whenever you set a TT, you consciously go for a certain effect on the sensitivity ratio
628 to be used *in that time window*.

629 Usually it will be stronger than the tweaking that the Activity Monitor would do.
630 So, you would not want the Activity Monitor interfere, and additionally tweak
631 things you just defined differently for a certain situation, and time window, by
632 setting a TT.

- 633 • During the set TT, your Activity Monitor **keeps collecting** the data on your activity/in-
634 activity. Immediately after the set TT ends, the Activity Monitor **automatically**
635 **resumes** its work

636 This is good news also for those who might use brief even/odd target switches
637 (e.g. when sneaking-in a small snack w/o triggering a SMB), but would hate to
638 see their Activity Monitor function go under for a while afterwards.

639

640 You can easy, in real-time, check the impact of your Activity Monitor on the sensitivity used
641 (to calculate your insulinRequired, or also to auto-adjust also iobTH) in your AAPS **SMB tab**.
642 From autoISF 3.0.1 onwards, this is super easy to retrieve in the 1st screen, on top of the
643 autoISF results.


644

645 *Skip next (last) page*

646

647 Idea for an eventual User Interface update for the Activity Monitor:

648 In this dialogue box (connected *in a future update* with the exercise button), the two scaling
649 parameters (set as default by the user during initial set-up in preferences) are displayed, and
650 can be temp. over written. (These settings will expire and revert to default as set in
651 /preferences, whenever the Activity Monitor closes (goes auto-off, or is pushed off)).



653 The resulting sensitivity effect is the roughly expected effect of requiring >100% insulin if
654 moving around a bit (activity), or needing a lesser %age when being very stationary.

655 It is displayed in the right side column of the dialogue box (*if already launched*) to give the
656 user a feeling for the expected effects from her/his „weight“ inputs.

657 The exact impact is calculated by the loop and shown on top of the autoISF results in the
658 SMB tab (every 5 minutes).