

## 5. Modulation of autoISF aggressiveness. V 2.4

**Please note that with autoISF 3.0 you are in an early-dev. environment,** where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



Once the initial tuning according to [section 4](#), is done, you are ready to use autoISF for your fully automated meal management.

You will have three major *other* challenges to manage:

- recognize and manage (partial) occlusions, or other technical (CGM or BT related) obstacles (see [section 2](#) on pre-requisites of FCL, and related case studies)
- deal with times when insulin given by the loop must be restricted (e.g. a snack could be “misinterpreted” as a meal)
- deal with times when the loop should be set “milder” as a precaution (e.g. nights; or in an exercise context).

How big the remaining challenge is depends very much on your individual lifestyle. [Sections 5](#) and [6](#) discuss this in more detail.

**In order to run the loop fully automatically around the clock, the times outside the meal blocks must also be precisely analyzed, and solutions to problems (if any) must be sought.**

It is up to every user to decide where to draw the line.

- With a technically well functioning system, moderate meals, moderate or no exercise, moderate %TIR expectations and a bit of mindfulness it should be possible to go into Full Closed Loop 24/7, after working through, and observing, [sections 1-4](#).
- Especially if you are a bit shy of using the emulator for really detailed analysis, it is likely that you will not hit *one* real good system calibration ([section 4](#)) for your *entire range* of diets.

In that case you will occasionally run out of range, and your options to prevent, react, or improve are

- 35 ○ accepting a few % higher time outside range for that day (and, if feasible, in the
- 36 future avoiding what seemed to have caused it)
- 37 ○ taking a snack (whenever you tend to go low from the “tails” of insulin activity that
- 38 was required to fight a peak)
- 39 ○ doing a manual override (if you can think of one in time, to manage the problem
- 40 manually)
- 41 ○ temporarily resorting to the well-known hybrid closed loop.

42 Instead of accepting such instances, you could launch “improvement projects”

- 43 • that refine your initial tuning ([section 4.](#))
- 44 • that make you and your FCL loop fit to manage an increasing number of disturbances
- 45 either automatically, or via a user intervention ([sections 5-6](#)).

46 To tailor the loop’s response to disturbances *other-than* your major meals probably will require

47 specific **modulation of the aggressiveness** that you set according to [section 4](#) for your meal

48 management.

49 There are many avenues to achieve this. The main ones, that are also easy accessible via

50 Automations in AAPS, are:

- 51 • temporary shut-off SMBs (odd-numbered target)
- 52 • temporary reduce bgAccel\_ISF-weight
- 53 • temporary reduce iobTH
- 54 • temporary reduce set %profile
- 55 • temporary set higher TT (especially in connection with exercise mode)t

56

57 In setting up your FCL, you now have another difficult and time-consuming job at hand, to define

58 solutions for any of your „other“ situations (outside of meal management) that tend to drive glucose

59 outside of the desirable range.

- 60
- 61 • In [section 5.1](#) we explore avenues towards fully automated management that in daily life
- 62 will require no user intervention at all.
- 63 • In [section 5.2](#) and [5.3](#) we will look at solutions that involve an easy user interaction like a
- 64 data entry or button push.

65  
66  
67  
68  
69  
70  
71  
72  
73  
74  
75  
76  
77  
78  
79  
80  
81  
82  
83  
84  
85  
86  
87  
88  
89  
90  
91  
92  
93  
94  
95  
96  
97  
98  
99  
100  
101

## 5.1 Fully automatic modulation of FCL aggressiveness

The following subchapters describe set-ups you may want to use for allowing completely hands-off FCL in as many daily situations as possible.

### 5.1.1 autoISF adaptation switched off outside of meal-time windows

If, aside from having to bolus for meals, your *hybrid* closed loop was running pretty well *without* other interventions from your side, you could continue to run in that mode, and just focus your FCL (= adaptation to glucose behaviour by your new autoISF) on management of meals.

In your initial transitioning phase this approach makes a lot of sense, even by focusing autoISF adaptation to glucose behavior on *just a sub-set* of them, like only dinners.

Note that SMB range extentions (and also the SMB delivery ratio) coming with your autoISF settings *remain in place 24/7*. Depending how well tuned your former hybrid closed loop was (with respect to *not* bouncing against the *then* lower SMB size limits) you might need further precautions, like setting a temp. lowered %profile (or an adjusted temp. iob\_TH\_percent) for your time windows outside of FCL meal management.

This *intermittent hybrid closed looping* avenue is taken by many FCL users for the night times, “hanging on” to their well performing hybrid closed loop with standard oref(1) SMB+UAM (A very good alternative is FCL with night time SMB shut-off, see next [section 5.1.2](#)).

For this, you define Automations

- that set meal time windows in which “autoISF / Enable ISF adaptation by glucose behaviour” is turned **ON** in AAPS preferences/OpenAPS SMB
- or: that turn *all* autoISF’s ISF modulations (or just *bgAccel\_ISF* ) **OFF** in time windows in which surely no meal occurs. For instance, you can go for all nights back into your Hybrid Closed Loop, as you had before.

Other early DEV AAPS variants (see [section 13.3](#)) all work with meal-time windows. The window is either set by time of day in the settings, or it always must be „set“ by the user via giving a mandatory small pre-bolus before any meal starts. **Outside** of these time windows, these loops then runs with less aggressive SMBs like oref(1) SMB+UAM in AAPS Master.

This mode is not really FCL, but an advance over traditional HCL that often achieves satisfying degrees of automation and performance.

The term **Meal Announcement** (MA) is often used to label this closed looping mode. Trigger to set a meal time window could also be a pre-bolus given by the user, a carb entry made, an EatingSoonTT set, or a meal announcement button pushed.

Note: Outside of the meal time windows you would be in hybrid closed loop. To the extent you rarely face disturbances (aside from meals), you could be looping in full automatic mode around the clock,

Your temp. shut-down of ISF adaptation by glucose behavior (core autoISF function) is meant to prevent problems from the loop over-reacting to bumps in the glucose curve in times of day (night) when standardoref(1) performance is sufficient.

### 5.1.2 Odd-numbered *profile* targets used to block SMBs

An alternative route of preventing the FCL loop from over-reacting to bumps in the glucose curve would be to make use of the option to temporarily shut down SMBs

Ensure the even/odd logic in the settings is toggled on in Preferences> openAPS SMB> autoISF settings> smb delivery settings>: "Enable alternative activation of SMB depending on profile target" ON.



In time blocks with an odd-numbered profile target you can prevent any SMBs being given by your loop. The (unchanged) aggressive settings (that still will do strong adaptation of ISF to glucose behavior) then can only translate within the limits set by %TBR possible.

This will very much slow down any more insulin being given, and is an excellent solution for night times, especially if you occasionally experience compression lows.

Alternatively, you could use the new included options for Automation Conditions and temporarily tune your bgAccel\_ISF\_weight much lower ([section 5.](#)).

Yet another alternative would be to go into hybrid closed loop for the night, with or without SMBs ([section 5.1.1](#)).

In case you occasionally do have nights that would benefit from a couple of SMBs (to treat temp. highs from a late fatty pizza, raclette and such): Define suitable Automations like the two „night“ ones in this *example*:



139  
140

141 Never underestimate the „trickyness“ of getting your Automations „right“.

142 With your thought-out Automations in place, night data need to be analyzed to see

- 143 • whether the bg and iob limits defined in the given example work sensibly four your data
- 144 pattern
- 145 • whether the TT duration is chosen appropriately
- 146 • how swapping the sequence in which the automations appear in the Automation list would
- 147 lead to different SMB impacts.

148

### 149 5.1.3 Odd-numbered *temp.* targets (TT) used to block SMBs

150

151 A widely used ACTION that strongly modifies how fast your FCL can add more iob is setting an  
152 **odd-numbered temp. glucose target** which makes the loop operate without giving any SMBs  
153 (%TBR modulation only).

154 Ensure the even/odd logic in the settings is toggled on in Preferences> openAPS SMB>

155 autoISF settings> smb delivery settings>: "Enable alternative activation of SMB depending  
156 on TempTarget" ON.

157

158 So, from patterns you find in YOUR data, at times where you want your loop act differently, you  
159 need to carve out *CONDITIONS* that describe the respective situations (and either for how long it  
160 typically lasts, or at which *other CONDITIONS* you want your loop get back to default FCL  
161 operation).

162

163 An odd TT is often set for an *anti-hypo* snack or *sports* snack. In both instances, you do not want  
164 SMBs to quickly counter act.

165 In case of sweet “fun” snacks, this is entirely different -> [section, 5.2.1](#) or for regular snacks  
166 (e.g. at school break) see next [section 5.1.4](#)

167

#### 168 5.1.4 Automatic differentiation of FCL aggressiveness using Automations

169

170 **Personalized Automations** tailor the loop exactly to **YOUR** data so **fully automated handling** of  
171 situations with different aggressiveness of the loop can be made.

172

173 From, autoISF 3.0 onwards, also the following parameters are provided as *CONDITION* and/or as  
174 *ACTION* for defining YOUR Automations:

- 175 • Enable ISF adaptations by glucose behaviour => Allows temp. ON/OFF for the key ISF  
176 modulation parts of autoISF (and, as a result, will usually decrease loop aggressiveness)
- 177 • Trigger/set iobTH percent => Keeps default aggressiveness, but only until a iob threshold  
178 (that your Automation modifies) is surpassed (which is when any further SMBs will be  
179 blocked blocked)
- 180 • Trigger/set bgAccel\_ISF\_weight => Modifies the default aggressiveness of just the  
181 acceleration component

182

183 To set up suitable Automations, you first must **analyze patterns** you find in **YOUR data**, at times  
184 (or geo-locations, or bg and iob patterns that point to a problem ...) **where you want your loop**  
185 **act differently**, to carve out *CONDITIONS* that describe the respective situations (and either for  
186 how long it typically lasts, or at which *other CONDITIONS* you want your loop get back to default  
187 FCL operation).

188

189 A variant of this mode is to define several meal time windows in which autoISF aggressiveness  
190 (bgAccel\_ISF\_weight) and/or iobTH are set differently

- 191 • for different meal time slots of your day –

192 (*Breakfast at home, school lunches, school intermission snacks, dinners at home* could for  
193 example all deserve special settings regarding ISF\_weights and iobTH).

- 194 • or even for a geo-location etc –

195 (*School lunches, or mother-in-law visits, would be examples*).

196

197 Unless your meals differ vastly in size and in fast carb content all this may not be needed.

198 Still, personalized Automations might help ease your initial job of setting the various ISF\_weights,  
199 and a best-suitable iob\_threshold\_percent that would work “always”.

200

#### 201 5.1.5 Automatic adjustment of FCL aggressiveness via the Activity Monitor

202

203 With the autoISF variant of AAPS you can make use of your smartphone's **stepcounter** and use it  
204 to fully automatically adjust insulin sensitivity ratio to activity level in the past minutes to one hour  
205 time frame.

206

207 This is another little tuning opportunity, in which you study your body's response to light exercise  
208 (like walking) or to not moving at all (like desk, couch), and select appropriate settings which, in the  
209 future, will automatically adjust insulin delivery to suit activity state of the past minutes (up to 1  
210 hour).(AAPS Preferences/OpenAPS SMB/Activity modifies sensitivity/ -> set **two scaling factors**.)

211

212 This autoISF feature (new since V.3.0) is much quicker responding than Autosens or dynamicISF  
213 to adjust insulin sensitivity to your current „lifestyle state“.

214 For loopers who do not have huge variations in exercise levels in their everyday lives, this feature  
215 might fairly much close the gap towards being able to do a 24/7 hands-off FCL.

216

217 More see [sections 3.5](#) and [6.5](#)

218

#### 219 5.1.6 Pro/con completely hands-off Full Closed Loop

220

221 To stay 24/7 in a completely „hands-off“ FCL can be a realistic goal with autoISF 3.0 if besides  
222 meals also some special challenges, as discussed in this [section 5.1.](#), were analyzed and could be  
223 addressed.

224

225 Clearly it depends very much on your lifestyle, and how interested, willing, and capable you are to  
226 recognize, deal with, (and in the future avoid?) situations that get you outside of your desired %TIR  
227 on occasion.

228 So, this is also about what %TIR you are aiming at, and can accept, as it averages out for  
229 the week, for instance.

230

231 Everybody must weigh for her/himself

- 232 • how much **upfront effort** to put into getting it all 100% automatic
- 233 • **or** whether to take an **easier start, with a couple of situations left to take care of when**
- 234 **and as they arise in daily life**

235

236 Also, even if a principal capability for a fully automatic running FCL is given, this still  
237 means that

- 238 • the user should be knowledgeable about what exactly is going on, and
- 239 • has a principal capability to „nudge“, or to take over in a manual mode.

240

241 In the sections that immediately follow, we present the options to nudge or temporarily take over  
242 from the AAPS home screen that will be serving as your **FCL cockpit**:

- 243 • [Section 5.2](#) describes how you can use available “buttons” from your AAPS home screen,  
244 and how to complete it towards a suitable DIY FCL cockpit, for an even better FCL  
245 experience.

- 246 • [Section 5.3](#) describes how you might be able to manage “disturbances” even better (with  
247 more convenience) with an **improved** FCL cockpit **in the future**.

248

## 249 5.2 Modulating aggressiveness manually, from the DIY-FCL-Cockpit\*

250

251 \* Like in the airplane cockpit: Cruising in full auto mode should involve having an eye on the  
252 instruments, and on potential disturbances ahead in the environment.

253

254 In [section 4](#), we dealt with major meals. In [section 5.1](#) we looked into fully automatable manage-  
255 ment of other situations. Life in Full Closed Loop can become extremely easy then...

256

257 However: Other **disturbances** might come up, that:

- 258 • are not noticeable in-time, or foreseeable, by the loop (*e.g. your plan to start exercise in an*  
259 *hour or two*), but **that influence sensitivity dramatically** and therefore require temporary  
260 non-default settings in order to remain in-range, and/or
- 261 • **require** a different “starting point” regarding iob and bg, which translates into **a different**  
262 **iobTH** that should **temporarily** be set much lower (*in case of exercise*) or noticeably higher  
263 (*e.g. with very fast absorbing carbs in a sweet snack “sin”*) .

264

265 In [section 5.1](#) we looked into ways to automate also a modified loop response to *foreseeable* situa-  
266 tions (tied to a time of day, geo-location etc), or to those *the loop could recognize* (with enough  
267 time to react).

268

269 Other „disturbances“ might come up, and you must find an easy way to



- 270 • call up a pre-programmed routine for automatic management, with adjusted  
271 aggressiveness, or:
- 272 • manually tweak a setting or two, to temporarily adjust the aggressiveness
- 273 • There may also arise a desire to just exit the FCL mode, and be your own captain for  
274 mastering a special situation.

275 For peace of mind, to learn, and to stay informed (especially so in your initial tuning phase, or  
276 when your glucose curve goes in unexpected ways) we also must be able to

- 277 • find the key parameters that frame and drive the recent and upcoming loop decisions.

278

279 All this is facilitated within seconds right from the AAPS home screen, serving as a **FCL cockpit**  
280 after you built a couple of DIY cockpit features via Automations (as described below and in [case](#)  
281 [studies 5.2](#) and [6.2](#)):

282

283 Thoughts went also into [how to improve the cockpit in future releases](#), see [section 5.3](#)

#### 284

#### 285 5.2.1 Status recognition

286

287 Before considering any manual interventions into the ongoing FCL, you should be aware what the  
288 current mode of action is, and hence which button eventually to fine-tune or lever to switch, in or-  
289 der to adjust to the disturbance you see coming up.

290

291 See [section 5.4](#)

#### 292

#### 293 5.2.2 Manual interventions from the (DIY-) FCL cockpit

294

295 Trouble with all these is, not to forget to set back manually, too (=> better solutions in 5.3)

#### 296

#### 297 5.2.2.1 Temporary tuning of FCL aggressiveness via temp. %profile or TT settings

298

299 The set **% profile** multiplies with both, the ISF resulting from autoISF, and also with the default  
300 iobTH you have set, so both are nicely modulated in a linear way with the % temporarily chosen  
301

302 Just taking profile e.g. to 110% for an afternoon might be an easy way to explore whether  
303 you might benefit from 10% more “aggressiveness” in your core settings for lunches (like  
304 bgAccel\_ISF\_weight). Make sure, though, that the extra 10% are not cut away by set safety  
305 limits.

306

307 A lowered (relative to profile glucose target) temporary **bg target (TT)** signals lowered sensitivity  
308 (more insulin need), and an elevated TT (as often used with exercise) increases sensitivity and  
309 hence works in the direction of a lowered % profile to also reduce insulin given by the loop.

310

311 Moreover, the **exercise button** (top center on your AAPS home screen) can be activated (turns  
312 yellow, then). This will **further boost** how your set TT elevates the resulting ISF, and sharply  
313 lowers iobTH, as often desired for sports. See [section 6.1](#)).

314

315 5.2.2.2 Making temporary changes in settings made in AAPS/preferences/Open APS SMB

316

317 Going into AAPS/preferences/Open APS SMB allows to:

- 318 - set milder or stronger ...\_ISF\_weights
- 319 - set different iob\_threshold\_percent (or iobMAX)
- 320 - elevate or lower the SMB\_delivers\_ratio
- 321 - limit or expand max. allowed SMB size
- 322 - change the the even <-> odd logic for SMB on/off

323

324 Doing temporary changes in AAPS/preferences should be the exception because

- 325 - they require multiple steps, including entering a password
- 326 - you will often forget to set everything back to default a couple of hours, or minutes, later

327

328 5.2.2.3 Triggered Automations: Grey extra DIY cockpit buttons for pre-programmed “responses”

329

330 Recognizing conditions for fully automatic handling by the loop may not be not possible, or come  
331 too late for the loop to act on. Examples would be

- 332 • *exercise*: Minimum an hour before starting “the loop should know” to be able to lower iob  
333 and elevate bg by the time exercise starts.
- 334 • *snacks*: High carb snacks, sweets, consuming ice cream or having a sweet drink comes  
335 with the problem of even steeper glucose rises, but overall a lesser insulin need, compared  
336 to major meals (for which we tuned our FCL according to [section 4](#)).

337

338 This not necessarily implies that snacks need different settings than a meal. After all, autoISF  
339 was designed to act to all available data, especially to where the developing glucose curve is  
340 headed. So, depending on your effort to set parameters for a broad variety of meals (notably:

341 how well you avoid to invariably bounce fast against your iobTH), you might be able to accom-  
342 modate low carb, snack, and major meals with *one* set of settings.

343

344 For that reason, or just for increased comfort and safety, you might want to differentiate, and make  
345 use of what follows for the *sweet snack* example.

346

347 Note that both, FCL and HCL autoISF use all autoISF features, including ISF adaptation by  
348 bg curve behaviour. Only, in FCL no user bolus is given, and that has big consequences on  
349 how aggressively dialed in the bgAccel\_ISF component is.

350 In FCL autoISF, dealing with snacks (or, generally, with a very wide range of different meal  
351 sizes) is more difficult than in HCL autoISF applications, because FCL involves revving up  
352 iob supply (largely via big bgAccel\_ISF-weights) often too much to be balanced by just a  
353 snack getting absorbed.

354

### 355 **Tuning aggressiveness**

356 Key is that a sweet snack, in absence of any pre-bolus, likely benefits from even *more*  
357 *aggressive initial* FCL performance than the meals in your normal spectrum of diets  
358 require.

359 Therefore, you could set

- 360 • a higher **temp. profile%** and/or
- 361 • a temp.elevated **bgAccel\_ISF-weight** (see screenshot of my Automation).
- 362 • a **low temp. target** (76 for instance; this additionally helps maximize the first SMBs  
363 that will automatically be triggered at detection of acceleration)..

364

365 When first defining and testing this Automation, also check:

- 366 • that the safety limits as discussed in [section 2](#) will not block the intended elevated  
367 aggressiveness
- 368 • SMBs will not get outrageously big and iobTH sometimes exceeded by too much


369 Note that “the last SMB” is allowed to overshoot the effective iobTH by 30%

370

### 371 **Limiting iob**

372 For “just a snack”, total insulin need will be lower than for a meal.

373 If you would just have your sweet drink, and your meal-oriented FCL would “attack”,  
 374 iob likely would become too high, and a glucose rollercoaster would start, with you  
 375 needing to consume more =>  
 376 If you just have a snack, or drink a glass of juice, you can lower the **iobTH\_percent**  
 377 accordingly.




„User action“ is always ticked-on

This will, in the defined time space ..

..offer the “DIY cockpit” button..

..which I must press any time (~90...30 minutes) before my snack.

Note: Make sure that iobTH and bgAccel\_ISF\_weight revert to default afterwards



378

379

380 So, this can be a little **extra “project” when setting up your FCL.**

381 You need to research your snack habits (if any), and over time find out which settings in  
 382 the snack-related Automations work well.

383

384 **In everyday life you then just must press the related button in your cockpit** (which is  
 385 not time critical at all, except it should be clicked *latest* a couple of minutes after you took  
 386 the drink or snack).

387

388 If you consume more, and also eat something with your sweet drink, this will more  
 389 resemble a full meal... however, with unusual amounts of fast carbs.

390

391 Note: Pressing your snack button *a second time* would **not** help because the  
 392 lowered iobTH then would not allow iob going high enough. So, if the snack keeps  
 393 growing in size (not an unlikely scenario, right?), then you are better off just letting  
 394 your normal FCL meal routine kick in and take over, after your snack mode expired.

395  
396 Other options when snacks keep extending would require a manual modulation  
397 regarding %profile and/or bgAccel\_ISF, but keeping the full default set  
398 iobTH\_percent, or even elevating it (refer to [section 5.2.3](#)). If that happens often,  
399 define for yourself an extra User action Automation for a bigger snack (= another  
400 grey DIY cockpit button).

#### 401 402 **Installing the DIY cockpit button**

403 In the related Automation, just keep the “User action” box clicked at all times, and define in  
404 the Conditions when you want to see that button available for cockpit use (see screenshot  
405 above) => you will see that button offered.

406  
407 Besides snacks, also any **other recurring special situations can be addressed via a**  
408 **DIY cockpit button, and receive different aggressiveness up to a suitable iobTH**  
409 **level.**

410  
411 Over time you can have a big number of User action Automations, and keep them  
412 “shelved” rather invisibly (clicked in-active, top left box) in your long list of potential  
413 Automations. Even when active, they only show in your cockpit (bottom grey field of your  
414 AAPS home screen) in the time slot you assigned as potentially relevant.

415  
416 In the future you might be able to set the stage for a snack and other “disturbances”  
417 also via an extended menu below the TT button on the AAPS home screen, see  
418 [section 5.3.3.1](#)

#### 419 420 **Discussion**

421 In case you *do have* a snack habit and

- 422 • *can not* find settings, as in [section 4](#) defined for your meals, that *also* suit your  
423 snacks
- 424 • *can not* define an Automation as in [section 5.1.4](#) that would pin the snack habit to a  
425 *CONDITION* (e.g. if in time window .... bg rises ... )(then temp. set ... e.g. as in  
426 example given above)

427 *then you will need* a “snack announcement”, preferably via the grey DIY button (e.g. as  
428 shown above).

429

#### 430 5.2.4 Temporarily exiting the FCL

431

432 The “last resort” alternative always is to **temporarily** leave the FCL mode, and handle any  
433 disturbance “the traditional way” in **hybrid closed loop**.

434

435 The suggested FCL cockpit user interface with an extra version of violet loop on the  
436 AAPS home screen ([section 5.3.1](#) ) would facilitate that, including automatic removal  
437 and re-appearance of the insulin button at the bottom of the APS home screen.

438

439 In case [this feature](#) is not yet available, you must:

440 Exit the FCL mode by going to AAPS/preferences/put in your password/OpenAPS SMB/scroll down  
441 to autoISF settings and switch “Enable ISF adaptation..” OFF.

442

443 Unfortunately, there is no way yet for it to come automatically back on, after a selected time for  
444 instance. So do not forget to switch your autoISF fully back on, later.

445

446 As this will often be forgotten, it may be worth doing a “User action” Automation, for a “temp.  
447 FCL OFF” grey button (see section 5.2.1).

448 Caution though, there is very limited experience with this brand new feature

449

450 To recognize whether autoISF currently runs with ISF adaptation or not, you must consult the  
451 profile\_sens -> actual\_sens indicator below the Autosens%. However, this gets also modified  
452 by %profile switches or TT +/- exercise mode. So it is not [as easy as it would be with the “violet](#)  
453 [loop” proposal mentioned already above.](#)

454 Ultimately, you can of course study the SMB tab to find out what is going on.

455

456 Note that three things get easily confused, especially when colloquially saying “FCL is off”:

457 • FCL really at its core means the user does not bolus (and no carb entries). Looping in a  
458 largely hands-off way, with exceptions when unusual disturbances are encountered.

459 • Using autoISF means, to utilize (enable) the ISF adaptation to glucose behaviour.

460 Shutting this core autoISF function temporarily off means also exiting the autoISF FCL,  
461 because to deal with meals AAPS then will require user boli.

462 There are also non-core ancillary functions that come with autoISF (some of them may  
463 become integrated also in future AAPS Master versions). Please be aware that settings like  
464 SMB\_delivery\_ratio and also your extended SMB size limits, will apply also in your hybrid  
465 closed loop, and even if you resort temporarily to standard oref(1) SMB+UAM algorithm.

466           ⇒ This could be a problem *in case* in your oref(1) hybrid closed loop your profile  
467           ISF was set a bit too low, and you often bounced into the 120 minute SMB size  
468           limit (...but now you must tolerate the wider limit “from autoISF” ... or temp. set  
469           a lowered %profile as a parallel precaution).

470       • Shutting SMBs (temporarily) off is an independent measure to slow down delivery of  
471       correction insulin by the loop. Temp. shut-off of SMBs is built into the autoISF FCL for  
472       whenever iob> iobTH.

473       •

### 474 5.3 Modulating aggressiveness manually from the improved FCL-cockpit

475  
476 autoISF 3.0 is an early dev variant of AAPS, and as user you are participating in an on-going  
477 development. Of note, autoISF 3.0 is first launched without many of the cockpit features that are  
478 presented below in this font color. (Only what is written in black is at this point of some relevance  
479 for using autoISF 3.0)

480

481 For the time being, multi-step work-arounds may become necessary

482       • In many cases, going into AAPS Preferences and changing settings would be needed  
483       (...plus not forgetting to change these settings back, afterwards).

484       • Automations allow a DIY FCL cockpit, see [section 5.2](#) and [case studies 5.2](#) and [6.2](#)

485

486       This is also an open invitation for you to contact us in case you could help program a  
487       module for one of the required user interface extras.

488       For future integration into AAPS Master, an eye should be kept also on the question which  
489       other modes (like FCL using Automations and others mentioned in [section 13.](#)) might  
490       benefit from some of the extra features.

491

492 Keep in mind, though, that the **goal should be to interfere with the loop as little as possible**.

493 Under the described conditions it can run **fully automatically** without any user interaction (= after  
494 the initial tuning phase, and related settings made in AAPS /preferences/SMB/autoISF. See [section](#)  
495 [4.](#) and [5.1](#)).

496

497 However, just like in the airplane cockpit: Cruising in full auto mode should involve having an eye  
498 on the instruments, and on potential disturbances ahead in the environment.

499       *E.g.: storm ahead => instruct your plane to climb to another flight height.*

500       *Analogy: exercise ahead => setting an exercise TT, or => pressing a button that activates a*  
501       *sequence of instructions (some of them probably hinging on conditions, like actual iob) how*  
502       *to manage through that exercise situation).*

503

504 So, for the occasional „disturbance“ coming up, you should find an easy way to

- 505 • call up a pre-programmed routine for automatic management, with auto-adjusted  
506 aggressiveness, or:
- 507 • tweak a setting or two, to temporarily adjust the aggressiveness
- 508 • There may also arise a desire to just exit the FCL mode, and be your own captain for  
509 mastering a special situation.

510 All this is facilitated within seconds right from the AAPS home screen's **cockpit features** to the  
511 **extent they are already incorporated, or** to the extent you can build alike DIY cockpit features via  
512 Automations, as described in [section 4.1.3](#) and [case studies 5.2](#) and [6.2](#)):

513

- 514 • The button that is integrated into the **violet FCL icon** serves as emergency off button, to  
515 quickly stop FCL, or to at least to immediately stop any more SMBs (...just for a couple of  
516 minutes, or for the remaining meal time: pick from the options offered with just one  
517 keystroke).

518 Via the violet FCL icon on your AAPS home screen, you also can access a temp. switch-off  
519 button for SMBs (see section that next follows below).

520

- 521 • The **three top fields** (%profile, exercise, TT) provide access to temp. tuning of core  
522 parameters, and/or to some pre-programmed routines.

523 Taken together with some **new indicator fields** about your loop state, and the **grey DIY cockpit**  
524 **buttons** this makes the AAPS home screen your **cockpit** for Full Closed Looping.

525

526 Let us look on each of these cockpit elements in some detail:



527  
528 5.3.1 Violet FCL icon and underlying buttons


529  
530 Novices to FCL, or really anyone running into a very special situation, may appreciate that the new  
531 closed loop icon on the AAPS home screen in pink (for FCL) has buttons to quickly shut off getting  
532 more SMBs (1st row), or to enter other loop modes (second row).

533  
534 It functions very much as the other ones that you know from HCL already, and in fact you  
535 get offered some of the same options (for instance, to switch the (full) closed loop off for 15  
536 minutes for going to take a shower)

537  
538 Note that in FCL you leave all BG regulation, notably against meal spikes, to the loop. So, try not to  
539 disconnect in phases when your FCL must ramp up your job.

540 The required insulin would still be supplied *after* you reconnect. However, without the user  
541 pre-bolussing, the delay would be more of an issue in FCL than it had been in HCL.

542  
543 Just pressing on the FCL icon, a dialogue box comes up:



<-add-arrow,·form-like-like-CL-Loop··¶

<-dotted+·arrow,·form-like-like-LGS·¶

**DEV:**·when·in·Hybrid-Closed-Loop·=·top·  
headline+·green·circle), there is no Suspend·  
SMB·part·but·it·starts·right·with·„Closed·  
Loop“. In that screen, There, the 1st element·  
must be option to pick·¶

·····pink·circle,·FULL-CL-LOOP·¶

where here the green HYBRID-CL-LOOP·  
stands·¶

For loopers who did not set up FCL, a·  
feedback must come up ~ „FCL not installed“·  
if they press on that button·¶

544

545 Pressing „**Suspend SMB**“ provides fast and easy „emergency braking“ regarding delivery of more  
546 SMBs:  
547 Select the one with the desired number of minutes: 5 or 12 for just blocking the potential next  
548 SMB(s), and up to 3 hours to manage the entire rest of this meal with %TBR from then on.  
549  
550 Whenever, and why-ever, your FCL is in „no SMBs allowed“ mode (e.g. automatically after  
551 surpassing an iobTH, or triggered by a set odd TT), the FCL icon will turn into a dotted one.  
552 Instead of remaining **duration to end time** it indicates in the middle „the condition“, „**iob**“ or „**TT**“  
553 Add an indication **if** suspend SMB comes from an Automation, e.g. add an „**(A)**“, **underneath** the  
554 #minutes, iob, or TT in the middle of the dotted violet field.  
555 So, as in other (already in HCL existing) cases, those icons show in the middle the minutes left that  
556 they will be running, or the condition which would have to go away for this temp. setting to stop.  
557 It always auto-reverts into the FCL state and FCL icon, when time (or other condition) has elapsed.  
558  
559 Pressing „**HYBRID CL. LOOP**“ or other buttons from the 2nd row provides fast and easy  
560 „emergency **exit**“ **into other modes**.  
561 This enables beginners an easy „temp. escape“ into their well-known HCL (green) at any  
562 point of time. bgAccel\_ISF\_weight is set to zero when going FCL->HCL. HCL can run with  
563 autoISF (for instance dura\_ISF) uninhibited otherwise. (check **implications** for HCL users of  
564 autoISF ?? ).  
565 Note: These options from row 2 have no time limit. Loop will **not** by itself go back to FCL. You see  
566 the different loop icon as a reminder to manually revert, when ready.  
567  
568 Note (re-iterated from end of section 5.2): Three things get easy confused, especially when  
569 colloquially saying “FCL is off”:  
570 • FCL really at its core means the user does not bolus (and no carb entries). Looping in a  
571 largely hands-off way, with exceptions when unusual disturbances are encountered.  
572 • Using autoISF means, to utilize (enable) the ISF adaptation to glucose behaviour.  
573 Shutting this core autoISF function temporarily off means also exiting the autoISF FCL,  
574 because to deal with meals AAPS then will require user boli.  
575 • Shutting SMBs (temporarily) off is an independent measure to slow down delivery of  
576 correction insulin by the loop. Temp. shut-off of SMBs is built into the autoISF FCL for  
577 whenever iob > iobTH.

578

### 579 5.3.2 Buttons „Insulin“, „Calculator“ etc at bottom of AAPS home screen

580

581 These buttons are **not useful any longer in FCL**, and automatically disappear whenever in FCL  
582 mode (also in Suspend SMB state), and re-appear when leaving FCL. This applies also when an  
583 Automation or technical system failure shut off FCL.

584 Users who, maybe in the beginning phase, feel better having those buttons, can override  
585 the removal (of the insulin button, or any other) by going into /preferences/overview/buttons  
586 and forcing them on. They only remain on until the next re-entry into FCL mode, when auto-  
587 off happens again.

588 The reason why we do this: It really is important to let the loop loop, and not interfere more  
589 than absolutely needed. Any bolus the user gives will sure distort the bg curve, on which  
590 autoISF, especially when aggressively tuned for FCL, builds a lot of its decisions!

591

592

593 5.3.3. Three top fields (%profile, exercise, TT)

594

595 Depending on the variedness of lifestyle, the desired %TIR, and the initial tuning effort put in, the  
596 user may want occasionally to „tweek“ the aggressiveness of her/his FCL.

597

598 The top 3 fields (grey in default mode, **yellow when temp. in mode with changed**  
599 **aggressiveness**) serve as quick and easy entry points to make temp. switches (as users will be  
600 used to for %profile switches, or for setting an EatingSoonTT in HCL, .. which they still can do in  
601 FCL ... but more:)

602

603 Expert FCL users might need this feature rarely, but probably at least to manage activity after  
604 meals: Each require opposite aggressiveness, and the switch has to come in a certain point in  
605 time that would be difficult to capture. (More see [section 6.4](#) )

606  
 607 5.3.3.1 TT dialogue field (( Currently not available in the pictured form and function ! )  
 608 *f extended design for FCL cockpit is already launched)*  
 609 The TT field (top right of AAPS home screen) is the primary daily interface, and a dialogue field  
 610 opens when pressing on it

Initially, the form can be entirely empty re. TT inputs and just show the default iobTH and bGAccel parameters.

Two lines appear for target and get labeled ES, AC or HY depending on what was pushed at the bottom buttons. (HY => red frame; evtl 1 empty => grey)

TT entries can be made or overridden.

iobTH calculates from Target AC and other settings shifts away from default set in /preferences (here 0.7) but could be overridden here 0!..200 % is allowed

bGAccel ISF\_weight can be modulated here, too. Note: it can change again if % profile is also changed. --0...200 % is allowed

CANCEL allows to start fresh (select one or two of the square buttons, ES, AC or HY)

OK needed to use the settings

Duration input is made in minutes. In the exceptional case that both, ES and AC targets are defined, the duration input is for AC and framed blue. (This is because the preceding AC mode is automatically determined in length by the loop observing when iobTH is exceeded)

611  
 612  
 613 This looks complicated but only because it allows 4 different modes of use. Each user will primarily  
 614 use her/his preferred one.  
 615 (1) Who is happy with the initially well tuned FCL and does not have huge variations in daily eating  
 616 and moving around, will **not use** the TT **at all**. FCL is possible without an intervention via the

617 TT button in your cockpit. Actually 4 of 8 modes (GGG ...YYY permutations, list see [section](#)  
618 [5.3.1](#)) are not making use of TT.

619 (2) Super easy is also, to just input **any odd-TT** (odd-numbered temporary target) that will shut out  
620 any SMBs for the set duration. *That can be a good idea when having a snack, for instance.*

621 *Super quick access to stop SMBs is possible also via the loop icon ([section 5.2.1](#)).*

622 Specifically, an **EatingSoon TT** can be activated here (limited relevance see [section 2.5](#)). It is  
623 time-un-critical, can be manually set, or come up via an Automation.

624 *The cockpit enables you to set the iobTH differently (override) for the current meal.*

625 *Alternatively, iobTH can be temporarily changed in /preferences or using an Automation.*  
626

627 *Temp. iobTH will always revert to default when the TT expires. If another TT immediately*  
628 *follows, like in the example of the screen above, it will calculate, (then) show and use a new*  
629 *temp. iobTH.*

630 (3) The third way is to **use the input mask** (*if already ncluded in your software version see picture*  
631 *above*) **to freely modulate the loop aggressiveness** for a declared number of minutes. Click  
632 the bottom big square(s): Either HYPO, or ACTIVITY, or EATING SOON, or ACTIVITY and  
633 EATING SOON (*example in the pictured screen above*). Make or override entries in the offered  
634 fields. Press OK.

635 (4) The fourth way is to exclusively use one of the 4+4+2 little buttons seen in the bottom part of  
636 the TT dialogue box (*if already included in your software version*). They provide a set of  
637 settings (as will immediately show in all input fields above) that the user has set up in  
638 Preferences/SMB/autoISF/FullLoop (refer to [section 6.3](#)), and can freely label there. *For*  
639 *instance „hiC“ at high carb EatingSoon, „piz“ for Pizza/fatty meals, „grd“ for garden work,*  
640 *„mtb“ for mountain biking ...*

641 **Capturing good settings for not-everyday situations in /preferences** (*if already included*)  
642 **allows calling them up within 1 second**, from your cockpit on the AAPS home screen (...and  
643 won't ruin the FCL experience at all , especially because in most cases it is not time-critical,  
644 how long before the intended exercise the buttons are pressed).

645 [Case study 6.2](#) demonstrates that nearly the same performance and comfort can be reached  
646 via the **DIY FCL cockpit** with the grey extra buttons appearing at the bottom of the AAPS home  
647 screen, based on Automations with User action (see also [section 4.2](#)).

648 *The example picture given above, and also [case study 6.2](#), is the most complicated (but also most*  
649 *useful) case, **when exercise follows after a sizeable meal**. It is then that you need (a) aggressive*

650 FCL initial performance at the meal, but, exactly when (!) a (for the intended sport already  
 651 temp.lowered) iobTH is exceeded, you need (b) to have SMBs automatically switched off and go  
 652 into the „milder“ mode, as defined for the exercise (with high instead of lowTT, that automatically  
 653 significantly reduces iobTH again, and insulin sensitivity(resistance) settings too).  
 654 Pressing exercise related buttons will automatically also light the **exercise button** on the main  
 655 screen yellow.

656  
 657 To summarize, the TT dialogue field offers easy but powerful ad-hoc modulation of loop  
 658 aggressiveness for FCL (if already included).

659  
 660  
 661 5.3.3.2 Exercise button (see more in [section 6.](#))

662  
 663 The exercise button automatically lights yellow when exercise related TTs are activated in the TT  
 664 dialogue box. 4 of 8 modes are making use of the exercise button.

665 If pressing on the exercise button, a dialogue box appears (*if extended design for FCL cockpit is*  
 666 *already launched*) with info on exercise setting first (and opportunity to override), plus below the  
 667 activity monitor (experimental for auto-tracking of lighter movement during the day, and effects on  
 668 sensitivity that may have. See [section 4.5](#)).

669  
 670 So, first the exercise settings (as set under TT) are there to read. Example :

activity	TT	dura	%-profile	iobTH	bgAccel	%overall
mtb	171	180	70	3.0-U	0.24	70

Mode set to run for 134 more minutes

Mode starting after meal when iob > 4.4 U or n/a

671  
 672  
 673 The exercise (here mtb) is selected in the dialogue box of the neighboring TT field, and there auto-  
 674 filled with settings made in the set-up and tuning stage by the user under preferences (see  
 675 above,...). They are reported also under the exercise button here, and TT, duration, and % sens  
 676 (which also shows active on the %-profile field on the left side of the exercise button) can be temp.  
 677 changed there. iobTH, bgAccel\_ISF and overall resulting sensitivity ratio is given in the other fields.  
 678 The **middle field** of the table, „% profile“ either picks up the % set under the %-profile button, or  
 679 an input can be made here, in the exercise button domain, which will:



- turn the neighboring %profile button on yellow and show that inputted % on it, too
- be multiplied with the result from the exercise mode settings per se, and change the % overall, accordingly.

So, if this middle field of above table (dialogue box of sports button) contains a figure other than 100, input field becomes yellow, and you are operating with a combination of traditional PLUS new exercise mode (with all three top buttons of your FCL cockpit yellow). This maximally will soften aggressiveness, for which you get an idea by the last calculated figure.

The mode is either running already (for another number of minutes, as probably also shown in the yellow TT field anyways). Or it is scheduled to run, after insulination for a started meal reaches iobTH (as in table). Or, no exercise is scheduled (both points red, no entries).

The lower part of the exercise dialogue box (not pictured above, but see in [section 6.5](#)) is dedicated to the Activity Monitor

#### 5.3.3.3 Profile button

The profile button can still be used to set a different profile, or profile%, for instance to adjust for days with sickness (as you are used to from hybrid closed looping). 4 of 8 modes are not making use of the profile button.

Any inputs made here will be used to modify profile\_ISF on which all further changes are made on (multiplied with).

The profile field remains grey if standard profile is applied.

It turns yellow, displaying a %figure relating to any altered loop overall aggressiveness:

- When no inputs (changes from 100% profile) are made here, but inputs in the TT field, e.g. for exercise, automatically lead to different insulin sensitivity ratio ~~that ratio is shown here~~
- when% is changed by input in the profile button itself, it will be multiplied with with profile\_ISF and be used in place of profile\_ISF *by the algorithm*.

However, for exercise (sports) you no longer must make an entry here, because reasonable %reductions should be automatically provided, driven by your set TT (and half-basal exercise target), see [section 6](#).

## 716 5.4 Recognizing your loop state in the AAPS home screen

717

### 718 5.4.1 Color scheme of top cockpit buttons tells kind of closed loop that is running

719 3 Buttons (%profile; exercise; TT) each in 2 states (yellow Y, or grey G) makes  $2 \times 2 \times 2 = 8$  possible  
720 combinations:

721 GYY = dynamic exercise mode

722 YGY = not-dynamic „traditional“ exercise mode (if <100%) or hypo mode (if >100%)

723 GYG = basic closed loop with Activity Monitor running

724 GGG = basic closed loop (FCL or HCL) without any altered sensitivities etc

725 YGG = basic closed loop but with a „long wave“ sensitivity shift (e.g. sickness)

726 GGY = temp. target like e.g. EatingSoonTT is set; or Hypo mode

727 YYG = closed loop with „long wave“ sensitivity adjustment and Activity Monitor running

728 YYY = dynamic exercise mode in time with additional „long-waved“ sensitivity shift

729

### 730 5.4.2 Information printed on the top buttons

731

732 The yellow TT field shows the currently valid TT (and further duration):

733 (profile) stands for the abbreviation you labeled your selected running profile



735 In the special case of settings for meal preceding sports, the field will look slightly differently:



...and ...

737 ... when iobTH is first time exceeded, this automatically switches to:



739

740 Likewise, if on the AAPS main screen just an **EatingSoonTT** is set (e.g. 72), this is entered with the  
741 desired duration. Afterwards, it automatically reverts to profile target and the display turns grey  
742 again there with e.g. 90 on it (and no time limit).

743 Without sports context, the middle field remains grey.





744

745

746 Independently from setting a TT, the user can choose to set a **%profile in the left top field**, for an  
 747 independent number of minutes, e.g. 70% in this screen example: Also, or additionally, this will  
 748 influence the resulting ISF and sensitivity%

749



750

751

752 The % might change and turn yellow also in context of making TT inputs in the related dialogue  
 753 box (see chapter TT dialogue field, above). Still, the % (or the length of time the profile switch shall  
 754 be active) can be independently overridden in the top left field, if so desired.

755

756 If an **Automation** sets a %profile, and/or a TT (e.g. *automatic detection of meal start at condition*  
 757 *e.g. when delta >10*), this would automatically show in respective field(s) turning yellow and  
 758 showing the temp. setting. To show the set parameter comes from an Automation, „ **(A)** „ is added  
 759 in the end of button text.

760 Note that an Automation is **usually / always (?)** only permitted to temp. change default profile  
 761 settings, not other pre-existing temp. settings. This is for a good reason : Why should a  
 762 sometimes in the past thought-out Automation supersede your - just for the occasion  
 763 specified – temp.settings that you consciously activated for the day?

764 Advice: Try to stay away from Automations that also aim at temp. modifying  
 765 aggressiveness. For the reason just given in above note, they often will not kick in anyways.  
 766 Generally, it also is no good idea to double up sub-algorithms for tweaking loop behaviours.

767

768 Try to **keep things as simple and clear as possible**.

769

770 That said, a limited number of Automations can be of help in distinct scenarios (that differ in  
 771 purpose and in applicable time of day).

772 A good one could be for night time, when your odd profile TT has SMBs shut off, but your  
 773 experience after pizza nights tells you that, under certain condition patterns (bg, iob), an  
 774 SMB or two should be „allowed in“. Another good example, if you go usually FCL without

775 any use of the TT button (which you could call a meal announcement of sorts), is to define  
776 an Automation that, after detecting a meal start, automatically sets a low TT to get  
777 maximally aggressive first SMBs.

778

#### 779 5.4.3 FCL related indicator fields in the AAPS home screen

780

781 In extra data fields of the AAPS main screen you can always see (not change) the key  
782 „aggressiveness“ parameters your loop currently operates with (see also home screen  
783 example below):

784 • how profile sensitivity (**ISF**) changes by the %profile input, by autoISF, and/or a set  
785 exerciseTT.

786 • next to current available iob number is an indication of your **valid iobTH** (the iob above  
787 which no more SMBs will be given)

788 • The AAPS home screen additionally shows, above the deltas, the current **acceleration**

789 Having a look at that can be valueable. For instance, when glucose is relatively low and still  
790 falling, a positive (and getting more positive) acceleration indicates that bg will swing back  
791 up, rather than crash low. This will give info about necessary snack size, and hence help  
792 avoid both, unnecessary calories, and going on a bg roller coaster.

793  
794  
795 5.4.4 Overall home screen:  
796

Overall home screen:



797  
798  
799