

12. Remarks for users of previous autoISF versions V 3.0

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



12.1 Main innovations in autoISF 3.1.0

Except for things relating to 1-minute Libre readings, the main point of migrating from the previous AAPS 3.2.0.4+autoISF 3.0.1 to the current newest AAPS 3.3 3.a+autoISF3.1.0 was really to update autoISF to the changes done in AAPS Master.

For instance, there is an extra AUTO ISF tab now, instead having different SMB tab designs depending on selected modes.

The FCL e-Book needed only very minor adaptations. Changes in autoISF are sketched below, more see in autoISF Quick Guide @ [Github/ga-zelle/autoISF](https://github.com/ga-zelle/autoISF).

From there:

The focus of this release was on merging AutoISF into recent AAPS3.3.3-dev-a as well as adding [calibration](#) and [smoothing](#) for Libre sensors without using xDrip. Other minor changes include:

- The mainstream AAPS3.3 logfiles were moved back to the AAPS install directory. Under Android14 and above this folder is heavily protected and the emulator on the phone was no longer able to access them. Therefore in AutoISF the logfiles get stored in "...emulated\0\Documents\aapsLogs".
- Now in AutoISF the logfiles contain data for one full hour and are named like "AndroidAPS._yyyy-mm-dd_hh.log.zip". This means it is easier to find the relevant logfile if the time of the event is known.
- Logfiles for one month will be kept before getting recycled.
- There is a PR for automation states not yet included in mainstream AAPS3.3. but included here with a pilot application related to [sleeping phase](#) in the activity monitor.
- As part of the Libre calibration and smoothing methods there is a trigger to disregard SMB and AutoISF weights for the next 20 minutes.
- When running *cleanup database* the period to keep can be [scaled down to just 1 week](#) of data instead of 3 months.
- The [half basal exercise target](#) can now be entered in mmol/L as well.
- Gaps in CGM readings up to 11 minutes are now accepted for the [parabola fit](#).

12.2 Main innovations in (previous) autoISF 3.0.1

- Integration of **Libre3** in **1-minute** mode
- **Cancellation of some parameters and options** (none of which were recommended to use in past FCL guidance documents):
 - *pp_ISF_hours* no longer required because ...
 - *enable_pp_ISF_always* is now always true which means ...
 - *delta_ISFrang_weight* is no longer used in favour of *pp_ISF*
 - *enable_dura_ISF_with_COB* is now always true
 - *enable_SMB_EvenOn_OddOff* was discontinued and unified with *enableSMB_EvenOn_OddOff_always*
- About **3 times as many logfiles** can be stored in the phone, making copying zip files into your PC less often necessary
- Clearer structure of the **SMB tab**
 - now with all autoISF effects detailed on top page
 - iobTH declared
 - language clearer:

Message	Condition	What does it affect?
Loop allows maximum power	even target < 100	increase in bg limited to 30%, otherwise no SMB; actual SMB delivery ratio is max of fixed smb_delivery_ratio and linearly growing ratio
Loop allows medium power	even target >= 100	increase in bg limited to 20%, the APS default, otherwise no SMB; actual SMB delivery ratio is either fixed smb_delivery_ratio or linearly growing ratio
Loop allows minimal power	odd target	no SMB, only TBR available for action
Loop power level temporarily capped	IOB > effective iobTH	temporarily no SMB, only TBR available for action; IOB is above user defined iobTH, potentially modulated by exercise mode, activity monitor and profile percent
Loop allows APS power level	no even/odd target option active	SMB enabled/disabled according to standard APS rules and settings; no iobTH threshold active

-
- AAPS home screen on smartphone now shows **autoISF result** (sens/profile.sens = amplification factor on the profile_ISF) **underneath** the **Autosens %**

See also preface in the Quick Guide: https://github.com/ga-zelle/autoISF/blob/A3.2.0.4_ai3.0.1/autoISF3.0.1_Quick_Guide.pdf

Also an error in the previous 3.0 version was rectified:

44 In 3.0 there was a bug related to the *iob_threshold_percentage* being inactive when setting it to 100%.
45 Once modulated by exercise or similar it became active unintentionally

46

47 12.2 Implications regarding your previous FCL settings

48

49 The changes in autoISF were minor.

50 **Please refer to AAPS readthedocs regarding the important changes in the migration to the**
51 **newest AAPS version.**

52

53 The following related to the previous autoISF update

54 Odd target => SMB-off

55 As the even/odd logic now uniformly applies to profile targets and to TT, users who had previously only used the odd
56 TT/SMB shut off MUST now extend this setting to profile targets.

57 This means, you need to look through your bg targets in your profile and may have to adjust the bg target for some hours
58 by 1 unit (or decimal point in mmol).

59 Cancellation of delta_ISF

60 If you had followed the suggestions – also in previous versions of this e-book - you would not have made use of any of
61 the cancelled feature.

62 Else you may have to switch to pp_ISF now for management of the near-linear bg rise phase in some meals. See [section](#)
63 [4.4](#)

64 Storing more logfiles in the internal phone memory

65 You can wait about 3 times longer than you were used to, for securing your logfiles in your PC.

66 Lower hurdles for “looking under the hood” of your autoISF

67 Looking the current aggressiveness your loop works under (ISF and iobTH actually in use) is made much easier in the reorganized
68 *SMB tab*.

69 The updated FCL e-book now comes with detailed guidance how to use the **emulator**, both on PC and on (Android) phone ([sections](#)
70 [10](#) and [11](#).)