

12. Remarks for users of previous autoISF versions V 3.2

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



12.1 Main innovations in autoISF 3.2.0

Currently (July 2026) we are in a transition phase where the main AAPS is upgraded to version 3.4., with major redesigns. At the same time there were some pod-related problems that got resolved.

So, we are working on getting autoISF upgraded and integrated, as well.

See the current status here: https://github.com/ga-zelle/autoISF/blob/A3.4.0.0_ai3.2.0/Features%20only%20available%20in%20this%20stand-alone.md and, specifically, the latest autoISF 2.0 guide here: https://github.com/ga-zelle/autoISF/blob/A3.4.0.0_ai3.2.0/autoISF3.2.0_Quick_Guide.pdf

Files

A3.4.0.0_ai3.2.0

Go to file

BerNie's guide for Full Closed Loop...

Calibration_proforma.ods

ExerciseMode.ods

ExerciseMode.xlsx

Features only available in this stand...

How-to-get-larger-SMBs.pdf

Libre_serpent.mp4

README.md

To prebolus or not to prebolus.pdf

autoISF3.2.0_Kurzanleitung_ENTWU...

autoISF3.2.0_Quick_Guide.pdf

autoISF / Features only available in this stand-alone.md

Preview Code Blame 29 lines (25 loc) · 1.26 KB

FUNCTIONALITIES

- graphical representation of the history of AutoISF factors and iobTH
- graphical overlay of fitted parabola in the main glucose graph
- smoothing of the full 1-minute Libre dataset when coming directly from Jugglucio
- calibration of data from Libre or G7 coming directly from Jugglucio
- pausing SMB and AutoISF weights during a phase of recalibration
- using Automation States
- using automation triggers and actions specific to some AutoISF variables
- display Script Debug section on follower phone
- fixed step trigger automation

EASE OF USE

- readability of AutoISF tab
- quick access to exercise mode and displaying its status
- easy access to logfiles on the phone
- keep 1 month of logfiles
- clean glucose records in database down to just last 7 days
- alternatively enter *half basal exercise target* in mmol/L
- automatic reminder to enable Physical Activity permission for phone step counter

21 12.2 Main innovations in (previous) autoISF 3.1.0

22

23 Except for things relating to 1-minute Libre readings, the main point of migrating from the previous
24 AAPS 3.2.0.4+autoISF 3.0.1 to the current newest AAPS 3.3 3.a+autoISF3.1.0 was really to
25 update autoISF to the changes done in AAPS Master.

26 For instance, there is an extra AUTO ISF tab now, instead having different SMB tab designs
27 depending on selected modes.

28 The FCL e-Book needed only very minor adaptations. Changes in autoISF are sketched below,
29 more see in autoISF Quick Guide @ Github/ga-zelle/autoISF.

30 From there:

The focus of this release was on merging AutoISF into recent AAPS3.3.3-dev-a as well as adding [calibration](#) and [smoothing](#) for Libre sensors without using xDrip. Other minor changes include:

- The mainstream AAPS3.3 logfiles were moved back to the AAPS install directory. Under Android14 and above this folder is heavily protected and the emulator on the phone was no longer able to access them. Therefore in AutoISF the logfiles get stored in "...emulated\0\Documents\aapsLogs".
- Now in AutoISF the logfiles contain data for one full hour and are named like "AndroidAPS._yyyy-mm-dd_hh.log.zip". This means it is easier to find the relevant logfile if the time of the event is known.
- Logfiles for one month will be kept before getting recycled.
- There is a PR for automation states not yet included in mainstream AAPS3.3. but included here with a pilot application related to [sleeping phase](#) in the activity monitor.
- As part of the Libre calibration and smoothing methods there is a trigger to disregard SMB and AutoISF weights for the next 20 minutes.
- When running *cleanup database* the period to keep can be [scaled down to just 1 week](#) of data instead of 3 months.
- The [half basal exercise target](#) can now be entered in mmol/L as well.
- Gaps in CGM readings up to 11 minutes are now accepted for the [parabola fit](#).

31

32

33 Main innovations in (previous) autoISF 3.0.1

34

- Integration of **Libre3** in **1-minute** mode

- **Cancellation of some parameters and options** (none of which were recommended to use in past FCL guidance documents):

37

- *pp_ISF_hours* no longer required because ...
- *enable_pp_ISF_always* is now always true which means ...
- *delta_ISFrange_weight* is no longer used in favour of *pp_ISF*
- *enable_dura_ISF_with_COB* is now always true
- *enable_SMB_EvenOn_OddOff* was discontinued and unified with *enableSMB_EvenOn_OddOff_always*

38

- About **3 times as many logfiles** can be stored in the phone, making copying zip files into your PC less often necessary

40

- Clearer structure of the **autoISF** tab
 - now with all autoISF effects detailed on top page
 - iobTH declared
 - language clearer:

Message	Condition	What does it affect?
Loop allows maximum power	even target < 100	increase in bg limited to 30%, otherwise no SMB; actual SMB delivery ratio is max of fixed smb_delivery_ratio and linearly growing ratio
Loop allows medium power	even target >= 100	increase in bg limited to 20%, the APS default, otherwise no SMB; actual SMB delivery ratio is either fixed smb_delivery_ratio or linearly growing ratio
Loop allows minimal power	odd target	no SMB, only TBR available for action
Loop power level temporarily capped	IOB > effective iobTH	temporarily no SMB, only TBR available for action; IOB is above user defined iobTH, potentially modulated by exercise mode, activity monitor and profile percent
Loop allows APS power level	no even/odd target option active	SMB enabled/disabled according to standard APS rules and settings; no iobTH threshold active

-
- AAPS home screen on smartphone now shows **autoISF result** (sens/profile.sens = amplification factor on the profile_ISF) **underneath** the **Autosens %**
- Note: More features for autoISF 2.0 and later. The above section will be updated. Meanwhile refer to the Quick Guide: https://github.com/ga-zelle/autoISF/blob/A3.4.0.0_ai3.2.0/autoISF3.2.0_Quick_Guide.pdf

12.3 Implications regarding your previous FCL settings

Please refer to autoISF Quick Guide (link given in [section 12.1](#)) regarding changes in the newest autoISF version.

Please refer to AAPS readthedocs regarding the important changes in the migration to the newest AAPS version.