

Case study 4.3: Hands-off FCL using autoISF on Xmas V.2.2



A pre-release tester reports here on trying, for a 5-day period in the 2023 Christmas season, a **completely hands-off** FCL utilization.

Method

FCL with (then not yet released) dev variant of AAPS 3.2.0.4 w/autoISF 3.0.1
Lyumjev 100 (DIA 7h) in Combo pump w/ 10mm Teflon cannula (0-48h)
2 x G6 overlapping (see case study 1.5; sensors used ~ d3 – d15; xDrip; no smoothing in AAPS)
profile basal ~ 14 U (0.41...0.75 U/h); profile_ISF 36...44 mg/dl/U (circadian);
TDD 43 U (which is above his normal TDD, see p.2).

Main settings for 24/7 adaptation of ISF *):

- SMB range extension and autoISFmax both = 2.9; SMB delivery ratio = 0.75 fixed
- bgAccel_ISF_weight = 0.24; pp_ISF_weight = 0.03; dura_ISF_weight 0.8
- iobTH_percent=60

***) Caution: Do not copy settings** from others,
not even for starting your tuning.

Why, see FCL e-book [section 4.1](#).

No user interference:


- no boli
- no carb inputs
- no setting Eating Soon
- no “manual nudging” via the top fields in AAPS main screen (%profile, exercise, TT).

Loop did use the following Automations (from the initial FCL set-up, as described in FCL e-book):

- *For meal management:* Automation that sets TT=74 mg/dl for 26m if delta>10 (FCL e-book section 2.5)
- *For nighttime management* “Adjunct Automations.” (discussed in FCL e-book section 5.1.2):
 - Operate with an odd profile target at 3am-10am for default SMB shut-off...
 - ...complemented by Automations to get temp. SMBs (if bg>160; iob<3.5) during nighttime
 - Note: If tester normally would do breakfasts, set odd profile TT would end earlier than 10 am
- **Activity monitor** 24/7 on. Scale factors 1.2 activity and 0.3 inactivity.

96% TIR (HbA1c ~ 5.6%) during 5 days at Christmas

Glucose distribution (5 days total, Friday 22.12.2023 - Tuesday 26.12.2023)




Negligible hypo tendency

The statistics for the “Low < 70 mg/dl” line of above table shows:

- Only 1.3 % of values were below 70 mg/dl
- Half of those (0.7%) were under the group-median of 66 mg/dl
- Virtually no data point below 62 mg/dl (= Median low minus 2 SD)(61 was lowest in Daily stats)

The percentile report shows, that briefly after midnight the low tendency is biggest (down regulation from dinner coincides with a last dogwalk):



Details for each of the 5 included days:

Due to some **excessive eating**, TDD was on average 42.8 U (+16% to normal TDD


averaging 37U). No big breakfasts. Dinners came relatively late; and once, on Dec 25/26th,

there was midnight chocolate fondue w/ 2 gl. sweet hot wine (Glühwein)


Two very big holiday lunches (25th and 26th), goose resp. turkey, with dumplings, vegetables

and small desserts; Christmas cookies. No settings were adjusted.


Daily stats report




Tuesday 26.12.2023




Bolus insulin: 42.0U
Base basal insulin: 13.7U
Positive temp basal insulin: 3.9U
Negative temp basal insulin: -7.7U
Total basal insulin: 9.8U
Total daily insulin: 51.8U
Total carbs: 0 g
Total protein: 0 g
Total fat: 0 g




Monday 25.12.2023




Bolus insulin: 34.2U
Base basal insulin: 13.7U
Positive temp basal insulin: 5.1U
Negative temp basal insulin: -7.7U
Total basal insulin: 11.1U
Total daily insulin: 45.3U
Total carbs: 0 g
Total protein: 0 g
Total fat: 0 g



Sunday 24.12.2023



Bolus insulin: 23.5U
Base basal insulin: 13.7U
Positive temp basal insulin: 5.9U
Negative temp basal insulin: -7.1U
Total basal insulin: 12.5U
Total daily insulin: 36.0U
Total carbs: 0 g
Total protein: 0 g
Total fat: 0 g



85 What the user should still do, is:

- 86 • Look occasionally after **BT connectivity** (especially after meal starts)
- 87 • Look occasionally into the bg and iob (or insulin activity, thin yellow curve)
- 88 development, and develop some “mindfulness” so to some extent, and only
- 89 sometimes, this might influence the eating/**snacking habit** a bit.
- 90 • With very special “disturbances” the hands-off FCL runs into limitations.
- 91 Notably if the need for extra snacks shall be kept low, ahead of **exercise** it can be
- 92 essential to take special precautions for limiting iob and to elevate the glucose target
- 93 (as known from hybrid closed loop). See example in [case study 6.2](#)
- 94 • If more “serious” exercise would have been included, the tester would have used some form of
- 95 “exercise announcement”, and/or would have required snacks to avert hypos.

96