

13. Other Avenues to FCL

V 2.2

Please note that with autoISF you are in an early-dev. environment, where the user interface is **not optimized for safety** of users who stray away from intended ways to use. Good safety features exist, but these are only as good as the development-oriented user understands and implements them. This is not a medical product, refer to disclaimer in [section 0](#)



- 13.1 FCL using AAPS Master and Automations
- 13.2 dynamicISF used for Full Closed Loop.
- 13.3 Methods involving simple Meal Announcement that might be stretched into a FCL
 - 13.3.1 Boost,
 - 13.3.2 AIMI,
 - 13.3.3 EatingNow
 - 13.3.4 Tsunami
- 13.4 No-bolus looping with precise carb inputs
- 13.5 Machine Learning (AI)
- 13.6 Dual Hormone Systems

Available related case studies:

- Case study 13.1: Comparison 1 month FCL Automation vs autoISF
- Case study 13.2: FCL using dynamicISF (call for an example so far un-answered)
- Case study 13.3: FCL using Boost

13.1 Full Closed Loop using AAPS Master and Automations

AndroidAPS 3.0 was (Sep.2023) the first DIY system to launch Full Closed Looping as an option to manage T1D, if a described set of pre-requisites apply.

Key pre-requisites were described in

<https://androidaps.readthedocs.io/en/latest/Usage/FullClosedLoop.html> , and are sketched also in [section 1](#), with [case studies 1.1 – 1.5](#) underscoring the importance.

You may (not) have noticed: There was no big „marketing fuzz“ made around that FCL option. Seeing how many AAPS users struggle with even getting their basal, ISF and SMB settings right, it would be foolish to allure everybody to a supposedly very easy way of looping. True, it can be easy. But only after doing a personalized set-up project. Setting up is easier than what autoISF and the methods we get to in [section 13.3](#). demand, but still a project. It also requires a well mastered hybrid closed loop, to start from.

With attention to the pre-requisites, and avoiding extreme high carb diets, many (mostly: adult) users achieve satisfactory %TIR after supplementing AAPS Master with personalized Automations that attempt to strongly elevate iob upon recognition of a meal-related bg rise.

42 See also Case Studies, and the randomized cross-over study involving AAPS FCL: PubMed [First](#)
43 [Use of Open-Source Automated Insulin Delivery AndroidAPS in Full Closed-Loop Scenario:](#)
44 [Pancreas4ALL Randomized Pilot Study](#);

45

46 This method is **highly recommended for an entry into FCL for those who do not have the**
47 **interest, or lack the time, to deal with the very much more sophisticated and demanding other**
48 **routes** towards FCL, like autoISF, or also like the methods briefly presented below in [section 13.3](#) .

49

50 Note that using the autoISF dev version of AAPS 3.2 (*with "Enable ISF adaptation.." OFF*) can
51 be a good idea, to make use of features like SMB_range_extention and SMB_delivery_ratios >
52 0.5. Compared to using AAPOS Master, this allows stronger boosting of SMB sizes, also when
53 *not making use of autoISF* but just Automations, for FCL.

54

55 13.2 FCL using dynamicISF with AAPS or with iAPS

56

57 As opposed to

- 58 • autoISF, with it's bgAccel_ISF component , or to....
- 59 • AAPS Master, with Automations strengthening ISF at meal-related bg rises ...

60 dynamicISF was **not** designed to help boost SMBs asap after an omitted user bolus.

61

62 Rather (as the name also suggests) it was designed to be used in hybrid closed looping to make
63 ISF react more dynamic to suspected swings in insulin sensitivity (which shows in bg values, and
64 in TDD trends). It does a similar job like Autosens, but can be much more amplified (by the users
65 tuning their dynamicISF adjustment factor (%)).

66

67 When using a fast insulin (and when some other pre-requisites discussed in section 1 are in place,
68 too), the dynamicISF method can be applied also to Full Closed Looping. (See [Case Study 13.X:](#)
69 [not available by time of publication => this is a call for a dynISF FCL user to provide a case study](#)
70 [that contains a 1 week 24h scatter plot as well as one analyzed meal where we can see when and](#)
71 [how dynISF helped build iob, after not having bolussed](#)).

72

73 It will have a principal timing-disadvantage because responses are more tied to high bg values
74 than to acceleration (in autoISF) or to delta (in the Automations route to FCL).

75

76 On the other hand, people who 1) do have strong sensitivity swings and 2) cannot pro-actively
77 deal with those (e.g. by making profile switches) might be satisfied with the automatic (although a
78 bit late) adjustments that dynamicISF automatically will provide.

79

80 dynamicISF therefore could be characterized, in the FCL context, as a potential solution to a rather
81 care-free approach for those who do not seek best-possible performance (or who take other
82 measures, like low carb diet, to still reach pretty acceptable performance in FCL mode).

83

84 **More info** (caution, both not focussed on FCL:)

85 AAPS / search term dynamicISF in: <https://discord.gg/DfvK5HnxXu>

86 iAPS / section dynamic-isf-cr: <https://discord.gg/gGKXW5uX3m>

87

88

89 **13.3 Methods involving simple Meal Announcement that might be stretched** 90 **into a Full Closed Loop**

91

92 See also [section 7](#) on using **autoISF** in “MA” mode, involving a pre-meal bolus.

93

94 **13.3.1 Boost**

95

96 All of the additional code outside of the standard SMB calculation requires a daily time period
97 („Boost window“) to be specified within which it is active.

98 A variation of dynamicISF is used in which also predicted bg will be considered in varying degrees
99 (40...75%) to mimic the effects of higher insulin sensitivity at lower glucose levels.

100 When using Boost without carb inputs (permanent cob=0) a special **boosting of SMBs** is provided
101 when an **initial bg rise** is detected with a meal:
102 delta, short_avgDelta and long_avgDelta are used to trigger an early bolus (assuming IOB is below
103 a user defined amount).

104 This procedure goes in the direction of the bgAccel_ISF route discussed for autoISF
105 ([section 4.1](#)). If used with an excellent CGM, autoISF acceleration detection should be a bit
106 earlier, and boosting can be made much stronger in autoISF

107 For safety, the user sets a value of 2.5% (up to 5%) of TDD for the max. Boost Bolus (Boost Bolus
108 Cap).

109 For stronger boost, the default AAPS 50% SMB_delivery_rate can be overwritten with a higher in-
110 sulin percentage determined by the user. The SMB_delivery_ratio is called „Boost insulin required
111 percent“ here, and suggested not to go over 75%. The % can be defined variable with bg value
112 (like also in autoISF).

113 The Boost function automatically shuts off as soon as delta and the average deltas are aligned,
114 i.e. when the accelerated rise goes over into a constant rise (compare pp_ISF in autoISF).

115 However, the boost function is only „dormant“ if the boost window lasts longer for more meal-
116 related accelerations.

117 Additional functions are a step-count modified dynamic_ISF, inactivity detection etc

118 A couple of safety feature are integrated. The user can define an iob limit for boosts, here called
119 UAM Boost max IOB . In Preferences/Treatments There is also a user adjustable Low Glucose
120 Suspend threshold. This allows the user to set a value higher than the system would normally use,
121 such that when predictions drop below this level (65...100), a zero TBR is set.

122 For an example, see [case study 13.3](#)

123 More info: <https://discord.gg/nYC4T9PgCR> ; <https://github.com/tim2000s/no-bolus-dev>
124 ; https://github.com/tim2000s/Boost_AAPS_3.2/blob/Boost-Master-3.2/README.md
125 Contact: Tim Street @ diabettech.com

126

127 13.3.2 AIMI

128

129 AIMI has a single goal: to minimize the decisions necessary to maintain the target range, simplify
130 the composition of the profile for the user or doctor accompanying the patient, and allow the patient
131 to live normally without having to count carbohydrates or even without signifying physical activity
132 (especially for brisk walking).

133 A key component of AIMI concept is to give a **small pre-bolus before each meal** (“Meal
134 Announcement” that also provides some pos. iob).

135 • A **simplified profile** composition (neutral ISF around 100, DIA 9, target 90-90, a single value for
136 basal, a ratio that is not used in AIMI, so not important) For a first basal estimate, you can use the
137 TDD / weight ratio.

138 • Some variables in preferences that are important (AIMI_UAM which allows AIMI to make
139 decisions, Max SMB size which is the highest value for an SMB, B30_duration (which is the
140 duration during which the **basal will be forced after a manual bolus**), B30_upperBG and
141 B30_Upperdelta (these last two variables represent the conditions for replacing smb with a
142 consistent TBR depending on the delta)

143 • The basal profile is calculated by a polynomial equation.

144 • The ISF is calculated from the TDD (**dynamicISF**) and is adjusted based on the evolution of TIR
145 throughout the day and the **detection of physical activity**.

146 • The detection of glycemic rise (or the opposite situation) is also calculated by a polynomial
147 equation, which will influence the change of target but also the replacement of SMB by a TBR
148 between 100% and 500% or by an SMB of the same equivalence.

149 • SMB calculation is done in several ways specific to AIMI depending on the evolution of the delta
150 and IOB, with a distribution that can be done in three parts depending on the conditions.

151 Example scenario of execution, on almost all existing variants:

- 152 1. Make a "standard" manual bolus. I usually do 1.5U or 2U with luymjev
- 153 2. Just after this bolus, AIMI will force the 500% TBR for a duration defined by the user. The
154 observation made is that the absorption of insulin such as humalog for example is acceler-
155 ated and will strongly limit the first wave.
- 156 3. Depending on the options chosen, it is possible to receive an SMB of the initial manual bo-
157 lus size after the duration of the 500% TBR
- 158 4. Then the rest of the calculations will depend on the result of a polynomial equation and its
159 evolution.
- 160 5. A few hours later, if the patient decides to take a walk to go shopping, or other activities re-
161 quiring movement, the phone sensor will send information on the number **of steps taken**.
162 This will result in a reduction of the profile to about 60%. The return of the profile to normal
163 will be done in stages, in the first half hour following the activity, the profile will be restored
164 to about 80%.

165 The AIMI developer has been working on using machine learning (using tensorflow lite).

166 More info <https://discord.gg/7ehczAfZ>

167 The developer hasn't kept the code public. AIMI can only be obtained as an apk via joining their
168 WhatsApp group or here:

169 https://github.com/MTR93600/OpenApsAIMI/tree/dev_mergemilos_addOAPSAIMI

170 Given the very high number of changes happening in this AAPS variant, it is probably deemed
171 important to keep it in a tight sub-community. But, caution: This can be seen as violation of the
172 Open Source principle

173 Contact: Mathieu Tellier @ AndroidAPS User FB / Twitter @MTR93600/, Discord: MTR

174

175

176 13.3.3 EatingNow (EN)

177

178 This version of AAPS has evolved over time using elements from AIMI and Boost. It includes a
179 modified dynamicISF which moves ISF modulation in the direction as pioneered by autoISF, and
180 also uses Automations for FCL.

181

182 "Eating Now" (EN) allows user definable SMB's when deltas are sufficient and accelerating.

183 The intent of this plugin is the same, to deliver insulin earlier using mostly the AAPS predictions.

184

185 As all other variants for FCL, also EatingNow requires to set glucose TT occasionally, to nudge the
186 loop in certain direction, notably to announce and be prepared for activity.

187

188 Operating Modes provide 3 levels of „aggressiveness“ in 3 time windows:

- 189 • Master AAPS w/up to 120 min basal per SMB when EN is off (usually set for night-time).
- 190 • EN (usually set for daytime) is when the modified algorithm is capable of boosting ISF and
191 insulin delivery. At BG level rises within the EN Window, a „UAM maxBolus“ is given as a
192 first SMB. Recommended Setting: 1h current basal in units (max allowed: 2).
- 193 • ENW: A further boosted SMB will be issued in this ENW time window (e.g. for breakfast, or
194 generally for the first meal of a day, after fasting, with higher insulin need). Upon detection
195 of rising glucose, a SMB called Breakfast COB maxBolus is given by the loop. Recom-
196 mended Setting: 25% of average breakfast total units

197

198 EN uses the dynamicISF concept, modified to making ISF stronger as and eventualBG predictions.
199 Increase.

200

201 Specifically for the ENW (usually: breakfast window), an additional boost factor called Breakfast
202 ISF/CR Percentage (e.g. 125 or 150%) can be applied

203 A setting „TIRS“ provides a very simple version of autoISF (dura_ISF) and sharpens ISF
204 temporarily when bg „seems stuck“ above a certain value.

205

206 Autosens sensitivityRatio will be overridden by EN sensitivity options.

207

208 SMB delivery ratio for insulinReq. Is set to 65% for when EN is disabled (overnight, usually).

209 It is recommended to set maxSMBBasalMinutes and maxUAMSMBBasalMinutes to 30 minutes
210 max as these will be used when EN is OFF or in SLEEP mode. Falling back on OAPS SMB
211 settings is considered as the safe mode should you experience any issues with sensitivity or EN
212 settings in general

213 It is set 85% for an active ENW, or 75% when EN is on but ENW not active

214

215 Furthermore, SMB optionally can be disabled day/night below defined bg level/s (SMB Disabled)

216 More info <https://discord.gg/XqhnPRChEP> (method description in pinned post)

217 <https://github.com/dicko72/AAPS-EatingNow> scroll down to README.md

218 Contact: dicko via Discord channel

219

220

221

222 13.3.4 Tsunami

223

224 The Tsunami loop algorithm analyses blood glucose and insulin activity developments to estimate
225 bolus requirements during meals, without the necessity of carb announcements.

226

227 Users must make a **meal announcement via a button** on AAPS main screen. It switches on the
228 main Tsunami algorithm for a finite amount of time.

229

230 In between meals (when Tsunami is inactive), users are given the choice between running a
231 weaker version of the Tsunami algorithm (called wave), or falling back tooref1.

232

233 A “historic” merit of this method was that it pioneered a BG smoothing algorithm that later
234 became included as a plugin in AAPS.

235 The insulin models dynamically readjust DIA based on bolus size so that a user-set, fixed
236 DIA value is no longer needed.

237

238 For best results, it is recommended to issue a **bolus** at the beginning of a meal to account for the
239 disadvantageous kinetics of subcutaneously administered insulin in a UAM setting.

240 More info <https://discord.gg/veRKcgwVUT> GitHub repository: <https://github.com/piecycle/tsunami>
241 official documentation: [https://cdn.discordapp.com/attach-](https://cdn.discordapp.com/attachments/969948954949189633/972852790739238992/tsunami_guide_3_2.pdf)
242 [ments/969948954949189633/972852790739238992/tsunami_guide_3_2.pdf](https://cdn.discordapp.com/attachments/969948954949189633/972852790739238992/tsunami_guide_3_2.pdf)

243 Contact: nichi#1391 on discord / piecycle on GitHub

244

245 13.4 No-Bolus Looping with Carb Entries

246

247 Some oreof(1) loopers attempting to go full closed loop reported that they do best when they (do not
248 bolus but) give their loop precise carb (and absorption time) information. This:

249 * announces a meal to follow (so it is not UAM, but might be called full closed looping if the
250 insulin management is left 100% to the loop)

251 * provides data on cob, and with the glucose and insulin activity info the loop has, it can
252 always calculate how much more carbs are to become absorbed (to the extent the carb-
253 related infos the user put in is correct)

254 * will display realistic cob info to the user, including cob info looking forward (rather than
255 only calculating carb deviations for the past minutes or hours, and making some coarse
256 assumptions for the upcoming hour). It gives the user better feeling of safety if she/he can
257 see cob info in addition to the available iob info, and insulin activity prediction.

258

259 With detailed carb (amounts + absorption times) inputs, the loop has best-possible info to provide
260 „the best expert fit“ of insulin activity and carb absorption.

261 It still rarely can come close to physiological values, because the time-delays inherent in
262 our „artificial pancreas“, notably the stretched out DIA, make it difficult still, compared to a
263 real pancreas.

264

265 So, carb inputs could help. However,

266 • only to the extent amounts and time pattern for absorption („eCarbs“) are correct ((which,
267 every day, is a mission impossible))

268 • the oref(1) loop still largely „waits for glucose to rise“, and there is no significant time ad-
269 vantage from inputting carb info

270 Only the **user**-bolussing for **expected** carb absorption in hybrid closed loop offers a
271 convincing time advantage (but with associated risks).

272 • inputs require actually more attention to detail than it is good practice even in AndroidAPS
273 hybrid closed loop, so in that respect a step back, not forward.

274 Entering **precise** carb information takes away a very large part of the attractivity of full closed
275 looping.

276 And entering *imprecise* carb info could easy be inferior to not doing *any* carb inputs = to letting the
277 UAM mode of oref(1) figure out further carbs that probably come to be absorbed in the next
278 minutes, judging from the pattern of the calculated past *carb deviations* (see [section 4.5](https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand-determine-basal.html#understanding-the-basic-logic-written-version) and
279 [https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand](https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand-determine-basal.html#understanding-the-basic-logic-written-version)
280 [-determine-basal.html#understanding-the-basic-logic-written-version](https://openaps.readthedocs.io/en/latest/docs/While%20You%20Wait%20For%20Gear/Understand-determine-basal.html#understanding-the-basic-logic-written-version)).

281 PS: Because that is so, also loopers who do carb inputs get the UAM predictions besides their
282 other predictions, and their algo makes a judgement (every 5 minutes) as to what the best
283 calculation might be for where glucose, underlying „real“ carb absorption, and estimated carb
284 deviation are headed.

285

286 13.5 Machine Learning

287

288 Involving machine learning (“artificial intelligence”) could help both in the learning/tuning phase, but
289 also in fine adjustments in daily utilization.

290 The study that was already referenced in [section 1](#) discusses on page 80 the application of
291 machine learning in some predictions of postprandial glucose response (IEEE Control Systems
292 Magazine, ResearchGate: The Artificial Pancreas and Meal Control. A. El Fathi et al, IEEE Control
293 Systems Magazine Feb.2018 p.67-85.).

294 So there is already a body of data and evidence. To which extent it lends itself to UAM remains to
295 be researched. For this, a body of data would have to be captured from UAM loopers, and I fear
296 many more data would be required than what could easily be captured in Clarity® or even in the
297 OPEN project database.

298

299 In the DIY universe, a prototype solution was already developed for AIMI ([section 13.3.3](#)).

300

301 We might see industry come up with a 1st generation solution that will probably be geared to folks
302 with miserable HbA1c and poor carb counting/meal handling, to offer a safe gradual improvement.

303

304 A top performing entirely self-learning system might be impossible to design:

305 For instance, if today you do something entirely different from yesterday (don't we all want this
306 freedom – even need it? Think about the fasting day following a feasting day...) there are
307 two problems:

- 308 • Such systems rely on information from the preceding day, or an average of several preced-
309 ing days
- 310 • The user does not know/learn much about how the system works, what it is calibrated for
311 today, how she/he might intelligently change something for the specific different situation
312 coming up. This seems like the opposite of the FCL solutions we discussed, for instance
313 self-defined Automations, combined with profile switches for to-be-expected temporary sen-
314 sitivity shifts (see [section 13.1](#) , or the more sophisticated options presented for exercise in
315 [section 6.1.3](#) and [6.2](#)).

316

317 13.6 Dual Hormone Systems

318

319 Many see a dual hormone “double full loop” as the ultimate system.

320

321 The beauty of this concept would be that the second pump could influence the glucose curve via
322 giving glucagon or an analogue, thus overcoming the strongest limitation our current systems
323 have:

324 Taking basal away (zero-temping) is only a severely limited course of action against impending
325 hypoglycemia, and therefore, to keep things safe at the back-end of each meal, fighting glucose
326 highs is more limited than we would like to see.

327

328 In conclusion, the glucagon component **not only helps stay out of hypos. It enables a more**
329 **aggressive treatment for preventing, or reducing, high glucose values, as well.**

330

331 While insulin and carbs have complex activity curves stretching over hours, glucagon has a
332 window of physiological activity starting 5-10 minutes after administration, and lasting only 30-40
333 minutes. Compared to insulin and carbs, that makes it a better component for rapid corrections
334 (without a lengthy “tail” of action).

335

336 As glucagon does not per se introduce more calories, but stimulates glucose release from the liver,
337 there should at least be no concern about gaining body weight from eventual roller-coasters the
338 dual loop might send us into. Actually there could be a nice side benefit of helping in body weight
339 control. Also, activity/sports management could become as easy as the meal management became
340 in the UAM step into full closed looping.

341

342 It will be interesting to see for which application(s) the dual loop will be developed and launched;
343 as part of a full closed loop with top performance, or as part of even only a hybrid closed loop for
344 problem patients?

345

346 It remains to be seen how well such systems work in day-to-day circumstances. And whether “real
347 people” will be able to handle all the involved technology, and use it in ways that truly could justify
348 the substantial extra cost.

349

350 The author currently is not really looking forward to become loaded with even more technology,
351 and quite happy with an aggressively tuned full UAM closed loop (...and an occasional nice post-
352 dinner or during- activity snack).

353

354 However, the dual hormone path holds enough promise to learn more about it, and to test it some
355 time in the near future.

356

357 This is an exciting time to be part of the open source T1D community. Anyone is welcome to
358 contribute ideas, help develop software or instructions how to use. Carefully weigh for yourself
359 what may be your entry point for eventually surmounting the initial hurdles, and **JUST EAT happily**
360 **ever after.**

361