

Case Study 6.2: Biking day with high carb lunch

V.2.3



This case illustrates how an entire exercise day can be managed by the loop nearly automatically, also in Full Closed Loop, when using the exercise mode plus some Automations (as outlined in FCL e-book [section 6.5.2](#)).

Method

Intended exercise and targeted sensitivity modulation

The meal challenge

Automations installed for meal management with exercise following

Overall result

Preserving settings; DIY cockpit

Addendum: Logfile case analysis

Method

Full Closed Loop with dev variant of AAPS 3.2.0.2 w/autolSF 3.0.

Lyumjev 100 (DIA 7h) in Combo pump w/ 10mm Teflon cannula (0-48h)

2 x G6 overlapping; xDrip; no smoothing in AAPS

profile basal ~ 14 U (0.41...0.75 U/h); profile_ISF 36...44 mg/dl/U (circadian);

TDD 37 U

Main settings for 24/7 adaptation of ISF *):

SMB range extention and autolSFmax both = 2.9; SMB delivery ratio = 0.75 fixed

bgAccel_ISF_weight = 0.24; pp_ISF_weight = 0.03; dura_ISF_weight 0.8

iobTH_percent=60

***) Caution: Do not copy settings** from others, not even for starting your tuning. Why, see FCL e-book [section 4.1](#).

Automations active: see under “Automations...”

User interference:

no boli ; no carb inputs ; no setting Eating Soon TT

Setting **% profile** and **TT** in top row of the AAPS main screen “to announce” exercise, before starting (just like done in Hybrid Closed Loop for exercise, too).

Because the case is complicated by also including a major meal in FCL, **one** additional (but time-uncritical) **manual step** is needed (that will start a “cascade of” Automations).

Occasionally *) having an eye on bg and on insulin activity curves = A few times, taking a small snack while biking

*) The curves, notably prediction curves, on the AAPS main screen (mounted on the bike handlebar) allow to judge at one glance, for how much longer there is sure safety to not need extra carbs against going low - and when, latest, to look again.

Intended exercise, and targeted sensitivity modulation

Leisure-time **biking** for about **4 hours** in total, with a sit-down longer lunch break.

I had only a minimal breakfast in order to keep **iob low when starting** to bike.

Rather I decided to carry along (in reach while biking) two protein bars (only 1 used, well before lunch, basically my breakfast replacement), and a couple of 2 grams glucose tablets.

*From past experience, **70% profile** suits well, set already at least 2 hours ahead of starting to bike, and continuing (or slightly elevated to 80%) also in the night that follows.*

Setting a 70% profile modulates my (average) 0.55 U profile basal to $0.7 * 0.55 = 0.39$ U (see Addendum, emulator-line 20)

*In the past, **glucose target** had been elevated to somewhere in the 120 to 170 range, out of a “gut feeling” looking at many factors, like current bg state, intensity of exercise, intended snacking.*

*This time, my setting had to consider the implications from the set TT on modulations that the (fairly new to me) exercise mode will do: In the **Exercise mode**, a **125 mg/dl TT** set translates ...*

... following a rather complicated calculation scheme, see autoISF Quick Guide p.7-9 at: <https://github.com/ga-zelle/autoISF..>

...into another (= getting multiplied) basal **reduction of 33 %**, to 67% of (70% of $0.55=0.39$ U= 0.26 U overall resulted (see Addendum, emulator-line 19)

Taking %profile, Exercise mode, and TT setting together, the loop operates at a **modulated sensitivity** of only **47%** of its profile value (or of any autoISF adaptations based on it).

*Note that, also for **other types of exercise**, I can, over time, establish sets of suitable % profile reduction and of TTs, that would result in any desired different modulated % sensitivity.*

Regarding preserving such settings for future use, see extra subchapter further below

The meal challenge

Duck with French Fries and a salad bowl that included fruit.

No carbs were entered into AAPS, no bolus was given by the loop

Going into a high carb meal with only 47% “power” seems no good idea, regarding the to-be-expected bg spike, when not bolussing at all (= that precedes exercise must be managed in FCL mode with these cornerstones in mind:

1. Before lunch time, ending the exercise TT and going to (the much lower) profile bg target provides
 - more iob
 - opportunity that an Automation can become activated (at bg delta > 10 mg/dl) that temporarily further lowers the TT to 74 mg/dl, which **maximizes first SMBs for the intended high carb lunch**
 - Overall, in this initial meal stage, the loop acts as aggressive “against” the high carb meal, as it would *outside of* an exercise context.
2. As exercise will (again) follow, it is here **very time critical** to end the 74 mg/dl TT again, so the exercise mode with its **125 mg/dl TT** can be immediately **re-installed**. That is of high importance, because the exercise mode also results in a **dynamic_iobTH** which, like basal, also goes lower and ends the aggressive SMB fire against the high carb meal earlier. That facilitates continuing the exercise at an **elevated bg** level, and with **significantly reduced iob** (only ~47% of otherwise “allowed” level, in our case).

Note, starting exercise with low iob level was already our receipt for success in the morning.

Rather than tightly watching the AAPS screen during lunch, to **manually** enact this important switching from aggressive meal management towards pre-exercise slow-down, in exactly the right point of time (before iob shoots too high), *I had* **automated this procedure** with a cascade of 3 Automations:


Automations for fully automatic meal management with exercise following

I used a sequence of 3 Automations, of which only the first one had to be manually triggered, in just one time-uncritical key stroke from the AAPS home screen.

The others come on automatically when the respective Conditions are met.

Automation 1

The key first task was, to approach a meal that precedes exercise with full loop aggressiveness, but to make sure that this aggressiveness stops immediately after a (reduced) iobTH is exceeded. The reduced iobTH ensures that not too much insulin is on board for exercise after the meal. Also it provides an elevated bg level at (re-)start of exercise.



„User action“ is always ticked-on

This will, in the defined time space *) ..

..offer the “DIY cockpit” button..

..which I must press any time (~90...30 minutes) before my lunch.



In this Automation, the box “User action” should be permanent. This will, in the defined time slot as specified under Conditions, automatically provide a grey button on the bottom of my AAPS home screen (DIY cockpit”) that I can freely name (= headline of my Automation).

To keep the AAPS home screen as clean (and relevant) as can be, that button will show only in the time slot as specified under Conditions.

*) It will be reduced to something realistic. Only for development and testing purposes it had to extend into the night.



Automation 2

When (my in Automation 1, for exercise following the meal, to 35 % reduced iobTH that translates for me into) $iob > 3.8$ is exceeded, I want two things:

(1) The loop shall now automatically run milder, on my reduced exercise %profile

(70%)(after the meal rise had been managed based on 100% profile, boosted by bgAccel_ISF driven full loop aggressiveness).

(2) I like also to set a exercise TT. This, however, is not possible. I first have to force an end to my EatingSoonTT of 84:



Automation 3

Now, that Automation 2 ended my TT of 84, Automation 3 can (max 5 minutes later) set the desired exercise TT=125 (which implies the exercise mode as in detail explained before).

Note that Automations 2 and 3 are fully automatic, no User Action involved.

139


140
141
142



143


144

145



< 4.4 U

The lunch iob hump in middle of the bottom graph (of the chart above, left) shows that the ~ 4 U iobTH was preventing higher iob as would be normal for a big lunch.



This chart from an emulator based analysis shows nicely how iobTH (dotted blue lines) changes with the TT set (red line, and scale with red numbers on the right):

Around 12:00 h TT 74 expired and profile target briefly was set, followed within a couple of minutes by switching back (as before 10:15 h in the morning biking time) to the 125 mg/dl exercise TT. This state allows only about 2 U as iobTH (see kink in the blue dotted line at 12:10 h), and my iob (blue curve) came nicely down to that low level, as desirable for sports - but sure unusual so shortly after eating a substantial meal.

Not astonishingly, few SMBs are in the picture (fat green line mostly at zero) and the loop can handle carbs in an exercise dominated period nearly by just elevating basal (%TBR up to 400%). (Actually, the very first graph shows for around 13:30 glucose from the meal had

risen to about 140 mg/dl and one additional SMB was necessary, and was also possible because from around 12:45 the iob (blue line) was below the (dotted blue line: lowered) iobTH again, as this graph shows.

Summary

Using the **top button row** on the AAPS main screen, I just **input** the **70%** (for the whole day) and **125 mg/dl** (for the first hours), which goes super fast and easy. It will immediately turn

- from all three fields grey,
- to „70%..“ on the profile field-turned-yellow, „125 .. „ in the TT field-turned-yellow, and also the exercise field lit yellow in the middle.

So, very easy to see on one glance, I am in the exercise mode, and which are the key settings (see picture with „95“ glucose).

Very roughly an hour before lunch time I activated the **User Action Automation #1**, reverting to profile target.

In just seconds I was able to „frame“ the upcoming exceptional situation for my loop, to manage me, from then on, fully automatically through the day!

And in case I would want to temporarily modify a parameter, or to prematurely exit: Same easy procedure, just within 1-2 seconds, right from my AAPS home screen.

The example demonstrated that using the exercise mode softened the loop response appropriately, while some Automations greatly helped to still get a major meal sufficiently treated under FCL conditions.

Preserving settings in **your** DIY cockpit

As I did in my initial set-up and tuning, you must first „learn“ from **your** past looping data (FCL and HCL), where about %profile and set TT should lie, for the kind of exercise that you occasionally were and are up to.

But, no need to make a science out of it. Unless you are competing in professional sports, it should be good enough to go by gut feeling, and by experience („what setting should I slightly alter the next time?“).

192 As the loop re-calculates every 5 minutes, it is *not* important to get things *exactly* right.
193 Automatic adjustments (every 5 minutes) allow the loop to still keep things under good-enough
194 control.

195 The time windows for doing the profile switch, and for setting a suitable TT can differ (and
196 they can also be automated, so not to require multiple inputs over the course of your
197 exercise day).

198 **Using all available tools allows a nearly surgical approach to what you want to achieve**
199 **for your favourite type(s) of exercise.**

200
201 Observe that your general setting (in AAPS Preferences) for **half-basal exercise target**
202 strongly influences where you land in your attempts to assign feasible TTs for various
203 exercise types. (More see autoISF Quick Guide p.7-9 at: <https://github.com/ga-zelle/autoISF>
204 or FCL e-book [section 6.1.3](#)).

205 Define one **User Action Automation** for each type of exercise, and name it as you would
206 like to see the according button labelled on your AAPS main screen.

207 If needed, define automatically triggered adjunct/tandem Automations (like we needed in our
208 very special case with a major meal inside the exercise activities)

209 To limit the number of buttons in your “cockpit”:

- 210 • Assign time spans in the Conditions (when, in any 24 h period, can this exercise happen?)
- 211 • Completely dis-activate (“shelve”) those Automations, that are for types of exercise that do not frequently
212 occur

213
214 **Caution:** As mentioned in other places, Automations can be tricky as to whether they actually
215 will ever work, because the loop goes through the exact **sequence of all your active**
216 **Automations**, and might be switched into a direction that no longer is compatible with the
217 conditions that must be a given, for the Automation you think that should kick in.

218

219

220

221

Addendum: Logfile analysis with the emulator

Skip this last section, unless you like to learn more about using the emulator.

Analysis of my biking day with the emulator required to load the logfiles from my phone into the PC, and to have the files for the emulator downloaded from Github.

- 1) Load an empty vdf file, and access to logfiles (erasing the end with UTZ time, and putting an asterisks after date_ => all of that date get loaded into the mask
- 2) Define the time window of interest, using UTZ (= MEZ summertime minus 2) in the last two lines (right hand side input fields) of the form:

Manage Inputs and Outputs for Emulating AAPS Settings

Your working directory
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autosf_UAMcodechanges/13_autosf_3.0/Sports day biking

Select Inputs | Select Graphics Options | Execute the Analysis | Inspect Results

Your variant definition file
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autosf_UAMcodechanges/13_autosf_3.0/Oct.12_worst.day/noChange-Bernd.vdf [Browse] [Edit]

Your AAPS logfile(s)
/Documents/B_Pump, CGM, Phone, Looping/Looping/autosf_UAMcodechanges/13_autosf_3.0/2023_Analyse_aISF3.0_SepOct/Logfiles_Sep22-Oct.22/AndroidAPS_2023-10-04_*.zip [Browse] [Show matches]

example date/time format ... 2019-11-06T12:30:00Z
2023-10-04T10:00:00Z
2023-10-04T16:00:00Z

☒ Use start time by entering UTC date/time
☒ Use final time by entering UTC date/time

- 3) Execute Analysis / Run Emulation yields results as table:

Select Inputs		Select Graphics Options				Execute the Analysis				Inspect Results																			
Messages from Emulation														Clear Messages				Run Emulation				Emulation finished ..							
--5% range--		--lin.fit--		---parabola		fit----		-----ISF factors-----						-----ISFs-----				insulin Req		---SMB---		--tmpBasal--							
B	dura	avg.	dura	rate	dura	last-Δ	next-Δ	auto	acce	bg	pp	delta	dura	orig	prof	emul	orig	emul	orig	emul	orig	emul							
0	10	103.3	20.0	2.3	15	-3	-8	1.0	-0.1	0.96	1	1	1	229.2	61.4	229.2	0	0	0	0	0								
0	15	103.3	25.0	1.7	25	-1.14	-2.57	1.0	1	0.96	1	1	1	95.7	61.4	95.7	0	0	0	0	0								
0	20	103.4	10.0	1.0	15	1.8	3.8	1.0	1	0.96	1	1	1	67	43	67.0	0	0	0	0	0								
0	25	104	10.0	2.0	15	2.6	3.6	1.0	1.1	0.96	1	1	1	60.7	43	60.7	0.1	0.1	0	0	0.57								
0	35	105	25.0	1.5	25	1.5	1.53	1.0	1.0	1.0	1.04	1	1.08	39.9	43	39.9	0.88	0.88	0.7	0.7	2.31								
0	0	117	30.0	2.2	15	7.2	11.2	1.0	1.78	1.01	1.16	1	1	24.2	43	24.2	3.02	3.02	2.4	2.4	5.5								
0	5	123.5	35.0	2.9	15	0.4	-4.1	1.0	0.55	1.02	1	1	1	76.4	43	76.4	0	0	0	0	0								
0	15	121.5	40.0	2.7	20	-2.37	-6.09	1.0	0.65	1.01	1	1	1.07	61.8	43	61.8	0	0	0	0	0								
0	0	116	70.0	1.8	25	-5.89	-9.53	1.0	0.65	1.01	1	1	1	65.3	43	65.3	0	0	0	0	0								
0	5	115.5	15.0	-3.0	20	-3.57	-4.28	1.0	0.98	1.01	1	1	1	43.4	43	43.4	0	0	0	0	0								
0	10	114.7	10.0	-1.5	15	-0.8	1.2	1.0	1.26	1.01	1	1	1.04	32.4	41	32.4	0	0	0	0	0								
0	15	115.3	25.0	-1.7	20	2.77	5.49	1.0	1.26	1.01	1.08	1	1.06	32.4	41	32.4	0	0	0	0	0								
0	20	115.8	10.0	2.5	25	3.39	5.35	1.0	1	1.01	1.02	1	1.08	38.1	41	38.1	0	0	0	0	0								
0	10	119.3	15.0	3.1	15	3.6	4.1	1.0	1.06	1.01	1.1	1	1.04	37.3	41	37.3	0.29	0.29	0	0.2	1.13								
0	0	130	10.0	6.0	15	7.4	10.4	1.0	1.65	1.02	1.14	1	1	24.9	41	24.9	2.89	2.29	0	1.8	5.5								
0	0	141	10.0	9.1	20	10.7	13.84	1.0	1.69	1.05	1.22	1	1	22.6	41	24.3	6.46	5.97	1.8	2.2	5.5								
0	0	141	10.0	9.1	20	10.7	13.84	1.0	1.69	1.05	1.22	1	1	22.6	41	24.3	6.46	6.01	1.8	2.2	5.5								
0	0	141	5.0	11.4	20	10.7	13.84	1.0	1.69	1.05	1.22	1	1	24.3	41	24.3	0	0	0	0	0								
0	0	149	5.0	8.1	20.1	10.03	11.39	1.0	1.28	1.06	1.16	1	1	32	41	32.0	0.69	0.69	0	0.4	1.99								
0	5	149	25.0	7.1	15	1.05	-4.39	1.0	0.38	1.06	1	1	1	102.2	41	102.2	0	0	0	0	0								
0	10	148.7	40.0	5.4	20	-2.41	-6.96	1.0	0.5	1.06	1	1	1.13	72.1	41	72.1	0	0	0	0	0								
0	20	146.6	10.0	-1.5	14.9	-2	-3.01	1.0	0.88	1.06	1	1	1.26	37	41	37.0	0	0	0	0	0								
0	25	145.7	10.0	-3.5	15	-4.6	-6.6	1.0	0.77	1.06	1	1	1.32	37.3	41	37.3	0	0	0	0	0								
0	5	138.5	10.0	-5.0	24.9	-5.53	-7	1.0	0.83	1.04	1	1	1	47.7	41	47.7	0	0	0	0	0								
0	10	136.7	15.0	-4.4	15	-3.4	-2.4	1.0	1.22	1.03	1	1	1.07	32.9	40	32.9	0	0	0	0	0								
0	10	136.7	5.0	-3.0	15	-3.4	-2.4	1.0	0.88	1.0	1	1	1	33.1	40	45.3	0	0	0	0	0								
0	15	136	5.0	1.0	15	0.6	3.6	1.0	1.63	0.99	1	1	1	49.1	50	49.1	0	0	0	0	0								
0	15	133.3	35.0	2.0	20	1.64	1.07	1.0	0.65	0.96	1	1	1	74.4	50	82.0	0	0	0	0	0								

In line 6 (12:34 CET) a 2.4 U big SMB is given, driven by bgAcel_ISF factor of 1.78

Highlighted is 13:29 CET (for some reason there is duplicated line for it). Here, a big insulinReq of 6,5 U is cut down to a SMB of only 1.8 U..

- 4) ... to see how our settings, TT and and sports button worked here, and what iobTH applied, we must look into the logs available from: ((alternatively, the SMB tab from exactly those 5 minutes would show, too))

Instead of Execute Analysis, press Inspect results...

Manage Inputs and Outputs for Emulating AAPS Settings

Your working directory
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Sports day biking

Select Inputs Select Graphics Options Execute the Analysis Inspect Results

Your variant definition file
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/noChange-Bernd.vdf Browse Edit

Your AAPS logfile(s)
/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/2023_Analyse_aISF3.0_SepOct/Logfiles_Sep22-Oct.22/AndroidAPS_2023-10-04_*.zip Browse Show matches

example date/time format ... 2019-11-06T12:30:00Z
2023-10-04T10:00:00Z
2023-10-04T16:00:00Z

☒ Use start time by entering UTC date/time
☒ Use final time by entering UTC date/time

Manage Inputs and Outputs for Emulating AAPS Settings

Your working directory
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Sports day biking

Select Inputs Select Graphics Options Execute the Analysis Inspect Results

*.log - Your file showing edits from the variant assignments
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.noChange-Bernd.log Browse Show

*.csv - Your table comparing key values of original vs emulation
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.noChange-Bernd.csv Browse Show

*.delta - Your table comparing bg deltas of original vs emulation
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.noChange-Bernd.delta Browse Show

*.orig.txt - Your short log of original analysis
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.orig.txt Browse Show

*.txt - Your short log of emulated analysis
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.noChange-Bernd.txt Browse Show

*.pdf etc. - Your graphic file comparing key values of original vs emulation
C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.noChange-Bernd.pdf Browse Show

If you press the 4th option you get the logs which include basically all SMB tab info in an exhaustively long list.

C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/Oct.12_worst.day/2023-10-04T10-00-00Z.orig.txt - Notepad++

Datei Bearbeiten Suchen Ansicht Kodierung Sprachen Einstellungen Werkzeuge Makro Ausführen Erweiterungen Fenster ?

AndroidAPS_2023-10-12_00-00-03_4.zip.orig.txt AndroidAPS_2023-10-12_00-00-03_4 AndroidAPS_2023-10-12_00-00-03_4.zip.noChange-Bernd.csv AndroidAPS_2023-10-12_00-00-03_4

```

16 ===== loop in row.13418 ===== of logfile C:/Users/Bernd/Documents/B_Pump, CGM, Phone,
17 Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/2023_Analyse_aISF3.0_SepOct/Logfiles_Sep22-Oct.22/AndroidAPS_2
18 ...created at= 2023-10-04T10:04:16.210ZCRIF
19 ----- Script Debug -----CRIF
20 Sensitivity ratio set to 0.67 based on temp target of 125;CRIF
21 Adjusting basal from 0.39 to 0.26130000000000003;CRIF
22 ISF from 61.4 to 91.7CRIF
23 CR: 10.714285714285715CRIF
24 -----CRIF
25 start autoISF 3.0CRIF
26 -----CRIF
27 SMB disabled; TempTarget 125 is odd numberCRIF
28 Loop at minimum powerCRIF
29 effective meal add on is 0CRIF
30 acce ISF adaptation is -0.1CRIF
31 bg ISF adaptation is 0.96CRIF
32 weakest autoISF factor -0.1 limited by autoISF min 0.4CRIF
33 final ISF factor is 0.27 including exercise mode impactCRIF
34 -----CRIF
35 end autoISFCRIF
36 -----CRIF
37 current temp: 0.27 lastTempAge: 0 m tempModulus: 25 mCRIF
38 profile.sens: 61.42857142857143.sens: 229.2.CSF: 21.39CRIF
39 Carb Impact: -8.3 mg/dL per 5m; CI Duration: 0 hours; remaining CI (~2h peak): 0 mg/dL per 5mCRIF
40 UAM Impact: -8.3 mg/dL per 5m; UAM Duration: 0 hoursCRIF
41 minPredBG: 130 minIOBPredBG: 130 minZTGuardBG: 107CRIF
42 minUAMPredBG: 113CRIF
43 avgPredBG: 154.COB: 0 / 0CRIF
44 BG projected to remain above 125 for 0 minutesCRIF
45 BG projected to remain above 83 for 240 minutesCRIF
46 naive_eventualBG: 128.bgUndershoot: -45.zeroTempDuration: 240.zeroTempEffect: 358.carbsReq: -19CRIF
47 ----- Reason -----CRIF
48 COB: 0.Dev: -25.BGI: 5.ISF: 229.CR: 10.71.Target: 125.minPredBG 130.minGuardBG 97.IOBPredBG 154.Eventual BG 1
49 AAPS scan from AAPS Logfile for SMB comparison created on Fri, 27 Oct 2023 18:09:40 +0200CRIF
50 FILE=C:/Users/Bernd/Documents/B_Pump, CGM, Phone, Looping/Looping/autoISF_UAMcodechanges/13_autoISF_3.0/2023_Analyse

```

Suchen Suchen Ersetzen In Dateien suchen


Suchen nach: autoISF

☐ Rückwärts suchen
☒ Nur ganze Wörter suchen
☐ Groß-/Kleinschreibung beachten
☒ Am Ende von vorne beginnen

Suchmodus
☒ Normal
☐ Erweitert (n, y, t, v, x...)
☐ Reguläre Ausdrücke ☐ findet \

The blue highlighted line shows we are at 10:04 UTZ (12:04 CET) . We can see how sensitivity ratio and also basal were lowered due to sports setting with a 125 mg/dl TT, or 129 used a bit later, which also disabled SMBs.

By using the search function you can jump, in that long list, to all places that e.g. have „autoISF“ in it, or „script debug“, or „SMB disabled“ (if you want to analyze when that happened...)



- 5) Now lets look at 11:29 UTZ (13:29 CET) which was when the first relevant size SMB was issued:

In the following, I copy and shorten texts from the log table as above, for some time points of high interest:

11:29 UTZ (13:29 CET), when SMB was delivered:

loop in row 14007created at= 2023-10-04T11:29:11.005Z

----- Script Debug -----

Sensitivity ratio set to 1.1 based on temp target of 74;

Adjusting basal from 0.55 to 0.6050000000000001;

ISF from 41 to 37.3

CR: 7.5

So, we see the meal TT of 74 (set by one of my Automations for cases in which I did not bother to set an EatingSoonTT) temp. overrides the exercise target and sharpens the applicable ISF: ((It also elevates basal ; note that TBRs can run several 100% of basal))

start autoISF 3.0

SMB enabled; TempTarget 74 is even number

Loop at full power

```

272         acce_ISF adaptation is 1.81
273         bg_ISF adaptation is 1.05
274         pp_ISF adaptation is 1.22
275         dura_ISF by-passed; bg is only 0m at level 141
276         final ISF factor is 1.81
277         -----
278         end autoISF
279         -----
280         profile.sens: 41 sens: 22.6 CSF: 3.01
281         minPredBG: 220 minIOBPredBG: 135 minZTGuardBG: 94 minUAMPredBG: 231
282         avgPredBG: 220 COB: 0 / 0
283         BG projected to remain above 74 for 240 minutes
284         IOB 2.033
285         Full loop capped SMB at 1.87 to not exceed 130 % of effective iobTH
286     Here we see that our, for the exercise day significantly lowered, iobTH not only blocked SMBs above
287     that lower threshold. Also the 30% “the last” SMB is max. allowed to go beyond the valid iobTH now
288     amounts to significantly less.

289     So, while, as intended for start of a high carb meal, my loop is at full power (line 209), the last SMB
290     (line 223) got capped to way below insulinReq (line 231).

291         ... maxBolus: 3.2
292         ----- Reason -----
293         Eventual BG 304 >= 74, insulinReq 6.46. Microbolusing 1.8U. adj. req. rate: 13.53 to
294         maxSafeBasal: 5.5, temp 2.75 < 5.5U/hr.

295
296     6) A look at the cake after 14:20 (16:20 CET):™
297         From 14:09 UTZ (16:09 CET) on, I discontinued the sports TT and the loop reverted to 90
298         mg/dl profile target.

299         As desired for starting to eat cake, this made SMBs possible, and the only reduction of FCL
300         aggressiveness was the 70% due to temp. profile set for this sports day.

301         However, 14:14 - .19 - .24 and :29 there was no insulinRequired yet (while ISF factor
302         gradually ramped up from 0.52 to 1,09 (line 1735, 1772, 1807) .

303         At 14:34 UTZ (16:34 CET) a first cake related SMB of 0.5 U was issued based on 0,66 U
304         ins.Requ (line 1897), and based on a bgAccel_ISF factor of 1,54 (and final ISF factor of 1,52).
305

```

306 As there is a **>10 mg/d delta** (glucose rise), **my Automation** kicks in and sets for the next
307 loop decisions (for 26 minutes is as my Automation defines it) a **TT=74** mg/dl which makes
308 the loop more aggressive from 14:37 on:

```
309         created at= 2023-10-04T14:37:15.371Z
310         ----- Script Debug -----
311         Sensitivity ratio set to 1.1 based on temp target of 74;
312         Adjusting basal from 0.48 to 0.528;
313         ISF from 39 to 35.5
314         CR: 6.6
315         -----
316         start autoISF 3.0
317         -----
318         SMB enabled; TempTarget 74 is even number
319         Loop at full power
320         acce_ISF adaptation is 1.64
321         bg_ISF adaptation is 1
322         pp_ISF adaptation is 1.08
323         dura_ISF adaptation is 1.11 because ISF 35.5 did not do it for 30 m
324         final ISF factor is 1.64
325         -----
326         end autoISF
327         -----
328         profile.sens: 39 sens: 23.8
329         avgPredBG: 94 , BG projected to remain above 74 for 240 minutes
330         IOB 0.438   Eventual BG 100 >= 74, insulinReq 0.84. Microbolusing 0.6U.
331
332         created at= 2023-10-04T14:39:16.537Z
333         ----- Script Debug -----
334         Sensitivity ratio set to 1.1 based on temp target of 74;
335         Adjusting basal from 0.48 to 0.528;
336         ISF from 39 to 35.5
337         CR: 6.6
338         -----
339         start autoISF 3.0
340         -----
341         SMB enabled; TempTarget 74 is even number
342         Loop at full power
```

343 effective meal add on is 0.04
344 acce_ISF adaptation is 1.25
345 bg_ISF adaptation is 1
346 pp_ISF adaptation is 1.06
347 dura_ISF by-passed; bg is only 5m at level 97.5
348 final ISF factor is 1.25
349 -----
350 end autoISF
351 -----
352 profile.sens: 39 sens: 31.2
353 avgPredBG: 78 BG projected to remain above 74 for 240 minutes
354 IOB 1.114
355 Eventual BG 84 >= 74, insulinReq 0.13; setting 30m low temp of 0.04U/h. Microbolusing 0.1U.