

Case Study 4.2: Low carb lunch

23 / 24.07.24

V.1.0



This paper is about putting the **settings** I had (about 2 years ago) established **for my entire usual meal spectrum** to the test, with two extremely low-carb lunches.

Method

FCL (no carb inputs, no user boli) with dev variant of AAPS 3.2.0.4 w/ autoISF 3.0.1:
Lyumjev 100 (DIA 7h) in Combo pump w/ 10mm Teflon cannula (0-48h)
2 x G6 overlapping (see case study 1.5; sensors used ~ d3 – d15; xDrip, no smoothing in AAPS)
TDD ~ 35 U; profile basal ~ 14 U (0.41...0.75 U/h); profile_ISF 36...44 mg/dl/U; iobTH% = 0.6
Key settings for entire meal spectrum (~ 20 ... 90 g carb per meal):
SMB size limited at ~ 3.5 U (=2.9 x 120 minutes basal)
autoISFmax = 2.9; SMB delivery ratio = 0.75 fixed
bgAccel_ISF_weight = 0.22; break_weight 0.12; lower_ISF-range_weight 0.7; higher_ISF-range_weight 0.1; pp_ISF_weight = 0.03; dura_ISF_weight 0.8

Fish + veggie lunches managed by autoISF FCL

Settings that were better proven to work for fairly high carb meals (see e.g. Xmas [case study 4.3](#)) were put to the test with two (for “my spectrum”) extremely low carb meals.


Two very similar meals were set up for two consecutive days, with an option to lower, on the second day, the iobTH or any of the ..._ISF_weights.

As the FCL e-book suggests (in [section 4](#)), it would be easy to “nudge” my FCL for an “outlier” of comparatively low insulin need

- by lowering iobTH% either directly, or indirectly (via low %profile, high TT + exercise button)
- by elevating the effective ISF via lower bgAccel_ISF_weight and / or pp_ISF_weight, and/or via a lower % profile



38 The first example shown was a cod fillet plus broccoli with herbal crème cheese (12:10)



39

40 As that went well without any manual intervention...

41 other than, 2.5 -3 hours after the meal, 2 small cookies before/after a dogwalk

42 ...no settings changes were made for the second meal, which had about same carbs and
43 calories.

44 Note that iob stayed under 4U all the time, and under 2U in first 2 hours after meal start. That
45 is well below my iobTH of 6 U.

46

47

48

49

50

51

52

53

54

55


56

57

The 2nd meal was (on 24.July'24)

2 small tuna steaks, and black beans (12:50).

Now the initial iob rose stronger (probably because black beans have more carb than broccoli) but remained in the 3 U magnitude in first 2 hours after meal start



Despite aggressive settings that would quickly deliver big SMBs and exceed iobTH, we experience very moderate SMBs, and no hypo tendency at all.

0:15

Statistics

19/07	33.6 U	24.1 U	9.5 U	28%	0 g
20/07	38.8 U	23.6 U	15.2 U	39%	0 g
21/07	37.4 U	24.2 U	13.2 U	35%	0 g
22/07	34.1 U	19.3 U	14.8 U	43%	0 g
23/07	25.5 U	12.4 U	13.1 U	51%	0 g
24/07	30.8 U	18.1 U	12.7 U	41%	0 g

Average

07 days	34.2 U	21.1 U	13.1 U	38%	0 g
---------	--------	--------	--------	-----	-----

Today

25/07	0.0 U	0.0 U	0.0 U	0%	0 g
-------	-------	-------	-------	----	-----

RESET

TIR (70-180)

Date	Below	In range	Above
18/07	1%	99%	0%
19/07	7%	93%	0%
20/07	7%	87%	6%
21/07	2%	98%	0%
22/07	0%	100%	0%
23/07	0%	100%	0%
24/07	0%	100%	0%

Average (70-180)

07 days	3%	96%	1%
30 days	3%	95%	2%

Average (70-140)


07 days	3%	82%	16%
30 days	3%	83%	15%

On both days, 100% TIR was achieved:

94 Analysis using the Emulator

95 Let us have a look how the autoISF factors contributed to provide appropriate iob (well
96 below iobTH in this low carb case).

97 Following [section 10.2.1 – 10.2.4](#) of the FCL e-book, analysis of the 2nd meal (24th July):



98

99 Execute analysis / run emulation / inspect results / csv =>

100 In column C: AAPS time = UTC time (B) + 2/24 ; format without seconds

101 Hide not important lines and columns =>

	A	C	E	F	G	L	M	N	Q	R	S	V	W	X	Y	AA	AB	AC	AD	AE	AF	AH	AI	AK	AO	AQ	AR
1	2	bg		targ			eff.	tol.	final	dura	min-	parab	parat	parat	parab	acce	bg	pp	delta	dura	final	ISF	ISF	ISF	ISF	ISF	ISF
3	id	AAPS	brake	brake		iob	emul	emul	orig	utes	avg.	correl	durat	last-	next-	emul	emul	emul	emul	emul	emul	prof	emul	emul	emul	emul	emul
10	6	11:51	92	92	90	0,02	6	7,8	0,89	75	96,1	1	45	-1,3	-1,5	0,97	0,89	1	1	1	0,89	44	49,5	0	0	0,39	
11	7	11:56	91	91	90	0	6	7,8	0,87	80	95,8	1	25	-1,3	-1,6	0,97	0,87	1	1	1	0,87	44	50,3	0	0	0,385	
12	8	12:06	93	90	-0,03	4,92	6,4	0,74	90	95,6	0	0	0	0	0	1	0,9	1	1	1	0,74	43	56,1	0,14	0,1	0,731	
13	9	12:11	90	90	0,08	3,84	4,99	0,55	90	95,1	0	0	0	0	0	1	0,86	1	1	1	0,55	43	78,1	0	0	0	
14	10	12:16	87	90	0,03	3,84	4,99	0,54	5	88,5	0	0	0	0	0	1	0,84	1	1	1	0,54	43	80,1	0	0	0	
15	11	12:21	85	90	-0,02	3,84	4,99	0,53	10	87,3	0	0	0	0	0	1	0,82	1	1	1	0,53	43	81,4	0	0	0	
16	12	12:26	84	84	90	-0,07	4,92	6,4	0,67	10	85,3	0,9995	20	-1,2	-0,5	0,92	0,82	1	1	1	0,67	43	64,1	0	0	0	
17	13	12:31	82	82	90	-0,11	6	7,8	0,8	15	84,5	0,9973	25	-1,2	-0,8	0,95	0,8	1	1	1	0,8	43	53,5	0	0	0	
18	14	12:51	84	90	-0,19	6	7,8	0,87	35	83,5	0,9929	45	1,1	1,7	1,06	0,82	1	1	1	1	0,87	43	49,5	0,1	0	0,75	
19	15	12:56	85	90	-0,17	6	7,8	0,87	40	83,7	0,9976	25	1,5	2	1,05	0,82	1	1	1	1	0,87	43	49,5	0,1	0	0,77	
20	16	13:01	84	90	-0,14	6	7,8	0,82	45	83,7	0,8103	45	1,4	1,7	1	0,82	1	1	1	1	0,82	41	50,1	0	0	0,47	
21	17	13:06	83	85	90	-0,14	6	7,8	0,83	50	83,8	0,9998	30	0,5	0,4	0,99	0,82	1	1	1	0,82	41	49,7	0,02	0	0,495	
22	18	13:11	91	90	-0,14	6	7,8	1,15	0	91	0,9826	20	4	5,7	1,31	0,87	1	1	1	1	1,15	41	35,8	1,01	0,7	2,57	
23	19	13:16	96	90	0,74	6	7,8	1,3	0	96	1	25	5,7	7,4	1,38	0,94	1	1	1	1	1,3	41	31,6	0,79	0,5	0,22	
24	20	13:21	100	90	1,2	6	7,8	1,18	5	98	1	40	5	5,8	1,18	1	1,12	1	1	1	1,18	41	34,8	0,06	0	0,67	
25	21	13:26	102	102	90	1,19	6	7,8	0,9	10	99,3	0,9997	20	2,4	1,1	0,85	1	1,06	1	1,01	0,9	41	45,7	0	0	0	
26	22	13:31	104	104	90	1,12	6	7,8	0,92	10	102	0,9997	25	1,5	0,3	0,86	1	1,06	1	1,02	0,92	41	44,8	0	0	0	
27	23	13:36	108	108	90	1,03	6	7,8	1,08	15	104	0,9992	35	2,6	2,3	0,96	1	1,12	1	1,03	1,08	41	38	0,08	0	0,71	
28	24	13:41	110	110	90	1,01	6	7,8	1,02	15	106	1	40	2,3	2	0,96	1	1,06	1	1,04	1,02	41	40,2	0	0	0,51	
29	25	13:46	108	108	90	0,96	6,36	8,27	1,02	15	108	0,9876	40	0,5	-0,3	0,93	1	1	1	1,04	1,02	41	40,2	0	0	0	
30	26	13:51	108	108	90	0,88	6,36	8,27	1,02	20	108	0,9886	45	-0,4	-1,2	0,91	1	1	1	1,05	1,02	41	40,3	0	0	0	
31	27	13:56	109	109	90	0,77	6,36	8,27	1,06	25	108	0,9853	45	-0,6	-1,3	0,93	1	1,03	1	1,07	1,06	41	38,8	0	0	0,29	
32	28	14:01	112	90	0,7	6,36	8,27	1,43	30	108	0,9999	20	2,9	4,4	1,35	1	1,09	1	1,08	1,43	40	28,1	0,71	0,5	2,003		
33	29	14:06	116	90	1,27	6,36	8,27	1,39	30	110	1	20	4,1	5,6	1,31	1,01	1,12	1	1,09	1,39	40	28,7	0,42	0,3	1,98		
34	30	14:11	118	90	1,64	6,36	8,27	1,14	10	115	0,9974	40	2,6	3	1,07	1,01	1,06	1	1,04	1,14	40	35,2	0	0	0,023		
35	31	14:16	120	120	90	1,53	6,36	8,27	1,05	10	118	1	20	1,9	1,4	0,93	1,01	1,06	1	1,04	1,05	40	38,2	-0,03	0	0,4	
36	32	14:21	123	90	1,41	6,36	8,27	1,16	15	119	1	30	2,7	2,7	1	1,02	1,09	1	1,07	1,16	40	34,6	0,66	0,4	0,16		
37	33	14:26	125	125	90	1,71	6,36	8,27	1,13	20	120	0,9999	25	2,1	1,9	0,97	1,02	1,06	1	1,09	1,13	40	35,5	0,34	0,2	0	
38	34	14:31	127	127	90	1,78	6,36	8,27	1,14	20	123	1	30	2	1,8	0,98	1,02	1,06	1	1,1	1,14	40	35,2	0,48	0,3	0	
39	35	14:36	128	128	90	1,96	6,36	8,27	1,09	20	125	0,9999	20	1,1	0,6	0,93	1,02	1,03	1	1,1	1,09	40	36,7	0,03	0	0,643	
40	36	14:41	130	130	90	1,87	6,36	8,27	1,15	20	127	1	35	1,6	1,4	0,98	1,02	1,06	1	1,11	1,15	40	34,7	0,89	0,6	0	

102

	A	C	E	F	G	L	M	N	Q	R	S	V	W	X	Y	AA	AB	AC	AD	AE	AF	AH	AI	AK	AO	AQ	AR
->	14:46	130	130	90	2,35	6,36	8,27	1,15	25	127	0,9981	30	0,7	0,3	0,95	1,02	1	1	1,14	1,15	40	34,8	0	0	0	0	
	14:51	133	133	90	2,21	6	7,8	1,15	30	128	0,9931	45	1,2	1	0,98	1,03	1,09	1	1,17	1,15	40	34,8	0,52	0,3	0	0	
	14:56	135		90	2,36	6	7,8	1,18	30	130	0,9999	20	2,3	2,8	1,09	1,03	1,06	1	1,18	1,18	40	34	0,91	0,6	0	0	
44	15:01	137		90	2,79	6	7,8	1,21	35	131	1	35	2,1	2,3	1,03	1,03	1,06	1	1,21	1,21	40	33	0,67	0,5	0	0	
45	15:06	136	136	90	3,12	6	7,8	1,09	35	132	0,9796	20	-0,1	-1,3	0,89	1,03	1	1	1,22	1,09	40	36,9	0	0	0	0	
46	15:11	135	135	90	2,94	6	7,8	1,08	40	132	1	25	-1,3	-2,4	0,86	1,03	1	1	1,25	1,08	40	37,1	0	0	0	0	
47	15:16	136	136	90	2,76	6	7,8	1,24	45	133	0,9991	45	0	-0,2	0,97	1,03	1,03	1	1,28	1,24	40	32,2	0	0	0	0	
48	15:26	133	133	90	2,39	6,36	8,27	1,35	55	133	0,9983	35	-1,5	-2	0,94	1,03	1	1	1,35	1,35	40	29,7	0	0	0	0	
49	15:30	132	132	90	2,21	6,36	8,27	1,43	60	133	0,9994	30	-1,3	-1,5	0,97	1,03	1	1	1,38	1,43	40	28	0	0	0,337	0	
50	15:36	131		90	2,05	6,36	8,27	1,15	65	133	0,9995	20	-1,1	-0,9	1,03	1,03	1	1	1,41	1,5	40	26,7	-0,04	0	0	0	
51	15:41	130		90	1,88	6,36	8,27	1,53	70	133	1	25	-0,9	-0,8	1,03	1,02	1	1	1,44	1,53	40	26,2	0,11	0	0,697	0	
->	15:46	130		90	1,82	6,36	8,27	1,56	75	132	0,9985	25	-0,3	0,1	1,08	1,02	1	1	1,47	1,56	40	25,7	0,27	0,2	1,017	0	
	15:51	130		90	1,94	6,36	8,27	1,59	80	132	1	30	0,1	0,4	1,08	1,02	1	1	1,5	1,59	40	25,2	0,28	0,2	0,08	0	
	16:01	127	127	90	1,83	6	7,8	1,5	90	132	1	25	-1,3	-1,6	0,97	1,02	1	1	1,56	1,5	39	26	-0,04	0	0,46	0	
55	16:06	130		90	1,71	6	7,8	1,59	95	132	0,8717	45	0,2	0,4	1	1,02	1,09	1	1,59	1,59	39	24,6	0,93	0,6	2,34	0	
56	16:11	131		90	2,45	6	7,8	1,62	100	132	0,9959	25	1,9	2,8	1,13	1,03	1,03	1	1,62	1,62	39	24,1	0,12	0	0,72	0	
57	16:16	132		90	2,35	6	7,8	1,65	105	132	1	40	1,4	1,8	1,09	1,03	1,03	1	1,65	1,65	39	23,7	0,34	0,2	0,29	0	
58	16:21	131	131	90	2,4	6	7,8	1,46	110	132	0,9949	20	-0,7	-1,9	0,87	1,03	1	1	1,68	1,46	39	26,7	0	0	0,1	0	
59	16:26	134		90	2,24	6	7,8	1,71	115	132	0,9998	40	1,7	2	1,07	1,03	1,09	1	1,71	1,71	39	22,8	1,14	0,8	2,76	4,9	
60	16:31	136		90	3,07	6	7,8	1,74	120	132	1	45	2	2,4	1,07	1,03	1,06	1	1,74	1,74	39	22,4	1,29	0,9	0	0	
61	16:36	136	136	90	3,78	6	7,8	1,77	125	132	0,9971	35	1,1	1,1	1	1,03	1	1	1,78	1,77	39	22	0	0	0	0	
->	16:41	136	136	90	3,58	6	7,8	1,59	130	132	0,9896	20	-0,5	-1,7	0,88	1,03	1	1	1,81	1,59	39	24,5	-0,41	0	0	0	
	16:51	132	132	90	3,15	6,36	8,27	1,5	140	132	0,9996	20	-2,3	-3,1	0,9	1,02	1	1	1,68	1,5	39	25,9	0	0	0	0	
	16:56	129	129	90	2,93	6,36	8,27	1,53	145	132	1	25	-3	-3,9	0,9	1,01	1	1	1,7	1,53	39	25,5	0	0	0	0	
65	17:01	123	123	90	2,72	6,36	8,27	0,9	15	130	0,9984	20	-5	-6,3	0,85	1,01	1	1	1,07	0,9	38	42	0	0	0	0	
66	17:06	114	114	90	2,5	6,36	8,27	0,71	0	114	0,9998	20	-8,5	-11	0,71	1	1	1	1	0,71	38	53,6	0	0	0	0	
67	17:11	108	108	90	2,29	6,36	8,27	0,85	0	108	0,9997	45	-7,9	-9,2	0,85	1	1	1	1	0,85	38	44,7	0	0	0	0	
68	17:16	101		90	2,08	6,36	8,27	0,94	0	101	1	20	-6,9	-6,7	1,03	0,92	1	1	1	0,94	38	40,2	0	0	0	0	
69	17:21	93	93	90	1,89	6,36	8,27	0,83	0	93	0,9999	25	-7,3	-7,4	0,98	0,83	1	1	1	0,83	38	45,7	0	0	0	0	
70	17:26	91	91	90	1,69	6,36	8,27	0,82	5	92	0,9966	25	-4,6	-3,7	0,89	0,82	1	1	1	0,82	38	46,5	0	0	0	0	
71	17:31	91		90	1,52	6,36	8,27	0,96	10	91,7	0,9997	20	-0,1	2,7	1,17	0,82	1	1	1	0,96	38	39,7	0	0	0	13:	
72		82							0,53						0,71	0,8	1	1	1		36	22	0	0	0	to	
93		137							1,77						1,57	1,03	1,27	1	1,81		44	81,4		1,1	16:	to	
94																								10,8	2,82		

105 Discussion of results

106 Optimization potential for low carb meals?

107 The investigated low carb meals were managed very well with the settings that also suit
108 higher carb meals.

109 *In case there were ideas* for potential further improvements:

- 110 • the difference that modified settings *would make* could be analyzed using the
111 Emulator with a “what-if” vdf. See FCL e-book [section 10.3](#) for analysis on your
112 PC, and (especially neat for AAPS:) the *real-time* emulator analysis of a *what-if*
113 question, with speech synthesis on your smartphone, telling you for every loop
114 decision how it would be impacted. See FCL e-book [section 11.4](#)
- 115 • it must be critically checked whether, after such optimization, *other* meals in your
116 spectrum then might suffer. It is problematic, to fine tune just for one kind of
117 meal.(See FCL e-book [section 8](#), and [case study 8.2](#)),
- 118 • the user could also choose to create 2 or more differentiated sets of parameters, suited
119 to different eating habits at lunch vs at dinner, for instance, and optimize both
120 independently.

121

122 Was I just lucky – Could there be trickier cases?

123 We know that super big meals do not pose much extra difficulty, if the system was set up for
124 a diet spectrum that included meals with max carb “burn rate” of around 30g/h (ref: Chapter 8
125 in: <https://github.com/danamlewis/artificialpancreasbook/> -)

126 The presented case study showed that low carb meals are also easily managed.

127 **Reduced-size meals (or high carb snacks) that trigger really aggressive first SMBs due**
128 **to** (a limited amount of) **fast resorbing carbs would be more problematic**, and might
129 require a manual intervention. - See [case study 5.2](#) re. sweet snacks.